

Keysight Technologies 8500 FE-SEM

Compact System for Low-Voltage, High-Performance Imaging

Data Sheet



Features and Benefits

- Resolution and imaging equal to that of conventional FE-SEMs
- Variable low voltage eliminates charging and the need for sample coating
- Programmable X, Y, Z stage allows user to set precise coordinates, scan, and save information
- Miniature electrostatic lens design ensures repeatable performance without constant re-tuning
- Compact size enables easy installation in any research laboratory and does not require special facilities

Overview

The new Keysight Technologies, Inc. 8500 FE-SEM offers researchers a field emission scanning electron microscope (FE-SEM) right in their own laboratory. This compact, innovative system has been optimized for low-voltage imaging, extremely high surface contrast, and resolution typically found only in much larger and more expensive field emission microscopes.

The 8500 is easy to install and use. No dedicated facilities are required, only an AC power outlet. About the size of a laser printer, the scientific-grade system provides capabilities to researchers in ordinary labs that previously were only available with conventional FE-SEMs installed in centralized facilities. The 8500 has been engineered to deliver consistent, reproducible performance and the industry's lowest total cost of ownership for a FE-SEM.

Enhanced Imaging Capabilities

The 8500 FE-SEM offers several imaging techniques for enhancing surface contrast and allowing nanoscale features to be observed on a wide variety of nanostructured materials, including polymers, thin films, biomaterials, and other energy-sensitive samples on any substrate, even glass.

The system's continuously variable imaging voltage is tunable from 500 to 2000 volts as an operational parameter, rather than a setup choice. The 8500 eliminates charging of nonconductive samples without the need to coat the samples, which can mask nanoscale features, or resort to increased pressure operation, which can degrade resolution.

In addition, the 8500 utilizes a four-segment microchannel plate (MCP) detector that provides topographic imaging along two orthogonal directions to enhance surface detail. This technique has been demonstrated to clearly resolve sub-nanometer atomic steps on the surface of crystalline substances such as polytype 6H-SiC

Applications

- Polymers
- Thin films
- Biomaterials
- Nonconductive samples
- Energy-sensitive materials
- Glass substrates

Novel Design

The core technology inside a scanning electron microscope is the electron beam column, which extracts, collimates, shapes, scans, and focuses the electron beam. A conventional electron beam column relies on a combination of precision-machined electromagnetic and electrostatic elements to control the electron beam. The coils that form the critical elements are typically hand-wrapped by technicians to achieve uniform electromagnetic fields. Moreover, closed-loop cooling and sophisticated vibration isolation are often required at the system level to manage the high currents in the lenses and other elements. The resultant high-resolution electron beam column is both large and expensive.

The system's thermal field emission electron source, meanwhile, provides high brightness, high stability, small virtual source size, and low energy spread and long-lasting consistent performance. A quad-segmented MCP detector is located just below the objective lens of the electron beam column directly and above the sample. This detector collects both backscattered and secondary electrons. The MCP may be operated either in a standard mode, where all the channels are added together, or in a differential mode (topographic mode), in which opposite sides of the detector are dynamically subtracted.

Performance Advantages


As described above, the 8500 FE-SEM embodies the successful miniaturization of the core technology in a scanning electron microscope. The 8500 is optimized for low-voltage imaging and sub-10nm resolution. Its thermal field emission electron source provides high signal-to-noise ratios and consistent, long-lasting performance, while secondary and backscatter electron detection capabilities provide a rich data set for each sample.

Ultimate Ease of Use

Just as the factory-calibrated 8500 FE-SEM is approximately the same size as a laser printer, it also offers similar plug-and-play performance. In fact, the compact size of the 8500 facilitates easy installation in practically any existing lab that has an AC power outlet. Sample preparation and loading is both simple and fast. An X, Y, Z programmable stage lets users set specific coordinates, scan and then store the locations to repeat experiments with precision and confidence. For optimum control, the system's powerful software package features an intuitive graphical user interface (GUI) designed for novice and expert users.

ECD Cartridge

The electron source, the electron beam column, and the electron detector of the 8500 FE-SEM are all combined in a field-replaceable ECD cartridge. When the long-lasting electron source is finally depleted, the entire ECD cartridge can be replaced on-site to provide the 8500 with not only a new source of electrons, but a new pre-aligned electron beam column and a new MCP detector – essentially a brand new FE-SEM.



Schematic of electrostatic electron beam column.




Figure 1. High resolution image of gold islands on carbon.




Figure 2. Low voltage image of uncoated electrode with CNTs on tip.




Figure 3. High resolution image of celtgard.




Figure 4. High resolution topographical surface information of de-processed chip metal conductor lines.

8500 Specifications

Performance

Beam voltage	500 to 2000V
Beam current	0.2 to 1 nA
Resolution	<10nm at 1000V
Magnification*	250 to 65,000 X
Digital zoom	10X
Scan field	1 x 1 mm (max)
Electron source	Schottky field emission
Detector modes	SE, BSE, Topo

Sample

Sample size	100 x 60mm (max)
Sample thickness	30mm (max)
Viewable area	50 x 30mm (max)
Sample mounts	Standard SEM stubs
Electrical activation	Vacuum feedthrough

System Control

PC	Windows 7
Basic user interface	Simple image capture and controls
Expert user interface	Expert image optimization

Image

Image formats	JPEG, TIFF, BMP, PNG
Image resolution	User selectable up to 2048 x 2048 pixels
Scan rate	Slow scan to video rate
Noise reduction	Frame and pixel averaging

Motorized Stage

Piezo stage	1 µm accuracy
X, Y, Z travel	50 x 50 x 10mm

Vacuum System

Chamber vacuum	1 e ⁻⁴ Torr
Pumpdown time	3 minutes
Turbo pump	80 liters per second
UHV pump	Ion pump with gettering

Dimensions

Microscope	584(W) x 470(D) x 584(H)mm; 72 kg
Pump unit	203(W) x 254(D) x 203(H)mm; 4 kg

Installation Requirements

Power	100/120/220-240VAC; 50/60Hz
Operating temperature	5 to 40°C
Humidity	20 to 80% RH
Compressed air	Not required
Dry nitrogen	Not required**
Water cooling	Not required

* Electron optical magnification relative to a 3nm pixel on a VGA monitor.

** A N₂ vent port is configured and available.

FE-SEM Instrumentation from Keysight Technologies

Keysight Technologies offers high-precision, compact FE-SEMs for research, industry, and education. Exceptional worldwide support is provided by experienced application scientists and technical service personnel. Keysight's leading-edge R&D laboratories are dedicated to the timely introduction and optimization of innovative and easy-to-use technologies.

www.keysight.com/find/fe-sem

Nanomeasurement Systems from Keysight Technologies

Keysight Technologies, the premier measurement company, offers high-precision, modular nanomeasurement solutions for research, industry, and education. Exceptional worldwide support is provided by experienced application scientists and technical service personnel. Keysight's leading-edge R&D laboratories ensure the continued, timely introduction and optimization of innovative, easy-to-use nanomeasure system technologies.

For more information on Keysight Technologies' products, applications or services, please contact your local Keysight office. The complete list is available at: www.keysight.com/find/contactus

Americas

Canada	(877) 894 4414
Brazil	55 11 3351 7010
Mexico	001 800 254 2440
United States	(800) 829 4444

Asia Pacific

Australia	1 800 629 485
China	800 810 0189
Hong Kong	800 938 693
India	1 800 112 929
Japan	0120 (421) 345
Korea	080 769 0800
Malaysia	1 800 888 848
Singapore	1 800 375 8100
Taiwan	0800 047 866
Other AP Countries	(65) 6375 8100

Europe & Middle East

Austria	0800 001122
Belgium	0800 58580
Finland	0800 523252
France	0805 980333
Germany	0800 6270999
Ireland	1800 832700
Israel	1 809 343051
Italy	800 599100
Luxembourg	+32 800 58580
Netherlands	0800 0233200
Russia	8800 5009286
Spain	0800 000154
Sweden	0200 882255
Switzerland	0800 805353
	Opt. 1 (DE)
	Opt. 2 (FR)
	Opt. 3 (IT)
United Kingdom	0800 0260637

For other unlisted countries:

www.keysight.com/find/contactus

(BP-05-19-14)

This information is subject to change without notice.

© Keysight Technologies, 2011 - 2014

Published in USA, August 1, 2014

5990-6183EN

www.keysight.com