

860 KF Thermoprep



Manual
8.860.8003EN



Metrohm AG

CH-9101 Herisau

Switzerland

Phone +41 71 353 85 85

Fax +41 71 353 89 01

info@metrohm.com

www.metrohm.com

860 KF Thermoprep

Manual

Teachware
Metrohm AG
CH-9101 Herisau
teachware@metrohm.com

This documentation is protected by copyright. All rights reserved.

Although all the information given in this documentation has been checked with great care, errors cannot be entirely excluded. Should you notice any mistakes please send us your comments using the address given above.

Table of contents

1	Introduction	1
1.1	Instrument description	1
1.1.1	Instrument components	1
1.1.2	Intended use	2
1.2	About the documentation	2
1.2.1	Symbols and conventions	2
1.3	Safety instructions	3
1.3.1	General notes on safety	3
1.3.2	Electrical safety	3
1.3.3	Personnel safety	4
1.3.4	Flammable samples and solvents	5
1.3.5	Recycling and disposal	5
2	Overview of the instrument	6
3	Installation	8
3.1	Setting up the instrument	8
3.1.1	Packaging	8
3.1.2	Checks	8
3.1.3	Location	8
3.2	Connecting the mains cable	8
3.3	Mounting the guidance rod	9
3.4	Mounting the sample insert	12
3.5	Mounting the needles	13
3.6	Mounting the heating tubing	14
3.7	Mounting the dust filter	16
3.8	Assembling the drying flask	16
3.9	Air/nitrogen connector	18
3.10	Inserting the heating tube into the KF titration cell	19
4	Operation	22
4.1	Switching the instrument on and off	22
4.2	The keypad	22
4.3	The display	23
4.3.1	Display elements	23
4.3.2	Status display	23
4.3.3	Error messages	24



10.3	Quality Management Principles	46
11	Accessories	48
11.1	Scope of delivery 2.860.0010	48
11.2	Optional accessories	55
	Index	56

1 Introduction

1.1 Instrument description

The 860 KF Thermoprep is used whenever the heating up of a sample and/or the thermal expulsion of moisture in solid substances or liquids is required. In combination with a coulometric or volumetric KF titrator, the 860 KF Thermoprep is the ideal analysis system for water determination in samples that contain disruptive components or from which moisture can be removed only with difficulty.

One of its decisive advantages is the reduction of sample preparation to a minimum. Thanks to the use of hermetically sealed sample vessels ("head-space vials"), the filling of the samples can be accomplished directly on-site. The PTFE-coated septa guarantee a constant, non-falsified water content, even after prolonged holding times.

The sample heated in the oven module releases its moisture in the form of water vapor, which is conveyed into a measuring cell with the aid of a gas flow. An air pump is installed for the purpose of generating the gas flow. An inlet valve is available for nitrogen or other inert gases. The determination of the moisture can be accomplished in the measuring cell either coulometrically or volumetrically in accordance with Karl Fischer.

1.1.1 Instrument components

The 860 KF Thermoprep has the following components:

- **Oven**
Oven module made of aluminum with software-operated temperature control for heating the sample vessel.
- **Fan**
Propeller fan for cooling the oven module.
- **Inlet valve**
Valve for switching over the source of the gas flow.
- **Air pump**
Pump for generating the gas flow.
- **Outlet heating**
Heating tube for preventing the condensation of moisture.
- **Stand rods**
Guidance device with needle adapter and tubing for the gas flow.
- **Operating unit**
Monochrome LCD display and keyboard.

1.1.2 Intended use

The 860 KF Thermoprep is designed for usage as an auxiliary device for sample preparation in analytical laboratories. Its main area of application is moisture determination according to Karl Fischer (coulometric or volumetric). The 860 KF Thermoprep enables the application of thermal gas extraction technology.

The present instrument is suitable for processing chemicals and flammable samples. The usage of the 860 KF Thermoprep therefore requires that the user has basic knowledge and experience in the handling of poisonous and caustic substances. Knowledge with respect to the application of the fire prevention measures prescribed for laboratories is also mandatory.

1.2 About the documentation









Caution

Please read through this documentation carefully before putting the instrument into operation. The documentation contains information and warnings which have to be followed by the user in order to ensure safe operation of the instrument.

1.2.1 Symbols and conventions

The following symbols and styles are used in this documentation:

(5-12)	Cross-reference to figure legend The first number refers to the figure number, the second to the instrument part in the figure.
1	Instruction step Carry out these steps in the sequence shown.
Method	Dialog text, parameter in the software
File ► New	Menu or menu item
[Next]	Button or key
	Warning This symbol draws attention to a possible life hazard or risk of injury.
	Warning This symbol draws attention to a possible hazard due to electrical current.

	<p>Warning</p> <p>This symbol draws attention to a possible hazard due to heat or hot instrument parts.</p>
	<p>Warning</p> <p>This symbol draws attention to a possible biological hazard.</p>
	<p>Caution</p> <p>This symbol draws attention to a possible damage of instruments or instrument parts.</p>
	<p>Note</p> <p>This symbol marks additional information and tips.</p>

1.3 Safety instructions

1.3.1 General notes on safety



Warning

This instrument may only be operated in accordance with the specifications in this documentation.

This instrument has left the factory in a flawless state in terms of technical safety. To maintain this state and ensure non-hazardous operation of the instrument, the following instructions must be observed carefully.

1.3.2 Electrical safety

The electrical safety when working with the instrument is ensured as part of the international standard IEC 61010.



Warning

Only personnel qualified by Metrohm are authorized to carry out service work on electronic components.



Warning

Never open the housing of the instrument. The instrument could be damaged by this. There is also a risk of serious injury if live components are touched.

There are no parts inside the housing which can be serviced or replaced by the user.

Mains voltage



Warning

An incorrect mains voltage can damage the instrument.

Only operate this instrument with a mains voltage specified for it (see rear panel of the instrument).

Protection against electrostatic charges



Warning

Electronic components are sensitive to electrostatic charges and can be destroyed by discharges.

Always pull the mains cable out of the mains connection socket before connecting or disconnecting electrical appliances on the rear panel of the instrument.

1.3.3 Personnel safety



Warning

Wear protective goggles and working clothes suitable for laboratory work while operating the 860 KF Thermoprep.



Warning

The **oven and the sample vessels** can be at a temperature of up to **250 °C**! Avoid direct skin contact. Use the septum-sealing tongs supplied or another suitable holding device for inserting the sample vessels in the oven or removing them. Wear heat-insulated gloves if necessary.



Warning

There is a **considerable risk of injury** connected with the needles. Grasp the handle of the needle adapter with both hands when you guide the needle into the sample vessel. Take care to ensure that no other person reaches into the danger zone during this procedure.

1.3.4 Flammable samples and solvents




Warning

All relevant safety measures are to be observed when working with flammable samples or when adding solvents.

- Use nitrogen or another inert gas for conveying the moisture.
- Set up the instrument in a well-ventilated location.
- Keep all sources of flame far from the workplace.
- Clean up spilled fluids and solids immediately.
- Follow the safety instructions of the chemical manufacturer.

1.3.5 Recycling and disposal



This product is covered by European Directive 2002/96/EC, WEEE – Waste from Electrical and Electronic Equipment.

The correct disposal of your old equipment will help to prevent negative effects on the environment and public health.

More details about the disposal of your old equipment can be obtained from your local authorities, from waste disposal companies or from your local dealer.

2 Overview of the instrument




Figure 1 Front 860 KF Thermoprep

1	Guide head	2	Needle holder (6.2049.040) With needle adapter.
3	Spring sleeve With vessel stripper.	4	Guidance rod
5	Sample vessel (6.2419.007) With 6 ml content.	6	Sample insert (6.2063.010)
7	Oven module With sheet metal housing.	8	Fixing lever
9	Operating unit With display and keyboard.	10	Drying flasks (6.1608.050) With 6.1602.145 drying flask insert.
11	Drying flask holder	12	Gas inlet With M6 thread.
13	Gas outlet With M6 thread.		




Figure 2 Rear 860 KF Thermoprep

1 Air pump inlet With 6.2724.010 dust filter.	2 Mains switch
3 Mains connection socket	4 Type plate Contains specifications concerning mains voltage and serial number.
5 Outlet heating connection For the heating tube.	6 Air/nitrogen connector With M6 interior thread. Inlet for external gassing.
7 Fan For cooling the oven module.	8 Cable holder (6.2060.010)

3.3 Mounting the guidance rod

The guidance rod with the guide head is enclosed separately with the 860 KF Thermoprep and must first be mounted.

Proceed as follows:

- 1 Move the guide head of the guidance rod a bit downwards by pressing the fixing lever.




Figure 4 Loosen the adapter

- 2** Loosen the adapter on the head of the guidance rod with a hexagon key.
- 3** Remove the adapter.
- 4** Insert the tubing jutting out of the housing from below into the guidance rod.



Note

This tubing conveys the gas from the flow controller upward through the guidance rod to the guide head. From there the gas is guided to the double hollow needle, see next chapter.



-
- Diagram illustrating the removal of the top cap (6) and the insertion of the top cap (7) into the top of the unit.

- 7** Fasten the tubing on the olive of the 4.860.4260 adapter.
- 8** Pull out the guidance rod as far as possible in order that the adapter rests close on the guidance rod.



-

11 Fix the guidance rod with the V.0434.008 hexagon screw according to the drawing. The screw and the associated hexagon key are packaged separately with the instrument.

3.4 Mounting the sample insert




Figure 9 Mounting the sample insert



Note

The sample insert may not be inserted or removed unless it has been cooled down.

The dimensions of the supplied sample insert are optimized for the usage of 6 mL sample vials 6.2419.007 from Metrohm. This ensures the best possible transfer of heat between oven and sample.

- 1 Place the **6.2627.000 oven insert** into the oven from above.
- 2 Place the **6.2063.010 sample insert** into the oven from above.

If vessels with other dimensions are to be used, then individually modified sample inserts can be ordered. The precise vessel dimensions (including tolerances) will be required when ordering. Non-standard sample vessels may not exhibit dimensions outside of the following limit values:

- Diameter 10.0...32.0 mm
- Immersion depth 20.0...45.0 mm

3.5 Mounting the needles

The length of the needle holder defines how deeply the injection needle penetrates the sample vessel. The **6.2049.040 needle holder**, which is 58 mm in length, ensures that the needle penetrates the liquid or powdery sample. The carrier gas can flow through the sample and effect an efficient expulsion of the moisture it contains.

If there is a danger that the heated sample could clog the needle, then use the **6.2049.050 needle holder** with 73 mm length. In this case the injection needle penetrates the sample vessel only slightly deeper than the outlet needle and has no contact with the sample itself. The 6.2049.050 needle holder can be ordered from Metrohm if required.

Needle holders with the dimensions required for situations calling for special sample vessels can be supplied by Metrohm upon request.




Figure 10 Mounting the needles

Mount the needles as follows

- 1** Screw the **6.2049.040 needle holder** into the distributor on the guide head.

- 2 Screw the **6.2816.080 outlet needle** onto the Luer connector of the needle holder.
- 3 Carefully guide the **6.2816.070 injection needle** into the opening of the distributor from above and allow it to drop down.



Note

Take care to ensure that the white PTFE seal is positioned securely on the needle.

- 4** Screw the **6.1805.100 FEP tubing** by hand onto the opening of the distributor. Screw the other end of the tubing into the opening on the upper side of the guidance rod.

3.6 Mounting the heating tubing





Figure 11 Mounting the heating tubing (view from rear)

Proceed as follows:

- 1 Screw the M6 connector of the **6.1830.030 heating tubing** into the side opening of the distributor on the guide head.

- 

- 
- Air/N₂
max. 0.3
Inlet heater
4

Rotate the plug in such a way that the three contact pins match the alignment of the corresponding openings on the socket. Press the plug against the socket and rotate the front knurled screw in clockwise direction.



- 2 Screw one **6.1821.040 filter tube** into each of the **6.1602.145 drying flask covers** from below. Tighten the filter tubes well by hand.
- 3 Screw the two drying flask covers with the filter tubes onto the drying flasks. Tighten the covers well by hand.



Note

If drying flask covers or filter tubes are not sufficiently tightly screwed on, then this will prevent a precise, regular flow of gas. The error message "Flow rate error" will appear as a rule when there are leaks in the threaded connections.




Figure 15 Mounting the tubings

Mount the FEP tubings as follows:

- 1 Load the two drying flasks that have been prepared into the holders, see illustration.
- 2 Screw one 6.1805.080 FEP tubing (25 cm length) to the gas outlet (at the front, labeled **1**). Screw the other tubing end on the front drying flask to the M6 connector **without** point marking (see left-hand arrow).

- 3 Screw the second 6.1805.080 FEP tubing (25 cm length) to the gas inlet (at the rear, labeled **2**). Screw the other tubing end on the rear drying flask to the M6 connector **with** point marking (see right-hand arrow).
- 4 Screw the 6.1805.010 FEP tubing (13 cm length) tightly onto the remaining M6 connectors of the drying flasks.

The marking **Drying flask** on the right-hand side of the instrument displays the diagram for the tubing.



Note

Tighten the screw connections well by hand.

3.9 Air/nitrogen connector

If compressed air, nitrogen or another gas is to be used for transferring or expelling moisture, then a separate connector is available on the rear of the instrument.

A tube with M6 thread can be connected directly to the connector **Air/N₂** in. Enclosed with the instrument is the **6.1808.040 M6/M8 tubing adapter** for tubing with an M8 thread. The **6.1808.050 M8/tubing olive** can also be put in place in order to connect a simple tubing.




Figure 16 External gas supply connection



Warning

If gas is supplied from a pressure line or a pressure vessel, then it is imperative that a pressure reduction valve be placed upstream. The gas pressure may not exceed a maximum overpressure level of 0.3 bar.

3.10 Inserting the heating tube into the KF titration cell

The 860 KF Thermoprep can be used for either coulometric or volumetric water determinations. The mounting of the tip of the heating tube is accomplished in different ways for the respective KF titration cells.

Coulometric KF titration cell

- Remove the protective cover from the tip of the heating tube and the E.3010.032 O-ring.
- Disassemble the accompanying **6.1446.170 heating tube stopper** into three parts.




Figure 17 Coulometric KF titration cell

Proceed as follows:

- 1 Guide the upper part of the heating tube stopper over the tip of the **6.1830.030 heating tube** as shown in the illustration.
- 2 Guide the O-ring of the heating tube stopper over the heating tube.
- 3 Guide the lower part of the heating tube stopper over the heating tube and screw the three parts together (not too tightly).
- 4 Insert the stopper with the tip of the heating tube into one of the two SGJ openings of the KF titration cell

- 5 Shift the tip vertically in such a way that the outlet opening of the tubing is immersed as deeply as possible. The tip must not however be permitted to get in the way of the stirring bar in the KF titration cell. Afterwards give the heating tube stopper its final tightening.

Volumetric KF titration cell

- Remove the protective cover from the tip of the heating tube and the E.3010.032 O-ring. The latter you will still need.
- A **6.2730.020 septum stopper** is enclosed with every KF titrator made by Metrohm. Disassemble this septum stopper into three parts and remove the septum. It will not be required.




Figure 18 Volumetric KF titration cell

Proceed as follows:

- 1 Guide the upper part of the septum stopper over the tip of the **6.1830.030 heating tube** as shown in the illustration.
- 2 Guide the **E.3010.032 O-ring** over the heating tube.

This O-ring is also part of the 6.1244.040 set of seals that is enclosed with each KF titrator. It cannot be reordered individually.
- 3 Guide the lower part of the septum stopper over the heating tube and screw the three parts together (not too tightly).
- 4 Insert the stopper with the tip of the heating tube into the front opening of the KF titration cell

- 5 Shift the height of the tip of the heating tube vertically in such a way that the outlet opening of the tubing is immersed as deeply as possible. The tip must not however be permitted to get in the way of the stirring bar in the KF titration cell. Afterwards give the septum stopper its final tightening.

- 【◀】 Reduces gas flow by 1 mL/min during ongoing operation. Keep the key pressed down briefly for this purpose.
- 【▶】 Increases gas flow by 1 mL/min during ongoing operation. Keep the key pressed down briefly for this purpose.
- 【▲】 Reduces the brightness of the display.
- 【▼】 Increases the brightness of the display.

4.3 The display

4.3.1 Display elements







Figure 21 Main dialog

1 Temperature display Shows the measured oven temperature.	2 Symbol of the status display
3 Gas flow display Displays the set gas supply and the flow rate.	4 Target temperature Displays the set specification (Temp. 1 or Temp. 2).
5 Status display Shows status and error messages.	

4.3.2 Status display

Table 1 Symbols of the status display

Symbol	Significance
	The oven is being heated.
	The oven is being cooled.
	The pre-selected temperature was reached.

4.4 Modifying numerical values and settings

The instrument settings can be modified in two different ways:

- Selection of specified values from a list. Parameters with list boxes are indicated with a colon (:).
- Modifying numerical values.

The different parameters can be selected with the arrow keys [▼] and [▲]. The current line is displayed in inverted fashion.

4.4.1 Selecting from a selection list

To modify parameters with list boxes (e.g. **Flow source:**), use the arrow keys [▶] or [◀]. The respective next value in the list box will be displayed.

Confirm a modified value with the [OK] key.

4.4.2 Modifying numerical values

Modify parameters with numerical values using the arrow keys:

- 1 Use the arrow keys [▶] or [◀] to select one of the digits of the numerical value.

The selected digit is displayed black on white. The rest of the line remains inverted.

Temp. 1 60 °C

- 2 Use the arrow keys [▲] or [▼] to modify the digits.

- 3 Select and modify the next digit, etc.

- 4 Apply the modified value with the [OK] key.

Press [QUIT] if the previous value is to be restored.

[OK]



Range	0.001 ... 2.000 (Increment: 0.001)
Default value	1.000

Gas type	Measuring factor
Argon	0.950
Helium	0.300
Oxygen	1.000
Air	1.000
Nitrogen	1.000

Selecting the flow source.

Selection	pump valve
Default value	pump

Built-in air pump.

Inlet valve for inert gas.

[TEMP]

The basic settings for the instrument can be made in the Setup menu. You can call it up during switch on by holding down the **[TEMP]** key. Press the **[QUIT]** key to exit the menu.



Note

The Setup menu can also be called up when the keylock is switched on.

Keylock. The **[OK]**, **[←]** and **[→]** keys can be locked. You can thus prevent changes from being made in the Parameters menu or of the gas flow rate during a determination. If the function is switched on (**on**), then the only switching that can be done is between Temperature 1 and Temperature 2,



in addition to switching the gas flow on and off. Furthermore it is still possible to modify the brightness of the display.

Selection	off on
Default value	off

Temp.corr

Temperature correction. The target value for the temperature control of the oven can be influenced with this correction value. A temperature differential between the oven module and the sample in the vessel can thus be brought into line as needed.

Range	-10 ... 10 °C (Increment: 1)
Default value	0 °C

The determination of the necessary temperature correction can be accomplished with a special oven insert and should be performed by a service technician. Contact your responsible Metrohm supplier.

5 Performing a determination

Moisture determinations according to the heating method require the conditioning of the entire tubing system and of the KF titration cell before a determination can be carried out. Because the sample vials and the septum seals may contain small amounts of moisture, three to five blank value determinations must be carried out with sealed, empty vessels.

Solid or liquid samples are weighed out into sample vials, which are then sealed. We recommend the use of 6.2419.007 sample vials with 6.1448.057 septum seals. The silicone septum of the aluminum cap resists temperatures of up to 250 °C and has proven itself to be an optimum solution.

5.1 Conditioning the system

The entire system must be conditioned prior to use, i. e. it must be absolutely water-free. The oven must be heated to the desired working temperature for this purpose. The needle must be inserted into a sealed sample vessel (conditioning vessel) in order to flush the entire tubing system with the carrier gas. The gas supply is thus connected with the heating tube.

5.1.1 Preparing the KF titration cell

Filling the KF titration cell

The KF titration cell must remain filled with conditioned working medium when not in use. The working medium must be replaced from time to time, even though the sample is not transferred into the working medium during the heating-out process.

Prepare the KF titration cell as follows:

- 1** Fill the KF titration cell with working medium. We recommend approximately 150 mL of working medium when using a Coulometer cell and approximately 35 mL with a volumetric KF titration cell. The shaft of the heating tube should, however, not be immersed.
- 2** Align the tip of the heating tube against the vessel wall.
- 3** Switch on the stirrer. Set the stirring rate to high.

It is advantageous to have air bubbles stirred into the working medium and forcefully broken up in the process. The air in the KF titration cell must also be dried.

Conditioning the working medium

- 1 Load a conditioning method on the KF titrator and begin the conditioning by pressing the **[Start]** key.

You will find the necessary settings in the KF titrator manual.



Note

It is imperative that an **extraction time of at least 300 seconds** be set for a determination using the heating method.

5.1.2 Preparing the 860 KF Thermoprep

Switch on the 860 KF Thermoprep and select the temperature

The necessary settings on the 860 KF Thermoprep:

- 1 Switch on the instrument.
- 2 Select one of the preset temperatures using the **[TEMP]** key.
If you would like to change **Temp. 1** or **Temp. 2**, press **[OK]**.

Inserting the conditioning vessel

- 1 Seal an empty sample vessel tightly with a septum seal. To do this, use the crimping tongs.
- 2 Insert the sample vessel into the oven.



Note

Caution! The oven is hot. Use the crimping tongs to insert the sample vessel.



- 1** Grip the guide head with both hands and press the fixing lever.

Figure 23 Move the guide head down

- 2 Move the guide head down and insert the needle into the sample vessel.




Figure 24 Insert the needle

Switching on the gas flow

- 1 Press the **[OK]** key and adjust the flow rate (**Gas flow**), e. g. 50 mL/min. Use the arrow keys **[▶]** or **[◀]** to select the number of decimal places, and then **[▲]** or **[▼]** to change the value.
- 2 Confirm the entry by pressing **[OK]** and select **Flow source**.
- 3 Using the arrow keys **[▶]** or **[◀]**, select **pump** or **valve** (if you are using an external carrier gas) and confirm by pressing **[OK]**.
- 4 Close the Parameters menu by pressing **[QUIT]**.
- 5 Press the **[FLOW]** key.

The gas flow is switched on. The flow rate gradually approaches the level set.



Note

The carrier gas should emerge from the tip of the heating tube into the working medium as tiny bubbles and be spun there forcefully. Alter the stirring rate if necessary. The tip of the tubing should be immersed until it reaches the bottom of the KF titration cell and be rotated against the vessel wall.

Allow the system to condition for a few minutes and swirl the KF titration cell now and again in order to remove the last traces of moisture. Check whether the set oven temperature has been reached.



5.2 Determining the blank value



Note

Caution! The vessel is hot.

- 4 Place the blank into the oven using the crimping tongs.
- 5 Use both hands to move the guide head downward and insert the needle into the blank.

Waiting for the determination

- 1 Wait until the determination is finished. Conditioning will then be resumed automatically.
- 2 Repeat the blank value determination with the next blank.

5.3 Determining samples

The blank value must be subtracted from the calculated consumption when determining a sample.

Preparing the KF titrator

- 1 Load a suitable method.
- 2 Start the conditioning by pressing the **[Start]** key.

Preparing the sample

- 1** Weigh the sample into a sample vessel and seal it tightly with the crimping tongs.
- 2** Check that the oven temperature and the gas flow rate correspond to the values which have been set.
- 3** Wait until the working medium in the KF titration vessel has been conditioned and the signal drift is stable.

- 1 Start the determination by pressing the **[Start]** key on the KF titrator.
- 2 Use both hands to move the guide head of the 860 KF Thermoprep upward while pressing the fixing lever.
Refer to the illustrations on *pages 31 ff.*
- 3 Remove the blank from the oven with the crimping tongs.



Caution! The vessel is hot.

- 4 Place the sample vessel into the oven using the crimping tongs.
- 5 Use both hands to move the guide head downward and insert the needle into the sample vessel.

- 1 Wait until the determination is finished. Conditioning will then be resumed automatically.
- 2 Repeat the determination with the next blank.

6 Handling and maintenance

6.1 General

The 860 KF Thermoprep requires appropriate care. Excess contamination of the instrument may result in functional disruptions and a reduction in the service life of the sturdy mechanics and electronics of the instrument.

Severe contamination can also have an influence on the measured results. Regular cleaning of exposed parts can prevent this to a large extent.

Spilled chemicals and solvents must be removed immediately. In particular, the mains plug should be protected from contamination.

6.2 Care

- Replace exhausted molecular sieves promptly. You should replace the molecular sieve as soon as increased drift values appear in the Karl Fischer cell.
- Check all tubing connections regularly for leaks.
- Flush out the tubing connections from time to time. Carefully dry the tubing afterward. The tubing must be replaced after prolonged usage.
- Replace the dust filter once a year.
- Clean the sample insert regularly.



Warning

Clean the oven only when it is switched off and cold.

6.3 Quality Management and validation with Metrohm

Quality Management

Metrohm offers you comprehensive support in implementing quality management measures for instruments and software. You can find information on this in the brochure available from your local Metrohm agent **«Quality Management with Metrohm»**.

Validation

Please contact your local Metrohm agent for support in validating instruments and software. Here you can also obtain validation documentation to provide help for carrying out the **Installation Qualification** (IQ = Installation Qualification) and the **Operational** (OQ = Operational Qualification). IQ and OQ are also offered as a service by the Metrohm agents. In

Maintenance



7 Troubleshooting

7.1 Problems and their solutions

Problem	Cause	Remedy
The drift is very high during conditioning.	<i>Molecular sieve of drying flasks and/or the titration cell exhausted.</i>	Replace molecular sieve.
	<i>The titration cell is leaking</i>	Check seals. If necessary, replace.
The titration time is too long.	<i>The sample is non-homogenous</i>	Reduce the size of the sample before weighing in.
	<i>Not the entirety of the moisture expelled.</i>	Select more stringent switch-off criteria on the KF titrator: lower stop drift, higher switch-off delay time.
	<i>Tubing connections leaking</i>	Inspect tubing and replace as necessary.
The results are spread widely.	<i>Molecular sieve of drying flasks exhausted.</i>	Replace molecular sieve.
	<i>Condensate in the heating tube.</i>	<ul style="list-style-type: none"> ▪ Dry the tubing. ▪ Reduce the gas flow. ▪ Perhaps lower the oven temperature.
	<i>Gas flow too high.</i>	Reduce the gas flow.
	<i>The sample is non-homogenous</i>	Reduce the size of the sample before weighing in.
Selected gas flow is not achieved.	<i>System is leaking</i>	Check all tubing connections and drying flasks for leaks.
	<i>Pressure from external gas supply either too low or too high.</i>	Adjust gas pressure (max. 0.3 bar).
	<i>Pump defective</i>	Contact a service technician.

8 Appendix

8.1 Practical instructions

Carrier gas selection

If the hot sample is sensitive to air or oxygen (decomposition) and releases substances that disrupt the KF reaction, then nitrogen (N_2) should be used as the carrier gas.

Temperature settings

The temperature selected should be as high as the sample allows (high temperature = shorter analysis time). The sample may not however decompose. It must not give off any oxidizable substances, only water.

The displayed temperature refers to the temperature in the heating block and not to the sample temperature. Depending on the size of the vial used, the gas flow and the temperature setting, the effective temperature of the sample can deviate up to 10%.

Gas flow

If the needle is inserted into the sample vessel until it stops, then the gas flow moves through the sample and then through the outlet needle into the titration vessel where the moisture is titrated. If the needle is not permitted to dip into the sample, then a longer needle holder can be used, see the chapter "Optional accessories".

The gas flow should be kept as low as possible. Particularly in the case of very moist samples, care must be taken to ensure that excessively large quantities of water are not released too quickly into the titration vessel. In addition, there is the danger of condensate formation in the heating tube. The working medium in the cell must be able to absorb the moisture that is expelled without delay. Normally, a flow rate of 40...60 mL/min is optimal.

Extraction time

An extraction time of a minimum of 5 minutes should be set on the titrator in order to prevent the titration from being stopped before the sample has released its water.

Conditioning the system

The system must be conditioned with an empty, sealed sample vial (conditioning vessel) before a determination is made.



8.2 Literature

- E. Scholz, *Hydranal®-Manual*, Riedel-de Haën, 1996
- P. Bruttel, R. Schlink, *Water determination by Karl Fischer Titration*. Monograph, Metrohm, 2006, Best. Nr. 8.026.5013

9 Technical data

9.1 Oven

<i>Temperature range</i>	50...250 °C
<i>Accuracy</i>	±3 °C
<i>Correction range</i>	-10...+10 °C
<i>Heating cartridge performance</i>	165 W typical Depending on mains voltage
<i>Heating rate</i>	Typically 15 °C/min (with 80...180 °C, 230 V) Dependent on temperature, mains voltage, sample amount and vessel dimensions
<i>Cooling rate</i>	Typical. 9 °C/min (at 180...80 °C) Dependent on temperature, sample amount and vessel dimensions

9.2 Gas flow

<i>Flow range</i>	10...150 mL/min Under normal conditions
-------------------	--

9.3 Outlet heating

<i>Socket connection</i>	U = 16 ±1 V I ≤ 0.8 A
<i>Typical tubing temperature</i>	approx. 50 °C

9.7 Ambient temperature

<i>Nominal function range</i>	5...45 °C Humidity < 80 %
<i>Storage</i>	−20...70 °C
<i>Transport</i>	−40...70 °C

9.8 Reference conditions

<i>Ambient temperature</i>	25 °C (±3 °C)
<i>Relative humidity</i>	≤60 %
<i>Mains voltage</i>	230 V

9.9 Dimensions

<i>Width</i>	0.44 m
<i>Height</i>	0.18 m (without support rod) 0.50 m (with support rod)
<i>Depth</i>	0.23 m
<i>Weight (without accessories)</i>	8.17 kg
<i>Material</i>	
<i>Housing</i>	Cover and base: steel sheet, stove-enameled
<i>Oven covering</i>	PTFE
<i>Operating unit</i>	Polycarbonate/Acrylonitrile butadiene styrene (PC/ABS)

EN 61010-2-010 Particular requirements for laboratory equipment
for the heating of materials

Manufacturer

Metrohm Ltd., CH-9101 Herisau/Switzerland

Metrohm Ltd. is holder of the SQS-certificate ISO 9001:2000 Quality management system for development, production and sales of instruments and accessories for ion analysis.

Herisau, 27 November, 2007



D. Strohm

Vice President, Head of R&D



A. Dellenbach

Head of Quality Management

10.2 Warranty (guarantee)

Metrohm guarantees that the deliveries and services it provides are free from material, design or manufacturing errors. The warranty period is 36 months from the day of delivery; for day and night operation it is 18 months. The warranty remains valid on condition that the service is provided by an authorized Metrohm service organization.

Glass breakage is excluded from the warranty for electrodes and other glassware. The warranty for the accuracy corresponds to the technical specifications given in this manual. For components from third parties that make up a considerable part of our instrument, the manufacturer's warranty provisions apply. Warranty claims cannot be pursued if the Customer has not complied with the obligations to make payment on time.

During the warranty period Metrohm undertakes, at its own choice, to either repair at its own premises, free of charge, any instruments that can be shown to be faulty or to replace them. Transport costs are to the Customer's account.

Faults arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc. are expressly excluded from the warranty.

tional Qualification and Performance Verification of the system components or in carrying out the System Validation for the quantitative determination of a substance in a given matrix.

11 Accessories






Note




Subject to change without notice.




11.1 Scope of delivery 2.860.0010

Qty.	Order no.	Description
1	1.860.0010	860 KF Thermoprep Instrument for the thermal extraction of humidity from solid and liquid samples.
1	6.1446.170	Heating tube stopper for coulometric KF vessels in connection with heatable tubing connection Material: PP Height (mm): 34 Outer diameter (mm): 13
1	6.1448.057	Aluminum septum caps / 100 pieces Including Al sealing. Material: Aluminum Material 2: Silicone Height (mm): 7.4 Outer diameter (mm): 20.7
2	6.1602.145	Cap for drying bottle Cap with GL45 thread for use with 6.1608.050





Qty.	Order no.	Description
2	6.1608.050	Drying bottle / 100 mL / GL 45
	Material:	Clear glass
	Height (mm):	100
	Outer diameter (mm):	56
	Volume (mL):	100
		
1	6.1805.010	FEP tubing / M6 / 13 cm
	With light and kink protection	
	Material:	FEP
	Inner diameter (mm):	2
	Length (mm):	130
		
2	6.1805.080	FEP tubing / M6 / 25 cm
	With light and kink protection	
	Material:	FEP
	Inner diameter (mm):	2
	Length (mm):	250
		



Qty.	Order no.	Description
1	6.1805.100	FEP tubing / M6 / 40 cm With light and kink protection Material: FEP Inner diameter (mm): 2 Length (mm): 400
		
1	6.1808.040	Adapter M6 outer / M8 inner Outer thread M6, inner thread M8. Material: PTCFE
		
1	6.1808.050	Connector tubing nozzle M8 1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of Exchange Units and stability measuring instruments. Material: PVDF Length (mm): 31.5
		

Qty.	Order no.	Description
2	6.1821.040	Filter tube Filter tube for 6.1608.050 Drying bottle. For Rancimats and Karl Fischer Ovens. Length (mm): 112
		
1	6.1830.030	Heatable outlet tubing Heatable outlet tubing for 860 KF Thermoprep.
		
1	6.2048.030	Guidance rod Guidance rod for 860 KF Thermoprep.
1	6.2049.040	Needle holder with Luer-lock for 860 KF Thermoprep Length (mm): 58
		






Qty.	Order no.	Description	
1	6.2060.010	Cable support Used with 860 KF Thermoprep.	
1	6.2063.010	Sample holder for 6.2419.007 sample vials (6 mL)	
1	6.2419.007	Sample vial 6 ml / 100 pieces 6.1448.050 septum caps. Used for KF-determinations with oven. Volume (mL): 6	
1	6.2621.110	Crimping tongs Crimping tongs for sealing 6 mL sample vials (6.2419.000, 6.2419.007). Length (mm): 210	

860 KF Thermoprep

Qty.	Order no.	Description	
1	6.2739.000	Wrench For tightening connectors Length (mm): 68	
1	6.2811.000	Molecular sieve Molecular sieve. Bottle containing 250 g. Pore size: 0.3 nm. Without moisture indicator. For Rancimats and Karl Fischer instruments.	
1	6.2816.070	Injection needle Used with 832 KF Thermoprep. Material: Stainless steel (AISI 304)	
1	6.2816.080	Outlet needle Used with KF Thermoprep and Oven Sample Processor. Material: Stainless steel (AISI 304)	
1	6.2122.0x0	Mains cable with C13 line socket IEC-60320-C13	

Qty.	Order no.	Description
		Cable plug according to customer requirements.
		Switzerland: Type SEV 12 6.2122.020
		Germany, ...: Type CEE(7), VII 6.2122.040
		USA, ...: Type NEMA/ASA 6.2122.070
1	8.860.8003EN	860 KF Thermoprep Manual

11.2 Optional accessories

Order no.	Description	
6.1448.050	Aluminum septum caps / 1000 pieces	
	Including Al sealing.	
	Material: Aluminum	
	Material 2: Silicone	
	Height (mm): 7.4	
	Outer diameter (mm): 20.7	
		
6.2049.050	Needle holder with Luer-lock for 860 KF Thermoprep	
	Length (mm): 73	
		
6.2419.000	Sample vial 6 mL / 1000 pieces	
	For 6.1448.050 septum caps.	
	Height (mm): 38	
	Outer diameter (mm): 22	
	Volume (mL): 6	
		

Numbers/Symbols

A

Assembling

B

C

D

E

F

G

H

1

Inlet

Inlet filter 16

K

L

Lock keys 28

M

N

O

P

Q

S

T

V

W