

862 Compact Titrosampler



Manual
8.862.8002EN



Metrohm AG
CH-9101 Herisau
Switzerland
Phone +41 71 353 85 85
Fax +41 71 353 89 01
info@metrohm.com
www.metrohm.com

862 Compact Titrosampler

Manual

Teachware
Metrohm AG
CH-9101 Herisau
teachware@metrohm.com

This documentation is protected by copyright. All rights reserved.

Although all the information given in this documentation has been checked with great care, errors cannot be entirely excluded. Should you notice any mistakes please send us your comments using the address given above.

Table of contents

1	Introduction	1
1.1	Instrument description	1
1.1.1	Instrument components	1
1.1.2	Titration and measuring modes	2
1.1.3	Intended use	2
1.2	About the documentation	3
1.2.1	Symbols and conventions	3
1.3	Safety instructions	4
1.3.1	General notes on safety	4
1.3.2	Electrical safety	4
1.3.3	Tubing and capillary connections	5
1.3.4	Personnel safety	5
1.3.5	Flammable solvents and chemicals	7
1.3.6	Recycling and disposal	7
2	Overview of the instrument	8
3	Installation	10
3.1	Setting up the instrument	10
3.1.1	Packaging	10
3.1.2	Checks	10
3.1.3	Location	10
3.2	Removing the safety shield and cable cover	11
3.3	Mounting the Dosino	11
3.4	Connecting the Dosino	13
3.5	Setting up the titration head	14
3.6	Connecting a stirrer	17
3.7	Connecting keyboard, printer and other USB devices	18
3.8	Connecting a balance	20
3.9	Connecting a sensor	21
3.10	Remote connections	23
3.10.1	Miscellaneous remote cables	23
3.10.2	Example systems	24
3.11	Mounting the cable cover and the safety shield	26
3.12	Connecting the mains cable	27
4	Titration and automation sequences	28
4.1	Dynamic equivalence point titration (DET)	28



4.2	Monotonic equivalence point titration (MET)	28
4.3	Endpoint titration (SET)	29
4.4	Automation sequences	29
4.4.1	Dipping in special	29
4.4.2	Dipping in special 2	30
4.4.3	Double dipping	31
4.4.4	Rinsing in sample	32
4.4.5	Rinsing in special	32
4.4.6	Pump control	33
4.4.7	Dosing auxiliary solutions	34
5	Operation	35
5.1	Switching the instrument on and off	35
5.2	Fundamentals of operation	36
5.2.1	The keypad	36
5.2.2	Structure of the dialog windows	37
5.2.3	Navigating in the dialog	37
5.2.4	Entering text and numbers	38
5.2.5	Selecting from a selection list	39
5.3	Formula editor	39
5.4	Methods	41
5.4.1	Creating a new method	41
5.4.2	Saving a method	42
5.4.3	Loading a method	43
5.4.4	Exporting a method	44
5.5	Sample data	45
5.5.1	Sample table	45
5.5.2	Entering sample data in the main dialog	48
5.6	Performing a sample series	49
5.6.1	Starting the sample series	49
5.6.2	Pausing a sample series and continuing	50
5.6.3	Status of the sample series	52
5.6.4	Special case: Calibrations	52
5.7	Live modifications	53
5.7.1	Editing the sample data of the running determination	53
5.7.2	Editing the sample table while a determination is running	54
5.7.3	Live parameters	55
5.8	Results	56
5.9	Statistics	58
5.10	Printing a report manually	59
5.11	Manual control	61
5.11.1	Rotating the sample rack	61
5.11.2	Moving the lift	62
5.11.3	Dosing	62

5.11.4	Measuring	66
5.11.5	Stirring	67
6	System settings	69
6.1	Basic settings	69
6.2	Lift settings (Lift)	72
6.3	Managing sensors	73
6.3.1	General	73
6.3.2	Editing the sensor data	74
6.4	Managing solutions	75
6.4.1	General	75
6.4.2	Editing the solution data	76
6.5	Managing common variables	77
6.5.1	General	77
6.5.2	Editing common variables	78
6.6	File management	79
6.7	Configuring external devices	80
6.8	Instrument diagnosis	83
6.8.1	Loading program versions and language files	83
6.8.2	Diagnosis functions	84
7	Parameters	85
7.1	Dynamic equivalence point titrations (DET)	85
7.1.1	Start conditions	85
7.1.2	Titration parameters	86
7.1.3	Stop conditions	91
7.1.4	Evaluation	93
7.1.5	Calculation	96
7.1.6	Statistics	98
7.1.7	Reports	98
7.2	Monotonic equivalence point titrations (MET)	100
7.2.1	Start conditions	100
7.2.2	Titration parameters	101
7.2.3	Stop conditions	106
7.2.4	Evaluation	107
7.2.5	Calculation	110
7.2.6	Statistics	112
7.2.7	Reports	112
7.3	Endpoint titrations (SET)	114
7.3.1	Conditioning	114
7.3.2	Start conditions	115
7.3.3	Titration parameters	117
7.3.4	Control parameters EP1	119
7.3.5	Control parameters EP2	122
7.3.6	Stop conditions	122



7.3.7	Calculation	123
7.3.8	Statistics	125
7.3.9	Reports	125
7.4	pH calibration (CAL)	127
7.4.1	Calibration parameters	127
7.4.2	Buffers	128
7.4.3	Reports	129
7.5	Automation: Dipping in special	130
7.6	Automation: Dipping in special 2	131
7.7	Automation: Double dipping	132
7.8	Automation: Rinsing in sample	133
7.9	Automation: Rinsing in special	134
8	Handling and maintenance	135
8.1	General	135
8.2	Quality Management and validation with Metrohm	135
9	Troubleshooting	137
9.1	SET titration	137
9.2	Miscellaneous	138
10	Appendix	139
10.1	Dosing unit	139
10.1.1	Maximum dosing and filling rate	139
10.1.2	Parameters for preparing (PREP) and emptying (EMPTY)	139
10.2	Stirring rate	140
10.3	Balance	140
10.4	USB devices	141
10.4.1	Numerical USB keypad 6.2147.000	141
10.4.2	Key assignment of a USB keyboard	142
10.4.3	PC mouse	142
10.4.4	Printer	143
10.5	System initialization	143
10.6	Remote interface	145
10.6.1	Pin assignment of the remote interface	145
10.6.2	Status diagram of the remote interface	146
11	Technical specifications	147
11.1	Measuring inputs	147
11.1.1	Potentiometry	147
11.1.2	Polarizer	147
11.1.3	Temperature	148

11.2	Dosing drive	149
11.3	Lift	149
11.4	Turntable	149
11.5	Interfaces and connectors	149
11.6	Mains connection	150
11.7	Safety specifications	150
11.8	Electromagnetic compatibility (EMC)	150
11.9	Ambient temperature	151
11.10	Reference conditions	151
11.11	Dimensions	151
12	Conformity and warranty	152
12.1	Declaration of Conformity	152
12.2	Warranty (guarantee)	153
12.3	Quality Management Principles	154
13	Accessories	156
13.1	Scope of delivery 2.862.0010	156
13.2	Optional accessories	160
	Index	163



Table of figures

Figure 1	Front 862 Compact Titrosampler	8
Figure 2	Rear 862 Compact Titrosampler	9
Figure 3	Removing shields	11
Figure 4	Mounting the Dosino	12
Figure 5	Connecting the Dosino	13
Figure 6	Equip the titration head	14
Figure 7	Installing rinsing nozzles and aspiration tip	15
Figure 8	Setting up the distributor	16
Figure 9	Rod stirrer 802 Stirrer	17
Figure 10	Connecting a stirrer	17
Figure 11	Connecting USB devices	18
Figure 12	Connecting the USB stick	19
Figure 13	Connecting the 6.2147.000 USB keyboard with USB stick and printer	20
Figure 14	Connecting the USB hub with USB stick, printer and the 6.2148.030 RS-232/USB Box (for connecting balances)	20
Figure 15	Connecting a balance	20
Figure 16	Connecting pH or redox electrodes	21
Figure 17	Connecting a reference electrode	22
Figure 18	Connecting a temperature sensor	22
Figure 19	Connecting a polarizable electrode	23
Figure 20	Remote cable with lettering	24
Figure 21	Remote connection 862 Compact Titrosampler - 843 Pump Station - Dosimat plus	24
Figure 22	Remote connection 862 Compact Titrosampler - Dosimat	25
Figure 23	Mounting the shields	26
Figure 24	Connecting the mains cable	27
Figure 25	Reagent dosing for DET	28
Figure 26	Reagent dosing for MET	28
Figure 27	Reagent dosing for SET	29
Figure 28	Keypad 862 Compact Titrosampler	36
Figure 29	Directory structure on the USB stick	80
Figure 30	Tubbs method for determining the equivalence point	96
Figure 31	Rotational speed depending on stirring rate	140
Figure 32	Pin assignment of remote socket and plug	145
Figure 33	Remote status diagram	146

1 Introduction

1.1 Instrument description

The 862 Compact Titrosampler is a versatile analysis instrument which combines a titrator and a compact sample changer in a single unit. It is the central control instrument in an automation system that can also include, in addition to a Dosimat (for adding auxiliary solutions), a pump for rinsing and aspirating the sample vessels.

Given automation sequences and titration modes can be individually parameterized and saved as sample-specific methods. Methods can be exported on a USB stick and then copied onto another instrument quickly and easily.

1.1.1 Instrument components

The 862 Compact Titrosampler has the following components:

- **Turntable**
Permanently mounted sample rack with 11 positions for sample beakers and one position for a rinsing beaker.
- **Lift with titration head**
For two electrodes, one rod stirrer, two dosing tips, one aspiration tip and three rinsing nozzles.
- **Sensor connectors**
Four connectors for the following sensor types:
 - pH or redox electrodes
 - reference electrodes
 - polarizable electrodes
 - temperature sensors (Pt1000 or NTC)
- **Stirrer connector**
For a rod stirrer with propeller stirrer.
- **MSB connector (Metrohm Serial Bus)**
For connecting a Dosino.
- **USB (OTG) connector**
The 6.2151.100 adapter can be used to connect, for example, a USB hub or a USB stick (for system backup or method export).
- **Remote connector**
For connecting a Dosimat and/or an 843 Pump Station as well as other instruments with a remote interface.



1.1.2 Titration and measuring modes

The following titration and measuring modes are supported:

- **DET**
Dynamic equivalence point titration. The reagent is added in variable volume steps.
Measuring modes:
 - **pH** (potentiometric pH measurement)
 - **U** (potentiometric voltage measurement)
 - **Ip_{ol}** (voltametric measurement with adjustable polarization current)
 - **Up_{ol}** (amperometric measurement with adjustable polarization voltage)
- **MET**
Monotonic equivalence point titration. The reagent is added in constant volume steps.
Measuring modes:
 - **pH** (potentiometric pH measurement)
 - **U** (potentiometric voltage measurement)
 - **Ip_{ol}** (voltametric measurement with adjustable polarization current)
 - **Up_{ol}** (amperometric measurement with adjustable polarization voltage)
- **SET**
Endpoint titration at one or two specified endpoints.
Measuring modes:
 - **pH** (potentiometric pH measurement)
 - **U** (potentiometric voltage measurement)
 - **Ip_{ol}** (voltametric measurement with adjustable polarization current)
 - **Up_{ol}** (amperometric measurement with adjustable polarization voltage)
- **CAL**
Electrode calibration.
Measuring mode:
 - **pH** (calibration of pH electrodes)

1.1.3 Intended use

The 862 Compact Titrosampler is designed for usage as an automation system in analytical laboratories. It is **not** suitable for usage in biochemical, biological or medical environments in its basic equipment version.

This instrument is suitable for dosing chemicals and flammable solvents. The usage of the 862 Compact Titrosampler therefore requires that the user has basic knowledge and experience in the handling of toxic and caustic substances. Knowledge with respect to the application of the fire

prevention measures prescribed for laboratories or production plants is also mandatory.

1.2 About the documentation







Caution

Please read through this documentation carefully before putting the instrument into operation. The documentation contains information and warnings which have to be followed by the user in order to ensure safe operation of the instrument.

1.2.1 Symbols and conventions

The following symbols and styles are used in this documentation:

(5-12)	Cross-reference to figure legend The first number refers to the figure number, the second to the instrument part in the figure.
1	Instruction step Carry out these steps in the sequence shown.
Method	Dialog text, parameter in the software
File ▶ New	Menu or menu item
[Next]	Button or key
	Warning This symbol draws attention to a possible life hazard or risk of injury.
	Warning This symbol draws attention to a possible hazard due to electrical current.
	Warning This symbol draws attention to a possible hazard due to heat or hot instrument parts.
	Warning This symbol draws attention to a possible biological hazard.



	<p>Caution</p> <p>This symbol draws attention to a possible damage of instruments or instrument parts.</p>
	<p>Note</p> <p>This symbol marks additional information and tips.</p>

1.3 Safety instructions

1.3.1 General notes on safety



Warning

This instrument may only be operated in accordance with the specifications in this documentation.

This instrument has left the factory in a flawless state in terms of technical safety. To maintain this state and ensure non-hazardous operation of the instrument, the following instructions must be observed carefully.

1.3.2 Electrical safety

The electrical safety when working with the instrument is ensured as part of the international standard IEC 61010.



Warning

Only personnel qualified by Metrohm are authorized to carry out service work on electronic components.



Warning

Never open the housing of the instrument. The instrument could be damaged by this. There is also a risk of serious injury if live components are touched.

There are no parts inside the housing which can be serviced or replaced by the user.

Mains voltage



Warning

An incorrect mains voltage can damage the instrument.

Only operate this instrument with a mains voltage specified for it (see rear panel of the instrument).

Protection against electrostatic charges



Warning

Electronic components are sensitive to electrostatic charges and can be destroyed by discharges.

Always pull the mains cable out of the mains connection socket before connecting or disconnecting electrical appliances on the rear panel of the instrument.

1.3.3 Tubing and capillary connections



Caution

Leaks in tubing and capillary connections are a safety risk. Tighten all connections well by hand. Avoid applying excessive force to tubing connections. Damaged tubing ends lead to leakage. Appropriate tools can be used to loosen connections.

Check the connections regularly for leakage. If the instrument is used mainly in unattended operation, then weekly inspections are mandatory.

1.3.4 Personnel safety



Warning

Wear protective goggles and working clothes suitable for laboratory work while operating the 862 Compact Titrosampler. It is also advisable to wear gloves when caustic liquids are used or in situations where glass vessels could break.



Warning

Always install the safety shield supplied with the equipment before using the instrument for the first time. Pre-installed safety shields are not allowed to be removed.

The 862 Compact Titrosampler may not be operated without a safety shield!



Warning

Personnel are not permitted to reach into the working area of the instrument while operations are running!

A **considerable risk of injury** exists for the user.



Warning

In the event of a possible blockage of a drive, the mains plug must be pulled out of the socket immediately. Do not attempt to free jammed sample vessels or other parts while the device is switched on. Blockages can only be cleared when the instrument is in a voltage-free status; this action generally involves a **considerable risk of injury**.



Warning

The 862 Compact Titrosampler is **not** suitable for utilization in biochemical, biological or medical environments in its basic equipment version.

Appropriate protective measures must be implemented in the event that potentially infectious samples or reagents are being processed.

1.3.5 Flammable solvents and chemicals




Warning

All relevant safety measures are to be observed when working with flammable solvents and chemicals.

- Set up the instrument in a well-ventilated location.
- Keep all sources of flame far from the workplace.
- Clean up spilled fluids and solids immediately.
- Follow the safety instructions of the chemical manufacturer.

1.3.6 Recycling and disposal



This product is covered by European Directive 2002/96/EC, WEEE – Waste from Electrical and Electronic Equipment.

The correct disposal of your old equipment will help to prevent negative effects on the environment and public health.

More details about the disposal of your old equipment can be obtained from your local authorities, from waste disposal companies or from your local dealer.

2 Overview of the instrument




Figure 1 Front 862 Compact Titrosampler

1 Safety shield (6.2751.130)

With knurled screws for mounting. The safety shield can be folded up.

3 Sample rack

For 11 sample beakers and one rinsing beaker (6.1459.300, 120 mL)

5 Keypad

2 Titration head

With two retracted dosing tips and four stoppers.

4 Display

6 Dosino holder (6.2057.110)

For fixing the dosing drive and the titrant bottle (1 L).

7 Special position

For a rinsing beaker




Figure 2 Rear 862 Compact Titrosampler

<p>1 Warning symbol (see Chapter 1.3.4, page 5)</p>	<p>2 Tubing and cable cover</p>
<p>3 Electrode connector (Ind.) For connecting pH or redox electrodes with integrated or separate reference electrode. Socket F.</p>	<p>4 Electrode connector (Ref.) For connecting reference electrodes. Socket B, 4 mm.</p>
<p>5 Temperature sensor connector For connecting temperature sensors of the Pt1000 or NTC types. Two B sockets, 2 mm.</p>	<p>6 Electrode connector (Pol.) For connecting polarizable electrodes, e.g. double Pt electrodes. Socket F.</p>
<p>7 Remote connector For connecting instruments with a remote interface. D-Sub, 9-pin.</p>	<p>8 MSB connector Metrohm Serial Bus. For connecting the 800 Dosino. Mini DIN, 9-pin.</p>
<p>9 USB (OTG) connector For connecting printers, USB sticks, USB hubs, etc.</p>	<p>10 Stirrer connection For 802 Stirrer (rod stirrer).</p>
<p>11 Mains connection socket</p>	<p>12 Type plate Contains specifications concerning mains voltage and serial number.</p>



3 Installation

3.1 Setting up the instrument

3.1.1 Packaging

The instrument is supplied in highly protective special packaging together with the separately packed accessories. Keep this packaging, as only this ensures safe transportation of the instrument.

3.1.2 Checks

Immediately after receipt, check whether the shipment has arrived complete and without damage by comparing it with the delivery note.

3.1.3 Location

The instrument has been developed for operation indoors and may not be used in explosive environments.

Place the instrument in a location of the laboratory which is suitable for operation, free of vibrations, protected from corrosive atmosphere, and contamination by chemicals.

The instrument should be protected against excessive temperature fluctuations and direct sunlight.

3.2 Removing the safety shield and cable cover




Figure 3 Removing shields

If you remove the safety shield and the cable cover before, the installation of the accessories is easier to carry out. Proceed as follows:

- 1** Loosen the knurled screws on the sides of the tower and remove the safety shield.
- 2** Loosen the knurled screws on the rear of the tower and remove the cable cover.

Do not forget to refasten these two shields after the installation of the accessories.

3.3 Mounting the Dosino

The 800 Dosino is used for adding titrant. The titrant bottle can be placed next to the tower of the 862 Compact Titrosampler.




Figure 4 Mounting the Dosino

1 Attach the dosing unit

Connect the filling tubing to port 2 on the underside of the dosing unit and screw the dosing unit onto the titrant bottle.

2 Place the titrant bottle

Open the bracket on the right-hand side of the tower with the aid of the knurled screw. Place the titrant bottle with the dosing unit on the support surface. Fix the bottle with the bracket and fasten the knurled screw.

3 Attach the dosing drive

Insert the 800 Dosino with the guide pins into the openings on the upper side of the dosing unit. Fix the 800 Dosino with a rotation to the left. Note the markings.

3.4 Connecting the Dosino

Connect the Dosino to the **MSB** socket on the rear of the instrument. The Dosino is automatically recognized when switching on the instrument.




Figure 5 Connecting the Dosino



Caution

Make sure that the flat side of the plug matches the marking on the socket.



3.5 Setting up the titration head

Equipping the titration head without rinsing equipment




Figure 6 Equip the titration head

- 1 Insert the rod stirrer (802 Stirrer) from above into the rear opening of the titration head.
- 2 Guide the 6.1909.050 stirring propeller over the drive shaft of the rod stirrer from below and press firmly.
- 3 Insert an electrode into the left-hand opening of the titration head.
- 4 Connect a 6.2104.020 electrode cable to the electrode. Connect the other end to the electrode connector **Ind.** (see "Connecting pH or redox electrodes", page 21).
- 5 Manually screw the enclosed 6.1805.100 FEP tubing to the dosing tip mounted on the titration head. Connect the other end of the tubing to the dosing unit on the Dosino.

The remaining openings of the titration head can be sealed with the enclosed stoppers provided.

Equipping the titration head with rinsing equipment

An 843 Pump Station (with two pumps) can be used when the sample processing requires the rinsing of the electrodes and the aspiration of the processed sample solution. The 843 Pump Station is available as model version with complete rinsing and aspiration equipment. The rinsing and aspiration equipment is installed as follows:




Figure 7 Installing rinsing nozzles and aspiration tip

- 1** Insert the three 6.2740.020 spray nozzles into the titration head according to the illustration. The position of the individual spray nozzles can also be adjusted as required in terms of height.
- 2** Manually screw the three 6.1805.420 FEP tubings (with M6 thread) firmly to the spray nozzles.
- 3** Insert the 6.1543.120 aspiration tip into the front opening of the titration head. It can be adjusted in terms of height and its tip can be cut as required to the necessary length.
- 4** Manually screw the 6.1805.200 aspiration tubing (with M8 thread) firmly to the aspiration tip.



Setting up the distributor

The 6.1818.240 distributor must be mounted on the rear side of the tower for complete installation of the rinsing and aspiration equipment. It is supplied with the 843 Pump Station.

First loosen the knurled screws of the cable and tubing cover and then remove it. Proceed as follows:




Figure 8 Setting up the distributor

- 1
 - Loosen both union nuts on the distributor and guide each of them over the end of one 6.1812.000 PTFE tubing.
 - Fasten the tubing ends to the distributor and fix in place with the union nuts.
 - Connect the free tubing ends with a disposal or rinsing canister.
- 2 Manually screw the 6.1805.420 rinsing tubings already mounted on the titration head firmly into the openings with M6 threaded bores on the distributor.
- 3 Manually screw the 6.1805.200 aspiration tubing with M8 thread firmly into the remaining opening on the distributor.
- 4 Loosen the two screws on the rear panel of the instrument with a hexagon key and use it to screw the distributor firmly.



Note

Enclosed with the 862 Compact Titrosampler is the 6.1815.010 spiral band. You can wrap cables and tubings with it. This will ensure that the cables and tubings are arranged in an organized manner.

- 5 Use the four knurled screws to remount the cable and tubing cover.



Caution

Close the safety shield again after the titration head has been equipped. The 862 Compact Titrosampler is not permitted to be operated unless the safety shield is correctly mounted.

3.6 Connecting a stirrer

A DIN socket for connecting a rod stirrer **802 Stirrer** is located on the rear of the instrument.




Figure 9 Rod stirrer 802 Stirrer

Take care to observe correct orientation of the contact pins when plugging in the connection cable. The rib on the outside of the plug must match the reference mark (above) on the socket.




Figure 10 Connecting a stirrer



3.7 Connecting keyboard, printer and other USB devices

The 862 Compact Titrosampler has a USB (OTG) connector. Use the provided 6.2151.100 adapter USB MINI (OTG) - USB A for connecting USB devices as e.g. printers, keyboards or USB sticks, see the following figure.




Figure 11 Connecting USB devices



Caution

Switch the instrument off before connecting or disconnecting a USB device or a USB stick.

The 862 Compact Titrosampler can recognize the device immediately after switching on.

The following devices can directly be operated on the **USB connector with the 6.2151.100 adapter**:

- USB sticks (for the backup or storing of methods)
- 6.2147.000 numerical USB keypad
- USB hub (with or without an own power supply)

The **6.2147.000 numerical USB keypad** serves for comfortable numerical input and for navigating in the dialog. In addition, it provides two USB connectors. Connect additional USB devices to the keypad.



Note

Most of the USB devices need a so-called hub in order to work correctly.

A USB hub is a distributor to which several USB devices can be connected. USB hubs are available in specialty stores in a number of different models.

The USB (OTG) connector of the 862 Compact Titrosampler has no such hub. The 6.2147.000 numerical USB keypad has a USB hub and two USB connectors.

The following devices can **only be connected to a 6.2147.000 numerical keypad or to a USB hub**:

- Printer (with USB connector, use the 6.2151.020 connecting cable)
- Barcode reader (with USB cable)
- Mouse (PC mouse with USB cable, for navigating in the dialog)

The following devices can **only be connected to a USB hub**:

- PC keyboard (with USB cable, for the comfortable input of letters and numbers)
- Keypad with numerical keypad (with USB cable)

If you wish to connect a **variety of instruments without own power supply**, then you must possibly use a USB hub with own power supply (*self powered*). The USB (OTG) connector of the 862 Compact Titrosampler is not designed for supplying power to several devices with elevated electricity requirements.

Examples:




Figure 12 Connecting the USB stick




Figure 13 Connecting the 6.2147.000 USB keyboard with USB stick and printer




Figure 14 Connecting the USB hub with USB stick, printer and the 6.2148.030 RS-232/USB Box (for connecting balances).

3.8 Connecting a balance

Balances are equipped with a serial RS-232 interface as a rule. To connect a balance, you require a 6.2148.030 RS-232/USB Box.




Figure 15 Connecting a balance

When a 6.2151.020 USB cable is used, then the 6.2148.030 RS-232/USB Box can be connected to the 862 Compact Titrosampler by means of a USB hub or a 6.2151.100 adapter (see *Chapter 3.7, page 18*).

Connect the 9-pin plug of the respective balance connecting cable to the **RS 232/1** connector. Consult the user manual of the balance in order to select the correct connecting cable.

The parameters for the RS-232 interface on the instrument must match those on the balance (see *"Editing the COM1 settings", page 82*). Additionally consult the user manual of the balance.

3.9 Connecting a sensor

Connecting pH or redox electrodes




Figure 16 Connecting pH or redox electrodes



Note

The electrode cable is protected against accidental disconnection of the cable by means of a pull-out protection feature. If you wish to remove the plug, then you must first retract the outer plug sleeve.



Connecting a reference electrode




Figure 17 Connecting a reference electrode

Connecting a temperature sensor or an electrode with integrated temperature sensor

Temperature sensors of the Pt1000 or NTC type can be connected to the **Temp.** connector.




Figure 18 Connecting a temperature sensor



Note

The red plug must always be plugged into the red socket at the temperature sensor for the purpose of shielding against disruptions.

If you use an electrode with an integrated NTC probe, then you must plug the red plug into the red socket.

Connecting a polarizable electrode




Figure 19 Connecting a polarizable electrode



Note

The electrode cable is protected against accidental disconnection of the cable by means of a pull-out protection feature. If you wish to remove the plug, then you must first retract the outer plug sleeve.

3.10 Remote connections

The 862 Compact Titrosampler can be used as a control instrument for a simple automation system with a large variety of different instruments. Even older Metrohm instruments can thus be integrated into an automated analysis system.

3.10.1 Miscellaneous remote cables

The following connecting cables can be used with the 862 Compact Titrosampler:

6.2136.010

- For connections with a Dosimat with dosing contact (banana plug socket).

The cable transmits only a starting pulse from the 862 Compact Titrosampler to the connected Dosimat.

6.2141.230

- For connections to an 843 Pump Station.

The cable transmits the control signals of the 862 Compact Titrosampler to the pump 1 and 2 of the 843 Pump Station.

6.2141.240

- For connections with a Dosimat plus.

The cable transmits start and stop pulses from the 862 Compact Titrosampler to the connected Dosimat plus.

In the event of an error with the connected Dosimat, the cable transmits a stop signal to the 862 Compact Titrosampler.



Caution

Some of these cables have asymmetrical wiring. You must connect the correct plug to a particular instrument in each case. Observe the lettering on the ends of the cables. The designation of the instrument is explicitly printed on the end of the cable wherever necessary.




Figure 20 Remote cable with lettering

3.10.2 Example systems

The following illustrations show typical automation systems with different instrument combinations.

862 — 843 Pump Station — Dosimat plus

The standard combination for titrations, with adding auxiliary solution by a Dosimat plus, with rinsing and aspirating.




Figure 21 Remote connection 862 Compact Titrosampler - 843 Pump Station - Dosimat plus

The Dosimat is operated in XDOS mode. The volume of the auxiliary solution is defined on the Dosimat plus. On the 843 Pump Station the 862 is connected to **Remote 1**, the Dosimat plus to **Remote 2**. Pump 1 is used for rinsing the electrode, pump 2 for aspirating the sample solution. The sample series is started on the 862 Compact Titrosampler



3.11 Mounting the cable cover and the safety shield

After having installed all accessories you can remount the shields. Proceed as follows:




Figure 23 Mounting the shields

1 Mount the cable cover

Fasten the cable cover with the four red knurled screws to the rear of the tower. Ensure that all cables and tubings are routed in order.

2 Mount the safety shield

Fasten the safety shield with the four black knurled screws to the sides of the tower. When each of the lower knurled screws is slightly loosened, the safety shield can be folded up if needed.



Warning

The 862 Compact Titrosampler may not be operated without a safety shield!

3.12 Connecting the mains cable



Warning

This instrument must not be operated except with the mains voltage specified for it (see rear of the instrument).

Protect the connection sockets against moisture.

The connector for the mains cable is located on the rear of the instrument.




Figure 24 Connecting the mains cable

4.3 Endpoint titration (SET)

Endpoint titration is a titration mode for rapid routine determinations to a preset endpoint (e.g. titrations in accordance with special norms) and titrations for which reagent overflow must be avoided. The titration termination at the endpoint takes place either drift-controlled or after a waiting period. The volume dosed until the endpoint is used for calculating the content of the sample.




Figure 27 Reagent dosing for SET

4.4 Automation sequences

4.4.1 Dipping in special

This automation sequence is suitable for simple determinations.

No pump is required for rinsing or for aspiration of the sample vessels. The titration head with electrode and burette tips is immersed in the filled rinsing beaker on the special beaker position after each determination. At the same time, the rinsing solution is stirred.

If required, auxiliary solution can be added prior to the beginning of a determination with a Dosimat/Dosimat plus.

A filled rinsing beaker should be placed on the **special beaker position**.

The individual steps:

- Move to the sample
- Lower the lift to work position
- Start the determination



- If necessary, initiate dosing (**Activation pulse**) and switch on the stirrer
- If necessary, wait for **Start delay time**
- Wait for end of determination
- Switch off the stirrer and move the lift upward
- Wait for **Dripping time** to be completed
- Moving to special beaker position
- Lower the lift to work position and switch on the stirrer
- Wait for **Rinsing time** to be completed
- Switch off the stirrer and move the lift upward
- Wait for **Dripping time** to be completed

After the last sample:

- Lower the lift in the rinsing beaker to work position

4.4.2 Dipping in special 2

This automation sequence is suitable for simple determinations.

No pump is required for rinsing or for aspiration of the sample vessels. The titration head with electrode and buret tips is immersed in the filled beaker on rack position 11 at the beginning of the sample series and after each determination. At the same time, the rinsing solution is stirred. The electrode is immersed on the special beaker position after the sample series.

If required, auxiliary solution can be added prior to the beginning of a determination with a Dosimat/Dosimat plus.

A filled rinsing beaker should be placed on **rack position 11** and on the **special beaker position**.

The individual steps:

Before the first sample:

- Move to rack position 11
- Lower the lift to work position and switch on the stirrer
- Wait for **Rinsing time** to be completed
- Switch off the stirrer and move the lift upward
- Wait for **Dripping time** to be completed

For each sample:

- Move to the sample
- Lower the lift to work position
- Start the determination
- If necessary, initiate dosing (**Activation pulse**) and switch on the stirrer
- If necessary, wait for **Start delay time**
- Wait for end of determination
- Switch off the stirrer and move the lift upward

- Wait for **Dripping time** to be completed
- Move to rack position 11
- Lower the lift to work position and switch on the stirrer
- Wait for **Rinsing time** to be completed
- Switch off the stirrer and move the lift upward
- Wait for **Dripping time** to be completed

After the last sample:

- Move to special beaker position
- Lower the lift in the rinsing beaker to work position

4.4.3 Double dipping

This automation sequence is suitable for simple determinations.

No pump is required for rinsing or for aspiration of the sample vessels. The titration head with electrode and burette tips is immersed in the filled rinsing beaker on rack position 11 and on the special beaker position after each determination. At the same time, the rinsing solution is stirred.

If required, auxiliary solution can be added prior to the beginning of a determination with a Dosimat/Dosimat plus.

A filled rinsing beaker should be placed on **rack position 11** and on the **special beaker position**.

The individual steps:

- Move to the sample
- Lower the lift to work position
- Start the determination
- If necessary, initiate dosing (**Activation pulse**) and switch on the stirrer
- If necessary, wait for **Start delay time**
- Wait for end of determination
- Switch off the stirrer and move the lift upward
- Wait for **Dripping time** to be completed
- Move to rack position 11
- Lower the lift to work position and switch on the stirrer
- Wait for **Rinsing time** to be completed
- Switch off the stirrer and move the lift upward
- Wait for **Dripping time** to be completed
- Moving to special beaker position
- Lower the lift to work position and switch on the stirrer
- Wait for **Rinsing time** to be completed
- Switch off the stirrer and move the lift upward
- Wait for **Dripping time** to be completed

After the last sample:

- Lower the lift in the rinsing beaker to work position



4.4.4 Rinsing in sample

This automation sequence requires an 843 Pump Station for rinsing and aspirating. The sample solution is aspirated after each determination. The titration head with electrode and buret tips is subsequently rinsed in the sample vessel. The rinsing solution is also aspirated.

If required, auxiliary solution can be added prior to the determination with a Dosimat/Dosimat plus.

A filled rinsing beaker should be placed on the **special beaker position**.

The individual steps:

- Move to the sample
- Lower the lift to work position
- Start the determination
- If necessary, initiate dosing (**Activation pulse**) and switch on the stirrer
- If necessary, wait for **Start delay time**
- Wait for end of determination
- Switch off the stirrer and switch on the aspiration pump
- Wait for **Aspiration time** to be completed, the aspiration pump remains switched on
- Switch on the rinsing pump and wait for **Rinsing time** to be completed.
- Switch off the rinsing pump and once again wait for **Aspiration time** to be completed
- Switch off the aspiration pump and move the lift upward
- Wait for **Dripping time** to be completed

After the last sample:

- Moving to special beaker position
- Lower the lift to work position

4.4.5 Rinsing in special

This automation sequence requires an 843 Pump Station for rinsing and aspirating. The titration head with electrode and buret tips is rinsed in the rinsing beaker after each determination. The rinsing solution is aspirated at the same time.

If required, auxiliary solution can be added prior to the determination with a Dosimat/Dosimat plus.

An empty rinsing beaker should be placed on the **special beaker position**.

The individual steps:

- Move to the sample
- Lower the lift to work position

- Start the determination
- If necessary, initiate dosing (**Activation pulse**) and switch on the stirrer
- If necessary, wait for **Start delay time**
- Wait for end of determination
- Switch off the stirrer and move the lift upward
- Wait for **Dripping time** to be completed
- Moving to special beaker position
- Lower the lift to work position
- Switch on the rinsing pump and the aspiration pump
- Wait for **Rinsing time** to be completed, the aspiration pump remains switched on
- Switch off the rinsing pump and wait for **Aspiration time** to be completed
- Switch off the aspiration pump and move the lift upward
- Wait for **Dripping time** to be completed

After the last sample:

- Lower the lift in the rinsing beaker to work position
- Switch on the rinsing pump and wait for **Rinsing time** to be completed.
- Switch off the rinsing pump

4.4.6 Pump control

Rinsing the electrode and aspirating the sample vessels is carried out with the aid of an 843 Pump Station with two membrane or peristaltic pumps. These are connected to the 862 Compact Titrosampler by means of a remote cable (see Chapter 3.10.2, page 24). The pumps can be operated manually by push-button or controlled by means of remote lines.

The method runs of the 862 Compact Titrosampler automatically switch the pumps on or off at a given time. The sequences cannot be modified.

You can define the duration of the rinsing and aspiration procedures under **Menu ▶ Parameters ▶ Automation**, see page 133ff.



Note

The pumps of the 843 Pump Station **cannot be stopped manually** on the 862 Compact Titrosampler. In the event of an **Emergency stop**, switch off the 843 Pump Station with the red **mains switch**.



4.4.7 Dosing auxiliary solutions

The addition of an auxiliary solution can be carried out with a 6xx/7xx Dosimat or a Dosimat plus. This is connected via remote cable to the 862 Compact Titrosampler .

The triggering of the dosing is accomplished by switching on the **Activation pulse** which is put out at the beginning of a titration. The dosing proceeds automatically and is not monitored by the 862 Compact Titrosampler. A waiting time must be observed in each case for the duration of the dosing. Define a sufficiently long **Start delay time**. Both settings can be found under **Menu ▶ Parameters ▶ Start conditions**, see *page 130 ff.* Select a start delay time that is sufficiently large so that the entire volume is dosed before the titration begins.



Note

The dosing of an auxiliary solution is parameterized on the Dosimat. An **8xx Dosimat plus** must be operated in **XDOS** mode, a **6xx** or **7xx Dosimat** in **DIS** mode. Enter the dosing volume on the Dosimat and select a dosing rate as high as possible.

5 Operation

5.1 Switching the instrument on and off

Switching on the instrument

Proceed as follows:



- 1 ■ Press the red **[STOP]** key.
The instrument is initialized and a system test performed. This process takes some time.
- If a buret unit has been attached, then a request appears to carry out the **PREP** function:

```

Buret unit                                010-126
-----
! The buret unit should be
  prepared.
  ◆ Execute a PREP function.
-----
OK
  
```

All tubings and the cylinder are rinsed with the **PREP** (Preparing) function. The preparing of the buret unit is described in chapter "Preparing the buret unit (PREP)", page 62.

- Confirm the message with **[OK]**.
The display of this message can be deactivated in the system settings (see "PREP warning", page 71).

The main dialog is displayed:

```

>Menu                                     ready
Method                                    DET
ID1
ID2
Sample size                               1.0
Unit                                       9
-----
Current sample 0 of table(0)
  
```



Switching off the instrument

The instrument is switched off with the **[STOP]** key. The fact that the key needs to be pressed down for an extended time prevents accidental switch off.

Proceed as follows:

- 1 ■ Keep the red **[STOP]** key pressed down for at least 3 s.
A progress bar is displayed. If the key is released during this time, then the instrument will not be switched off.

5.2 Fundamentals of operation

5.2.1 The keypad




Figure 28 Keypad 862 Compact Titrosampler

BACK	Apply the input and exit the dialog.
↑ ↓	Move the selection bar either up or down by one line at a time. Select the character to be entered in the text editor.
← →	Select the character to be entered in the text and number editor. Select the individual functions in the function bar.
OK	Confirm the selection.
STOP	Stop an ongoing method run or a manual function. Switch the instrument on/off.
START	Start a method run or a manual function.

5.2.2 Structure of the dialog windows

```

System                               ready
>Settings
>Lift
>Sensors
>Solutions
>Common variables
>File management
>External devices

```

The current dialog title is displayed on the left-hand side of the title line. The current status of the system is displayed in the upper right-hand corner.

ready The instrument is in normal state.

busy A method has been started.

hold A method has been paused.

Some dialogs have a so-called function bar on the bottom line. The functions contained therein can be selected with the arrow keys [←] or [→] and executed with [OK].

```

Solution list                         ready
Reagent 1                             #IDU
Reagent 2                              DU
Edit New Delete

```

5.2.3 Navigating in the dialog

The selection bar is displayed in inverted style. Use the arrow keys [↑] and [↓] to move the selection bar upward or downward one line at a time. If a dialog text is marked with ">", then additional settings are available in a subordinate dialog. Use [OK] to access this dialog.

Example: System settings

```


System                               ready
>Settings
>Lift
>Sensors
>Solutions
>Common variables
>File management
>External devices

```

Use the [BACK] key to return to the next higher level.



5.2.4 Entering text and numbers



In the editing dialog for text or numerical input you can select the individual characters with the arrow keys. Use **[OK]** to apply the character in the input field. The following functions are available:

Editing function	Description
Accept	The modification is applied and the editing dialog is exited.
Cancel	The editing dialog is exited without applying the modification.
Delete	The contents of the input field is deleted completely.
[+-]	The character left of the cursor is deleted (back-space).
+ - 	Text editor only The cursor within the input field is shifted to the left by one character each time that [OK] is pressed.
 ->	Text editor only The cursor within the input field is shifted to the right by one character each time that [OK] is pressed.
[BACK]	The modification is applied and the editing dialog is exited.

The **[BACK]** key has the same function as **Accept**.

A commercially available USB keyboard can be connected to make it easier to enter text and numbers. The assignment of the keys on the PC keyboard is described in *chapter 10.4.2, page 142*.

5.2.5 Selecting from a selection list

```
Unit
█
mg
µg
mL
µL
pieces
>User-defined
```

In a selection list, select the individual entries with the arrow keys [↑] and [↓]. Accept the selection with [OK] or [BACK].

5.3 Formula editor

The formulas for the calculations are entered with the formula editor. The formula editor is equipped with an automatic syntax check, which is activated when applying the formula. The generally valid rules of priority apply for the calculation operations.

```
R1=
█
0123456789
.+-*/()
C00 EP# CI# R# FP# CV0# SMN#
TITER CONC Var Templates
Accept Cancel Clear [+ -]
+-| | ->
```

Variable	Description
C00	Sample size
EP#	Volume of endpoint EP# (# = 1...9)
CI#	Sample identification (# = 1...2)
R#	Result (# = 1...5)
FP#	Volume of fixed endpoint FP# (# = 1...9)
CV0#	Common variable (# = 1...5)
SMN#	Mean value of result R# (# = 1...5)
TITER	Titer of selected solution
CONC	Concentration of selected solution
Var	List of additional variables (see "Variables", page 40)
Templates	List of predefined calculation formulas (see "Calculation templates", page 40)



"#" stands for a sequential number that you must enter manually. Example: when applying the variable **EP#** in the formula, only **EP** is entered. You will still need to enter the number yourself.

The meanings of the editing functions are explained in *chapter 5.2.4, page 38*.

Variables

Pressing **Var** causes a list to be displayed containing additional variables. You can enter these variables either directly into the formula or also by selecting them from the list and applying them with **[OK]**.

Variable	Description
MIM	Initial measured value, i.e. measured value prior to the processing of the start conditions
MSM	Start measured value, i.e. measured value after the processing of the start conditions
MCV	End volume, i.e. total dosed volume at the end of the titration
ET#	Temperature at endpoint EP# (# = 1...9)
EM#	Measured value of endpoint EP# (# = 1...9)
ED#	Time at endpoint EP# (# = 1...9)
MSV	Start volume
MEN	Electrode zero point pH(0)
MSL	Electrode slope
DD	Duration of the entire determination
MST	Start temperature
MCT	End temperature
FT#	Temperature at fixed endpoint FP# (# = 1...9)
FM#	Measured value of fixed endpoint FP# (# = 1...9)
FD#	Time at fixed endpoint FP# (# = 1...9)

For Molw, see the following section.

Calculation templates

Pressing **Templates** displays a list with calculation templates. You can apply these templates directly with **[OK]**.



Note

Some templates contain the wildcard **Molw**, which stands for the molar mass of the sample. You must replace this wildcard with the correct value in the calculation formula.

The templates available:

Template	Description
Content %	Content in % Unit of the sample size = g
Content mmol/L	Content in mmol/L Unit of the sample size = mL
Content mol/L	Content in mol/L Unit of the sample size = mL
Content g/L	Content in g/L Unit of the sample size = mL
Content ppm	Content in ppm Unit of the sample size = g
Titer	Titer calculation Unit of the sample size = g
Blank mean value	Blank value as mean value of single results
Blank single value	Blank value as single value

5.4 Methods

5.4.1 Creating a new method

Proceed as follows to create a new method:

1 Open the method table


- In the main dialog, select **Method** and press **[OK]**.

The method table opens:



2 Select a titration and measuring mode

- In the function bar, select **New** and press **[OK]**.



- Select **Mode** and press **[OK]**.
- Select the desired titration mode in the selection list and apply with **[OK]**.
- Select **Measured quantity** and press **[OK]**.
- Select the desired measuring mode in the selection list and apply with **[OK]**.

3 Select an automation sequence

- Select **Automation** and press **[OK]**.
- Select the desired sequence in the selection list and apply with **[OK]**.
- Press **[BACK]**.

The method is now loaded and is displayed in the main dialog under **Method**.

If a new method has been created, then the individual parameters can be modified under **Menu ▶ Parameters**.

5.4.2 Saving a method

If you modify method parameters, then you can save these as a your own method. A maximum of 100 methods can be saved.

To save a method, proceed as follows:

1 Open the method table


- In the main dialog, select **Method** and press **[OK]**.

The method table opens:



2 Modify/apply the method name

- In the function bar, select **Store** and press **[OK]**.
A method name will be suggested for new methods. If the method has already been saved once, then the method name will be displayed:



Apply the name:

- Press **[BACK]**.

The method will be saved and the method table is displayed.

Enter a new name:

- Press **[OK]**.
The text editor opens.
- Enter a method name (max. 12 characters) and apply with **Accept** or **[BACK]**.
- Press **[BACK]**.

The method will be saved and the method table is displayed.

5.4.3 Loading a method

To load a method, proceed as follows:

1 Open the method table

- In the main dialog, select **Method** and press **[OK]**.

The method table with the stored methods opens:



2 Select a method

- Select the desired method.

3 Load the method

- In the function bar, select **Load** and press **[OK]**.

The method is now loaded and is displayed in the main dialog under **Method**.

5.4.4 Exporting a method

Methods can be exported to a connected USB stick.



Note

This function is possible only if a USB stick is connected as an external storage medium.

To export a method, proceed as follows:

1 Open the method table

- In the main dialog, select **Method** and press **[OK]**.

The method table with the stored methods opens:



2 Select a method

- Select the desired method.

3 Export the method

- In the function bar, select **Export** and press **[OK]**.

The method is exported. The directory structure on the USB stick is listed in *chapter 6.6, page 79*.

5.5 Sample data

You can enter the sample data (identification, sample size, etc.) in a variety of ways:

- Use of the sample table. This is particularly useful with sample series.
- Directly in the main dialog, if the same sample data is to be used for an entire sample series.
- Automatic request immediately after the start of the determination. This is only useful with single determinations.

You can also send the sample size and the unit from a connected balance in any case. With some balances, the sample identification and the method can be sent in addition to the sample size.

5.5.1 Sample table

The sample table is a table in which the sample data for up to 99 samples can be entered. The sample data can also be entered while a determination is running (*see Chapter 5.7.2, page 54*).

The sample table contains numbered lines. The identification (**ID1**) and the sample size of each sample are displayed.

Sample table		ready
1	#8805923 1.0 g	
2	#8805923 1.0 g	
3	#8805924 1.0 g	
4	#8805924 1.0 g	
5	...	
Edit Delete Insert New		

Edit

Edit data of the selected line

Delete

Delete the selected line from the sample table.

Insert

Insert a new line above the line selected.



New

Delete the sample table completely. This function is visible only if the instrument is in **ready** status.

Editing the sample data

Sample data		ready
Method	Me4155	
ID1	#8805923	
ID2		
Sample size	1.0	
Unit	g	

←- Line 1 of 4 -→

You will see at the very bottom the line number of the selected line and the line number of the last line containing data. In this example, the first line is opened and the sample table contains four lines.

One can scroll between the individual data sets with the keys [←] and [→].

Inserting a new line

If you find yourself on the last line (i.e. **Line 4 of 4** in the above example), you can add a new line to the sample table by pressing [→] again. The sample data of the previous sample will be applied thereby.

Method

Method used for processing the sample.

Selection	Selection of stored methods empty
Default value	empty

empty

The currently loaded method is used.

ID1

Sample identification. The sample identification can be used in calculations as the variable **CI1**.

Input	max. 10 characters
Default value	empty

ID2

Sample identification. The sample identification can be used in calculations as the variable **CI2**.

Input	max. 10 characters
Default value	empty

Sample size

Sample size. The value of the sample size can be used in calculations as the variable **C00**.

Range	-999999999 ... 999999999
Default value	1.0

Unit

Unit of sample size.

Selection	g mg µg mL µL pieces User-defined
Default value	g

User-defined

A user-defined unit can be created. This will be added to the selection list. The previous entry will be overwritten as soon as the new unit has been defined.

Sending the sample size from a balance



Note

In order for a balance to be able to send a sample size to the sample table, the sample table has to be switched on in the sample series dialog (see "Starting a sample series", page 49).

Sample table		ready
1	#8805923 1.0 g	
2	#8805923 1.0 g	
3	#8805924 1.0 g	
4	#8805924 1.0 g	
5	...	

Edit Delete Insert New

If the sample size is sent directly from the balance, then it will always be entered in a new line at the end of the sample table. It does not matter which line is highlighted or whether the sample table is even opened. The sample size is entered in line 5 in the above example.



Note

If you would like to enter the sample size in a particular line, then you must open the corresponding editing dialog (i.e. the **Sample data** dialog is displayed).

If the editing dialog for the sample size is opened, then the sent value will be ignored.



5.5.2 Entering sample data in the main dialog

In the main dialog you can enter the sample data even while a determination is running (see Chapter 5.7, page 53). It will be used for the running determination.

Menu	ready
Method	DET
ID1	
ID2	
Sample size	1.0
Unit	g
Current sample 0 of table(0)	

ID1

Sample identification. The sample identification can be used in calculations as the variable **CI1**.

Input	max. 10 characters
Default value	empty

ID2

Sample identification. The sample identification can be used in calculations as the variable **CI2**.

Input	max. 10 characters
Default value	empty

Sample size

Sample size. The value of the sample size can be used in calculations as the variable **CO0**.

Range	-999999999 ... 999999999
Default value	1.0

Unit

Unit of sample size.

Selection	g mg µg mL µL pieces User-defined
Default value	g

User-defined

A user-defined unit can be created. This will be added to the selection list. The previous entry will be overwritten as soon as the new unit has been defined.

5.6 Performing a sample series

Samples can be placed anywhere on the rack. They are processed in accordance with ascending rack position.

The following is to be observed:

- In addition to the sample vessels, a rinsing beaker has to be placed on the last rack position, marked with the sign ▲. This vessel must be either empty or filled with a rinsing solution, depending on the automation sequence, *see chapter 4.4, page 29ff*.

5.6.1 Starting the sample series

Starting a sample series

A suitable method must be loaded before a sample series is started (*see Chapter 5.4.3, page 43*). The necessary parameters (*see chapter 7.1, page 85ff*) can then be modified.

If the sample table is used and if it contains methods defined, these methods will be used. In this case, previously loading a certain method is not necessary.



1 Define the sample series

Press the **[START]** key.

```

Sample series          ready
Number of samples     table
Next sample pos.      1
Sample table           on

Press [START] key to continue
  
```

You can now select the quantity and the first rack position of the samples to be processed as well as the location of the sample data.

2 Enter the number of samples

- Select **Number of samples** and press **[OK]**.
- Enter the number of samples. **Table** means that all samples in the sample table will be processed until it will be empty.
- Close the entry dialog with **[BACK]** or **Accept**.

3 Enter the rack position of the first sample

- Select **Next sample pos.** and press **[OK]**.



- Enter the starting position of the sample series.
- Close the entry dialog with **[BACK]** or **Accept**.

The value for the number of samples remains saved for the next sample series. The position of the first sample is increased with each method run.

You can still cancel the start of the sample series at this time with **[BACK]** or **[STOP]**.

4 Activate or deactivate the sample table

If the sample table is activated, the sample data of the sample table is used. If the sample table is deactivated the sample data of the main dialog is used.

5 Close the sample series dialog

Close the dialog with the **[BACK]** key.



6 Start the sample series

Press the **[START]** key.

Stopping a sample series

A sample series can be canceled at any time.

A sample series can be canceled at any time. When this is done, instruments connected via remote connections, such as a Dosimat plus or an 843 Pump Station will also be stopped.



1 Press the **[STOP]** key.

The method run is stopped. The sample series cannot be resumed.

5.6.2 Pausing a sample series and continuing

Interrupting a sample series

A method run of the 862 Compact Titrosampler can be paused and then continued again. The connected instruments are however **not** stopped.

```
Me4155          busw
Move to sample
Hold Stirrer
```

A function bar with the entry "**Hold**" is displayed during the run of a sample series in the so-called "Live" dialog.

- 1 Press the **[OK]** key.

```
Me4155 hold
Move to sample
Continue Stirrer
```

The method run is paused. However, currently running movements of the sample rack or the lift will be finished.

Instead of the "**Hold**" function, "**Continue**" is displayed in the function bar.

Continuing a sample series

If a method run is paused, then the "**Hold**" status is displayed in the title bar, see previous illustration. The sequence can be continued with the "**Continue**" function.

In the "**Hold**" status, a method run can be stopped completely, and with it the entire sample series, by pressing the **[STOP]** key.

- 1 Press the **[OK]** key.

As is also the case at the start of a sample series, a request dialog appears here in which the number of samples to be processed can still be changed. It is thus possible to shorten a sample series or to extend it, without stopping it.

```
Sample series hold
Number of samples 99
Press [START] key to continue
```

- 2 Press the **[OK]** key and enter the number of samples that still need to be processed. The current sample must be taken into account.
- 3 Press the **[START]** key.





The sample series continues.

5.6.3 Status of the sample series

The status of the sample series is continuously displayed in the main dialog. In addition to the name of the method loaded and the sample data of the current sample, information concerning the running or the previous sample series will be displayed on the lowest line. This includes:

- The number of samples already processed (including the current sample)
- Total number of samples in the sample series
- Number of assigned lines in the sample table (only when sample table is switched on)

Examples:

Current sample 2 of 5(3)

The second sample in a sample series with a total of 5 samples is currently being processed. 3 unprocessed samples remain in the sample table.

Current sample 1 from table(10)

The first sample in a sample series which contains all of the samples in the sample table is currently being processed. 10 unprocessed samples remain in the sample table.

One can use the key **[BACK]** to switch back and forth between live display and main dialog while a sample series is running. This makes it possible to check the current status of the sample series at any time. If changes are made in the sample table, then this status bar is updated without delay in the main dialog.

5.6.4 Special case: Calibrations

The electrodes must be calibrated in advance for a correct SET titration. This can be accomplished in the same sample series.

- Create a calibration method (CAL mode) in which you define the number of buffer solutions to be processed. Save the method.
- Add a line to the sample table which calls up the calibration method. It should be positioned immediately ahead of the lines with the data for the samples to be titrated.
- Place the required buffer solutions on the rack positions ahead of the sample solutions on the sample rack.

The calibration method is run first when the sample series is being executed. The buffer solutions are approached automatically.



Note

Enter the rack position of the first buffer solution as the **Next sample pos.** at the start of the sample series.

When specifying the **Number of samples**, enter (for the calibration) one sample more than is to be titrated. The calibration is considered to be one single sample processing, even if several buffer solutions are being processed. The number of lines in the sample table is the determining factor. You can also select **Table** as the **Number of samples**.

5.7 Live modifications

5.7.1 Editing the sample data of the running determination

The sample data can be entered or modified in the main dialog while a determination is running. In calculations always the sample data entered at the end of the titration in the main dialog is used.

Proceed as follows to edit the sample data:

1 Display the main dialog

- Press **[BACK]**.

The main dialog is displayed. The determination continues to run in the background.

2 Edit the sample data

- Edit the sample data and accept with **Accept** or **[BACK]**.

3 Display the live dialog

- Press **[BACK]**.

or

- Select **Menu** and press **[OK]**.

```

Menu                                     busy
>Results
>Live dialog
>Parameters
>Sample table
  
```

- Select the menu item **Live dialog** and press **[OK]**.

The live dialog is displayed once again.



Note

If the determination is finished while the editing dialog is opened (e.g. of the sample size), then this will be closed automatically and the result dialog will be displayed. The value entered must be entered once more and the determination must be recalculated.

Make sure that the editing dialogs are closed before the determination is finished.

5.7.2 Editing the sample table while a determination is running

You can insert new lines or delete existing ones or edit sample data while a determination is running.



Note

We recommend that the editing dialogs always be closed in order to ensure that no problems occur during the run and that the current data is always available for calculation purposes.

Editing the sample table

Proceed as follows to edit the sample table:

1 Display the main dialog

- Press **[BACK]**.

The main dialog is displayed. The determination continues to run in the background.

2 Open the main menu

- Select **Menu** and press **[OK]**.

```

Menu                                     busy
>Results
>Live dialog
>Parameters
>Sample table
  
```

3 Select the sample data

- Select the menu item **Sample table** and press **[OK]**.
- Select the desired line.

- In the function bar, select **Edit** and press **[OK]**.

4 Edit the sample data

- Edit the sample data and accept with **Accept** or **[BACK]**.



Note

In addition to the sample data, the method can also be modified, except in cases where the determination is running.

5 Display the live dialog

- Select the menu item **Live dialog** in the main menu and press **[OK]**.

or

- Press **[BACK]** in the main dialog.

The live dialog is displayed once again.

Editing the sample data of the running determination

When you use the sample table, the editing of the sample data of the running determination proceeds as described in *chapter 5.7.1, page 53*. In addition, you have the option of editing these in the sample table. The first line always contains the sample data of the running determination. Simply select the menu item **Sample table** (see "*Editing the sample table*", page 54) for this purpose.

5.7.3 Live parameters

Certain method parameters can be edited while a determination is being carried out. The only parameters that can be modified are those that can be selected. Nevertheless, all of the parameters are visible. The modified parameters are taken into account at once. If you modify e.g. the start conditions after the start volume has been dosed, then these modifications will not be taken into account until the next determination.

Proceed as follows to edit the parameters:

1 Display the main dialog


- Press **[BACK]**.

The main dialog is displayed. The determination continues to run in the background.



2 Open the main menu

- Select **Menu** and press **[OK]**.



3 Edit the method parameters

- Select the menu item **Parameters** and press **[OK]**.
- Change the desired parameters accordingly.

4 Display the live dialog

- Select the menu item **Live dialog** in the main menu and press **[OK]**.

or

- Press **[BACK]** in the main dialog.

The live dialog is displayed once again.

5.8 Results

Menu ► Results

After the completion of the sample series, the results dialog is displayed.



Note

It is also possible to display the results dialog during a running sample series. The statistics results of the previous determinations are also available. However, the results of the current determination can only be displayed between the end of one titration and the beginning of the next titration. The method run can be interrupted with the **Hold** function in the live dialog.

Recalculations and the display of the titration curves in the results window are not possible during a sample series.

Results	ready
Content	10.3 %
EP1	pH 7.499 10.0000 mL
	72.0 s ERC 85.7
Stop volume reached	
Curve Recalc Statistics	

The calculated result of the last determination and details concerning the endpoint are shown in the overview.

Curve

Display the curve of the current determination.

Recalc


Recalculate the current determination. The procedure will be executed immediately.

Statistics

Display the statistical overview of a determination series (*see Chapter 5.9, page 58*).

Displaying the curve

The curve of the current determination can be displayed with the **Curve** function.



The arrow keys [←] and [→] can be used to move to the individual measuring points. A cross hair is used to show the current position on the curve. The data (volume, measured value, time, etc.) for the respective measuring point is indicated on the right-hand side.

Recalculating



Note

Recalculation cannot be undone.

All of the results are recalculated with the **Recalc** function. This is necessary if, for example, the calculation, the titer or the sample size has been modified.

The result and the sample size of each determination are shown.

On/Off

Remove the selected determination from the statistics. The line will then be marked with an asterisk (*), the statistics will be recalculated automatically. If several calculations are defined in the method, then all the results will be removed from the statistics.

Deleting statistical data

All statistical data is deleted with the **Reset** function. The statistics data is deleted automatically in following cases:

- when all of the determinations of the determination series have been carried out and a new determination has been started afterwards.
- when a new method is loaded.

Adding a determination to a determination series

In case that, for example, a determination was faulty and had to be removed from the statistics, an additional sample can be added to a determination series with the **Increase** function. The second number in the **Statistics** line will be increased automatically by one.

5.10 Printing a report manually

Menu ► Print reports

To print a report manually, proceed as follows:

1 Open the main menu

- In the main dialog, select **Menu** and press **[OK]**.

```


Menu                                ready
>Manual control
>Results
>Parameters
>Sample table
>System
>Print reports

```

2 Open the print dialog

- Select the menu item **Print reports** and press **[OK]**.

The dialog window with the available reports opens:



3 Select a report

- Select the desired report and press **[OK]**.

The report is printed out.

The following reports can be printed out manually:

Results	Result report with determination properties, sample data, calculated results, etc.
Curve	Curve report. The width of the curve is defined in the system settings (see "Graphics width", page 81).
Measuring point list	Measuring point list report.
Parameters	Report with all method parameters of the loaded method.
System	System report with system settings, solution list, external devices, etc.
Calculations/Statistics	Calculation report. The statistics are also printed out when there are multiple determinations. The individual determinations with the respective sample size will be printed out for each result, along with the mean value and the absolute and the relative standard deviation.
Report as in method	The reports that are defined in the method will be printed out.
PC/LIMS	Machine-readable report with all of the data for a determination. This report can be saved as a TXT file on a connected USB stick or sent to a terminal program or to a LIMS via an RS-232 interface. The definition is made in the system settings (see "PC/LIMS report", page 80).

5.11 Manual control

Menu ► Manual control

The following functions are available in the manual control:

- Move the sample rack (**Rack position**)
- Move the lift (**Lift position**)
- Dosing (**Dosing**)
- Measure (**Measure**)
- Stir (**Stir**)

Manual control		ready
Rack position		1
Lift position		work pos.
Dosing		
Measure		
Stirrer	off	Rate 8
PREP DOS ADD EMPTY		

The available subfunctions are listed for each function in the function bar.

5.11.1 Rotating the sample rack

If the **Rack position** line is selected, then the arrow keys [⇒] and [⇐] can be used to select one of the following functions, which can then be executed by pressing [OK]:

- Next** The lift is moved upward and the next-higher rack position is placed in front of the lift.
- If the [OK] key remains pressed, the rack automatically moves to the next position.
- Previous** The lift is moved upward and the next-lower rack position is placed in front of the lift.
- If the [OK] key remains pressed, the rack automatically moves to the next position.
- Reset** The rack is initialized. The lift is moved upward and the sample rack is rotated to the starting position. At the same time, the starting position (**Next sample pos.**) is reset to **1** for the start of the next sample series.

The rack position display is always updated as soon as the rack is in the new position.



5.11.2 Moving the lift

If the **Lift position** line is selected, then the lift can be moved to the position suggested in the function bar by pressing **[OK]**. Only two positions are possible:

Work pos.	The working height. It can be set under Menu ▶ System ▶ Lift (see page 72).
Shift pos.	The rotating height. The lift moves all the way to the top.

The current lift position is displayed. The respective other possible position is offered in the function bar.

5.11.3 Dosing

The following dosing functions are available in the manual control:

Prepare buret unit (PREP)	Rinse the cylinder and tubings of the buret unit (see Chapter 10.1.2, page 139).
Continuous dosing (DOS)	Dose while the [START] key is pressed.
Dose fixed volume (ADD)	Dose a specified volume.
Empty dosing device (EMPTY)	Empty the cylinder and tubings of the buret unit (see Chapter 10.1.2, page 139).

Preparing the buret unit (PREP)

The **PREP** function is used to rinse the cylinder and tubings of the buret unit and fill them air bubble-free. You should carry out this function before the first determination or once per day.

Proceed as follows:

1 Open the manual control

- In the main dialog, select **Menu** and press **[OK]**.
The main menu opens.
- Select the menu item **Manual control** and press **[OK]**.

Manual control opens.

2 Select the dosing function

- Select the **Dosing** item.

```
Manual control      ready
Rack position      1
Lift position      work pos.
Dosing
Measure
Stirrer      off      Rate 8
PREP DOS ADD EMPTY
```

- In the function bar, select **PREP** and press **[OK]**.

The following message is displayed:

```
Splash warning      010-132
! Check the buret tip. It
  should point into a vessel.
  Do you want to continue?
Yes No
```

3 Start the preparing



Caution

Make sure that the buret tip is directed into a vessel that can accommodate the volume of your buret unit several times over.

- Select **Yes** and confirm the message with **[OK]**.

Preparing is carried out.

Continuous dosing (DOS)

Continuous dosing will be carried out with the **DOS** function for as long as you keep the **[START]** key pressed down.

Proceed as follows:

1 Open the manual control

- In the main dialog, select **Menu** and press **[OK]**.
The main menu opens.
- Select the menu item **Manual control** and press **[OK]**.

Manual control opens.



2 Select the dosing function

- Select the **Dosing** item.

```
Manual control          ready
Rack position          1
Lift position          work pos.
Dosing
Measure
Stirrer      off      Rate 8
PREP DOS ADD EMPTY
```

- In the function bar, select **DOS** and press **[OK]**.

```
Dosing          ready
Dosing rate     max. mL/min
Filling rate    max. mL/min

Press [START] key
```

3 Configure the dosing function



Note

- The dosing and filling rate should be decreased for viscous liquids.
- The maximum dosing and filling rate depends on the cylinder volume (see Chapter 10.1.1, page 139).

- Enter the dosing rate.
- Enter the filling rate.

4 Start dosing

- Press **[START]**.

The status changes to **busy**, the dosed volume is displayed. When a cylinder volume has been added, the dosing cylinder will be refilled automatically.

5 Fill the cylinder

- Press **[STOP]** or **[BACK]**.

The dosing cylinder is filled. If you start the filling with **[BACK]**, then the dialog will also be exited.

Dosing a particular volume (ADD)

You can dose a particular volume with the **ADD** function.

Proceed as follows:

1 Open the manual control

- In the main dialog, select **Menu** and press **[OK]**.
The main menu opens.
- Select the menu item **Manual control** and press **[OK]**.

Manual control opens.

2 Select the dosing function

- Select the **Dosing** item.

```
Manual control      ready
Rack position      1
Lift position      work pos.
Dosing
Measure
Stirrer           off           Rate 8
PREP DOS [ADD] EMPTY
```

- In the function bar, select **ADD** and press **[OK]**.

```
Dosing      ready
Volume      10 mL
Dosing rate max. mL/min
Filling rate max. mL/min
Press [START] key
```

3 Configure the dosing function



Note

- The dosing and filling rate should be decreased for viscous liquids.
- The maximum dosing and filling rate depends on the cylinder volume (see Chapter 10.1.1, page 139).
- Enter the desired volume.
- Enter the dosing rate.
- Enter the filling rate.



4 Start dosing

- Press **[START]**.

The status changes to **busy**, the dosed volume is displayed. When a cylinder volume has been added, the dosing cylinder will be refilled automatically.

5 Fill the cylinder

- Press **[STOP]** or **[BACK]**.

The dosing cylinder is filled. If you start the filling with **[BACK]**, then the dialog will also be exited.

5.11.4 Measuring

Open the dialog for manual measurement as follows:

1 Open the manual control

- In the main dialog, select **Menu** and press **[OK]**.
The main menu opens.
- Select the menu item **Manual control** and press **[OK]**.

Manual control opens.

2 Select a measuring mode

- Select the **Measure** item.

```
Manual control          ready
Rack position          1
Lift position          work pos.
Dosing
Measure
Stirrer      off      Rate 8
pH U
```

- In the function bar, select the measuring mode and press **[OK]**.

```
Measure          ready
Electrode        pH electrode
Temperature      25.0 °C


Press [START] key
```

3 Configure the measuring mode

- Select the desired electrode from the sensor list.
The selection depends on the measuring mode. Sensors are defined under **System ▶ Sensors**.
- Enter the measuring temperature if no temperature sensor is connected. If a temperature sensor is connected then the temperature will be measured automatically.
This temperature is used for automatic temperature compensation with pH measurements.

4 Start the measurement

- Press **[START]**.



```

Measure                                     busy
*** pH      6.362 ***
***          20.5 °C ***
  
```

The status changes to **busy**. The current measured value and the measuring temperature are displayed.

5 Stop the measurement

- Press **[STOP]** or **[BACK]**.

The measurement is stopped. The status changes back again to **ready**. If you stop the measurement with **[BACK]**, then the dialog will also be exited.

5.11.5 Stirring

You can control a connected stirrer manually.

Proceed as follows:

1 Open the manual control

- In the main dialog, select **Menu** and press **[OK]**.
The main menu opens.
- Select the menu item **Manual control** and press **[OK]**.

Manual control opens.



2 Set the stirring rate

- Select the **Stirrer** item.

Manual control		ready
Rack position		1
Lift position		work pos.
Dosing		
Measure		
Stirrer	off	Rate 8
On Stir- Stir+		

- In the function bar, select **Stir-** or **Stir+**.
The stirring rate will be increased or decreased by one step each time the **[OK]** key is pressed.
The algebraic sign changes the direction in which the stirring is done. When the stirrer is viewed from above, this means:
 - "+": counterclockwise rotation
 - "-": clockwise rotation

3 Switch on the stirrer

- In the function bar, select **On** and press **[OK]**.

The stirrer is started and stirs at the rate which has been set. **Off** is now displayed in the function bar.

4 Switch off the stirrer

- In the function bar, select **Off** and press **[OK]**.

The stirrer is stopped.

6 System settings

6.1 Basic settings

Menu ▶ System ▶ Settings

This chapter contains a description of general instrument settings.

User name

A user name can be entered here for the report. This parameter will only be printed if a user has been defined.

Input	max. 12 characters
Default value	empty

Instrument name

An instrument name for the report can be entered here. This parameter will only be printed if a designation has been defined.

Input	max. 10 characters
Default value	empty

Serial number

Serial number of the instrument. This is printed as a component of the instrument identification in the report header.

Program version

Version number of the instrument software. This is printed as a component of the instrument identification in the report header.

Time

Current time. Only valid numbers can be entered.

Format: hh:mm:ss

Date

Current date. Only valid numbers can be entered.

Format: YYYY:MM:DD

Language

Setting the dialog language. In addition to English a further language can be selected.



Note

A second language must be installed in advance in order to be able to select it here. The installation may only be carried out by competent personnel.

Dialog type

The user dialog can be limited for routine operations. One can operate normally with methods in the limited dialog. However, no settings can be made or methods deleted.

The resetting of the dialog does not take effect until the main menu is exited.

The limitation of the dialog has the following effects:

- The menu items **System** and **Parameters** are not shown in the main menu.
- Methods can only be loaded, but not deleted, exported or created.



Note

If the limited dialog is activated for routine operations, then the expert dialog cannot be switched on during running operations. To change the dialog type, the 862 Compact Titrosampler must be switched off and then back on again. The expert dialog can be forced at the time the instrument is started. Then it is possible to enter whatever settings one wishes, e.g. the changing of the dialog type. If the instrument is switched off again without changing the dialog type, then the routine dialog will remain activated.

Forcing the expert dialog:

- Switch on the instrument.
- Wait for the display of the instrument logo with the lettering **easy, safe, precise**.
- Press the **[STOP]** key once again and hold it down while also briefly pressing the **[BACK]** key.
- Release both keys once again.

Selection	Expert Routine
Default value	Expert

Expert

Complete dialog.

Routine

Limited dialog for routine operations.

Contrast

The contrast of the display can be adjusted with the arrow keys [←] and [→].

- [←]: the contrast will be decreased by one step each time the key is pressed.
- [→]: the contrast will be increased by one step each time the key is pressed.

Range	150 ... 240
Default value	212



Note

Alternatively, the contrast can also be modified in the following manner:

Keep the red **[STOP]** key pressed down. As soon as the progress bar appears, also press the arrow key [↓] or [↑] repeatedly.

This method will however cause the contrast to be modified by several steps.

Beep

If this parameter is switched on, then a short beep will be heard in the following cases:

- When a key is pressed.
- At the end of the determination.

Selection	on off
Default value	on

PREP warning

If this parameter is switched on, then the recommendation to carry out the **PREP** (preparing) function will be displayed in the following cases:

- After the instrument has been switched on.
- Each time a buret unit has been attached.

This function causes all tubings and the cylinder to be rinsed (*see Chapter 10.1.2, page 139*).

Selection	on off
Default value	on



Temperature sensor

The instrument supports the use of two different temperature measurement techniques:

- NTC (Negative Temperature Coefficient)
- Pt1000 (Platinum resistance)

Select here the type that has been connected to the instrument. If an NTC sensor is used, then it is also necessary to enter two sensor characteristics. These characteristics are listed in the specifications of the sensor.

Selection	Pt1000 NTC
Default value	Pt1000

R (25 °C)

This parameter is visible only when **Temperature sensor = NTC**.

Nominal resistance of the NTC sensor at 25 °C.

Range	1000 ... 99999 ohm
Default value	30000 ohm

B value

This parameter is visible only when **Temperature sensor = NTC**.

Material constant of the NTC sensor. B values of NTC sensors are frequently based on different reference temperatures (usually 25 °C and 50...100 °C).

Range	1000 ... 9999 K
Default value	4100 K

6.2 Lift settings (Lift)

Menu ▶ System ▶ Lift

Lift	ready
Work position	60 mm
Initial lift pos.	shift pos.
Work pos.	Up Down

Work position

The working height of the lift can be set to the desired value. This is accomplished by means of the direct operation of the lift.

Three functions can be selected from the function bar with [←] and [⇒] and then run by pressing [OK]:

- **Work pos.** moves the lift to the current working height.
- **Up** moves the lift 6 mm upward.
- **Down** moves the lift 6 mm downward.

When this dialog page is exited, the respective current lift position will be applied as **Work pos.**.

Range	0 ... 132 mm (Increment: 6)
Default value	60 mm

Initial lift pos.

After the 862 Compact Titrosampler has been switched on, the lift moves all the way to the top to the home position for the initialization of the drive. It can then be moved back down to the working height if desired.

Selection	shift pos. work pos.
Default value	shift pos.

shift pos.

Resting position (0 mm) all the way up

work pos.

The set working height

6.3 Managing sensors

6.3.1 General

Menu ▶ System ▶ Sensors

Sensor list	ready
pH electrode	
Metal electrode	
Temperature sensor	
Edit New Delete	

Three standard sensors are defined in the sensor list: **pH electrode**, **Metal electrode** and **Temperature sensor**. These sensors cannot be deleted or renamed. The sensor list can contain a maximum of 10 sensors.

Every sensor is identified with an unambiguous name. This means that it is not possible to use the same name twice, e.g. for a pH electrode and for a metal electrode.

Edit

Edit the data of the selected sensor, see following chapter.

**New**

Add a new sensor to the list, see following chapter.

The following sensor types can be selected:

- pH electrode
- Metal electrode
- Temperature sensor
- Other sensor, e.g. Spectrosense

Delete

Delete the selected sensor from the list.

6.3.2 Editing the sensor data**Name**

The designation of the sensor is used for unambiguous identification.

Input	max. 24 characters
Default value	empty

Type

The sensor type is displayed.

Slope

This parameter only appears for pH electrodes.

Slope of the pH electrode. With a 1-point calibration, only pH(0) can be calculated, 100.0 % is used as the slope.

Range	-999.9 ... 999.9 %
Default value	100.0 %

pH(0)

This parameter only appears for pH electrodes.

pH value of the pH electrode at 0 mV. Apart from the slope, pH (0) is the second characteristic of the calibration curve.

Range	-20.000 ... 20.000
Default value	7.000

Calibration temp.

This parameter only appears for pH electrodes.

Temperature at which the last calibration was carried out.

Range	-20.0 ... 150.0 °C
Default value	25.0 °C

Calibration date

This parameter only appears for pH electrodes.

Date of the last calibration.

Monitoring

This parameter only appears for pH electrodes.

Activating and deactivating the calibration monitoring.

Selection	on off
Default value	off

Time interval

This parameter is visible only when **Monitoring = on**.

You will be notified that this time interval (in days) has elapsed when starting a method. You can then select whether or not you would still like to start the method.

Range	1 ... 999 d
Default value	999 d

6.4 Managing solutions

6.4.1 General

Menu ► System ► Solutions

Solutions can be used in intelligent buret units or in non-intelligent buret units. Intelligent buret units have a built-in data chip on which the data for the reagent is stored. This data is automatically read out during attachment and entered in the solution list.

Solution list	ready
Reagent 1	*IDU
Reagent 2	DU

Edit New Delete

The name and the type are specified for each solution in the solution list. The asterisk (*) on the right-hand side indicates that this buret unit is attached (only for intelligent buret units). An unlimited number of solutions in buret units with data chip can be added to the solution list. The number of solutions in buret units without data chip is limited to 10 items.

Meaning of the type:

- **DU**: dosing unit without data chip



- **IDU**: dosing unit with integrated data chip

Edit

Edit the data of the selected solution, see following chapter.

New

Add a new solution to the list, see following chapter.

Delete

Delete the selected solution from the list.

6.4.2 Editing the solution data**Name**

The designation of the solution is used for unambiguous identification.

Input	max. 24 characters
Default value	empty

Type

The model of the buret unit is displayed.

Cylinder volume

Cylinder volume of the buret unit in mL. The cylinder volume is automatically read out with intelligent buret units.

Selection	2 5 10 20 50
Default value	20

Concentration

Concentration of the solution.

Range	-999999999 ... 999999999
Default value	1.000

Concentration unit

Units used to specify the concentration.

Selection	μmol/mL mmol/L mol/L g/L mg/L mg/mL μg/L ppm % mEq/L User-defined
Default value	mol/L

User-defined

A user-defined unit can be created. This will be added to the selection list. The previous entry will be overwritten as soon as the new unit has been defined. A blank entry can be generated this way as well.

Titer

Titer of the solution.

Range	-999999999 ... 999999999
Default value	1.000

Titer unit

Unit of the titer.

Selection	$\mu\text{mol/L}$ mmol/L mol/L g/L mg/L mg/mL $\mu\text{g/L}$ ppm $\%$ mEq/L empty User-defined
Default value	empty

User-defined

A user-defined unit can be created. This will be added to the selection list. The previous entry will be overwritten as soon as the new unit has been defined. A blank entry can be generated this way as well.

Date titer det.

Date of the last titer determination.

Monitoring

Switching the titer monitoring on and off.

Selection	on off
Default value	off

Time interval

This parameter is visible only when **Monitoring = on**.

You will be notified that this time interval (in days) has elapsed when starting a method. You can then select whether or not you would still like to start the method.

Range	1 ... 999 d
Default value	999 d

6.5 Managing common variables

6.5.1 General

Menu ▶ System ▶ Common variables

The instrument offers the possibility of saving five **method-independent variables**, so-called common variables. These variables remain saved in the instrument and can be used in future calculations. Common variables are useful, e.g. for the following applications:

- Determination of a blank value which will be taken into account during the content determination of the sample.
- Determination of the content of a standard solution, which will be taken into account during the content determination of the sample.



Common variables		ready
CV01	1.0472	
CV02	0.9638	
CV03	0.0	
CV04	0.0	
CV05	0.0	

Edit Delete

The common variables have the non-changeable designation **CV01...CV05**. The value is displayed for every variable. No unit can be assigned to the common variables.

Edit

Edit the data of the selected common variable, see following chapter.

Delete

Delete the elected common variable from the list.

6.5.2 Editing common variables

The common variables can be modified as follows:

- Manually in this dialog.
- Automatic assignment from the determination sequence. A calculation result must be configured accordingly for this purpose (see below).

Assigning a result automatically to a common variable

Proceed as follows:

1 Open the editing dialog of the result

- Select the menu item **Parameters ► Calculation** and press **[OK]**.
- Select the result whose value is to be assigned to a common variable.
- In the function bar, select **Edit** and press **[OK]**.

Edit result		ready
Result name	Blank	
R1=	EP1	
Decimal places	2	
Result unit	%	
Save as titer	off	
Save as CV	off	

2 Adjust the result properties

- Select the parameter **Save as CV** and press **[OK]**.
- Select the entry **on** in the selection list and apply with **[OK]**.

The assignment of the result to a common variable occurs automatically in accordance with the following scheme:

- Result **R1** ⇒ Common variable **CV01**
- Result **R2** ⇒ Common variable **CV02**
- etc.

6.6 File management

Menu ► System ► File management



Note

This menu item is visible only when a USB stick has been connected as an external storage medium.

Methods can be imported and deleted from a USB stick in this dialog. Only methods located in the **Files** directory are displayed in the list (see *"Directory structure on the USB stick"*, page 79).

A backup copy can be made of the system (all data and settings). Similarly, an existing backup copy can be stored.

Import

Import the selected method.

Delete

Delete the selected method.

Backup

Create a backup copy of all data and settings on the USB stick.



Note

Only **one** backup copy can be created on the same USB stick.

If a backup copy has already been stored on the stick, then it will be overwritten when this function is carried out again.

Restore

Load the backup copy from a connected USB stick.

Directory structure on the USB stick

A directory with the instrument number is generated on the USB stick. The structure within the directory appears as follows:



Figure 29 Directory structure on the USB stick

Backup	All of the files of the backup copy are stored in this directory. The directory will be created the first time a backup copy is saved.
Files	Exported methods will be stored in this directory. The directory will be created the first time a method is exported. Only methods being located in this directory can be imported.
pc_lims_report	PC/LIMS reports are stored in this directory as a TXT file. The directory will be created the first time a PC/LIMS report is printed out.

6.7 Configuring external devices

Menu ► System ► External devices

PC/LIMS report

Specification of the memory location for the PC/LIMS report. The PC/LIMS report is a machine-readable report with all of the data important for a determination. It can be saved as follows:

- as a TXT file on a USB stick.
- to a LIMS via an RS-232 interface. The 6.2148.030 RS-232/USB Box is required for this purpose.

Selection	COM2 USB Stick
Default value	USB Stick

COM2

The report is sent via the serial COM2 interface. The interface parameters set in the dialog **COM2 settings** are used (see "Editing the COM2 settings", page 82).

USB Stick

The report will be saved as a TXT file on the USB stick in the folder **pc_lims_report**.

Printer

If a printer is connected, then the printer type needs to be defined here in order for the reports to be printed out correctly.

The printers that have the designation **ESC-POS** are so-called POS printers (point-of-sale printers), i.e. they print on continuous paper.

Selection	Citizen (ESC-POS) Custom (ESC-POS) Epson Epson (ESC-POS) HP DeskJet HP LaserJet Seiko (ESC-POS)
Default value	HP DeskJet

Graphics width

Adjust the width of the curve to be printed out to the paper width of the printer to be used. The standard value depends on the selected printer. The height of the curve is 2/3 of the width.

Range	100 ... 3000 Pixels
-------	----------------------------

Keyboard layout

A commercially available USB keyboard can be connected to make it easier to enter text and numbers. Define the country-specific keyboard layout here.

Selection	English US French FR German CH German DE Spanish ES
Default value	English US

Balance

If a balance has been connected, then the balance type needs to be defined here.

Selection	AND Mettler Mettler AT Mettler AX Ohaus Precisa Sartorius Shimadzu
Default value	Sartorius

The following table indicates the balance type that needs to be selected for the balance model:

Balance	Type of balance
AND	AND
Mettler AB, AG, AM, PM, XP, XS	Mettler
Mettler AT	Mettler AT
Mettler AX, MX, UMX, PG, AB-S	Mettler AX
Ohaus Voyager, Explorer, Analytical Plus	Ohaus
Precisa	Precisa
Sartorius	Sartorius
Shimadzu BX, BW	Shimadzu



Editing the COM1 settings

Under **COM1 settings** the interface parameters for the connected balance are set.

Baud rate

Transfer rate in characters per second.

Selection	1200 2400 4800 9600 19200 38400 57600 115200
Default value	9600

Data bits

Number of data bits.

Selection	7 8
Default value	8

Stop bits

Number of stop bits.

Selection	1 2
Default value	1

Parity

Type of parity testing.

Selection	even none odd
Default value	none

Handshake

Type of data transfer protocol.

Selection	hardware software none
Default value	hardware



Note

In case of communication problems, try the software handshake (**software**).

Editing the COM2 settings

Under **COM2 settings** the interface parameters for devices connected to the **RS-232/2** connector of the RS-232 USB Box are set (e.g. PC). The parameters and input ranges are the same as for the COM1 interface.

6.8 Instrument diagnosis

6.8.1 Loading program versions and language files

Menu ► System ► Diagnosis

New program versions or language files can be loaded from a USB stick. The corresponding file must be saved on the USB stick in a directory with the instrument number (e.g. 848 or 863).

You can distinguish between language files and program files by noting how the file name is constructed.

Program files

They are instrument-specific. The file name has the following structure:

5XXXyyyy.bin whereas

XXX = Instrument type (e.g. 848 for the 848 Titrino plus)

yyyy = Program version

Language files

They can be recognized by means of the two-digit language code in the file name. A language file contains the dialog texts for various instrument types. It is not instrument-specific. The file name has the following structure:

5848xxxxYY.bin whereas

xxxx = Version number

YY = Language, e.g. DE (German), FR (French), ES (Spanish)

Loading a file

Proceed as follows:

1 Connect the USB stick

- Plug in the USB stick with the 6.2151.100 adapter (USB MINI (OTG) - USB A) at the USB port on the instrument.
- Switch on the instrument.

2 Open the update dialog

- Under **Menu ► System ► Diagnosis**, select the menu item **Software update**.
- Press **[OK]**.



```
Software update ready
Program version 58480011

Press [START] key to continue
```

3 Open the file selection

- Press [OK].

The selection list with the program and language files available on the USB stick is opened.

4 Select the file

- Select the required file with the arrow keys.
- Press [OK].

5 Start the update

- Press [START].

The update process is started, it runs automatically. At the end of the process, the instrument is automatically switched off and then back on again. No user intervention is required.

6.8.2 Diagnosis functions

Electronic and mechanical functional groups in Metrohm instruments can and should be checked as part of regular maintenance by specialist personnel from Metrohm. Please ask your local Metrohm agent regarding the precise terms and conditions involved in concluding a corresponding maintenance agreement.

7 Parameters

7.1 Dynamic equivalence point titrations (DET)

7.1.1 Start conditions

Menu ► Parameters ► Start conditions

The parameters that are carried out before the start of titration are defined under **Start conditions**.

Activation pulse

Output of an activation pulse on a remote line. This activation pulse starts a connected Dosimat.

Selection	on off
Default value	off

Start delay time

Waiting time after the start of the determination, before titration takes place. An auxiliary solution can, e.g. be added with a Dosimat during this time (parameterization at the Dosimat). The **Activation pulse** parameter must however be switched on for this purpose.

Range	0 ... 999999 s
Default value	0 s

Start volume

Volume that is dosed prior to the start of the titration.

Range	0.00000 ... 9999.99 mL
Default value	0.00000 mL

Dosing rate

Rate at which the start volume is dosed. The maximum dosing rate depends on the cylinder volume (*see Chapter 10.1.1, page 139*).

Range	0.01 ... 166.00 mL/min
Selection	max.
Default value	max.

Pause

Waiting time, e.g. in case of a slow response of the electrode after the start or reaction time after the dosing of a start volume.

Range	0 ... 999999 s
Default value	0 s



Request sample ID

Selection of the sample identification that is queried in the method run.

Selection	off ID1 ID2 ID1&ID2
Default value	off

Request sample size

If this parameter is switched on, then the value for the sample size will be queried.

Selection	on off
Default value	off

Request sample unit

If this parameter is switched on, then the unit for the sample size will be queried.

Selection	on off
Default value	off

Hold at request

If this parameter is switched on, then the method run will be paused during the query. If the parameter is switched off, the titration will be started in the background.

Selection	on off
Default value	on

7.1.2 Titration parameters

Menu ▶ Parameters ▶ Titration parameters

The parameters that influence the sequence of the entire titration are defined under **Titration parameters**.

Titration rate

Three predefined sets of parameters can be selected for the titration rate.

Selection	slow optimal fast user
Default value	optimal

slow

For titrations in which the finest details are also to be visible. This could however also lead to an increase in noise, which could result in unwanted equivalence points.

optimal

For all standard titrations. The parameters have been optimized for the most frequent applications.

fast

For fast and less critical titrations.

user

The individual titration parameters can be modified.

**Note**

Select the **optimal** titration rate when you are developing a new titration method. This parameter is suitable for virtually all titrations and only needs adaptation in special cases.

The settings of the individual titration rates are listed in *table 1, page 89*.

Meas. point density

This parameter is visible only when **Titration rate = user**.

A small value means small volume increments, i.e. a high measuring point density. The curve then shows all the finest details which also include noise; this could cause unwanted equivalence points to be found. A larger value, i.e. a smaller measuring point density, permits quicker titrations. If you are using a dosing device with a small cylinder volume then a smaller measuring point density value may be beneficial. However, you should also set a smaller signal drift and a higher EP criterion at the same time.

Range	0 ... 9
Default value	4

Min. increment

This parameter is visible only when **Titration rate = user**.

This smallest permitted volume increment is added at the start of the titration and with steep curves in the region of the equivalence point. Very small values should only be used if a low titrant consumption is expected; otherwise unwanted equivalence points could be evaluated.

Range	0.05 ... 999.90 µL
Default value	10.00 µL

Max. increment

This parameter is visible only when **Titration rate = user**.

A maximum volume increment should be selected in the following cases:

- when titration consumption is very low up until the equivalence point is reached.
- when a start volume is dosed up until shortly before the equivalence point is reached.



- when the change of direction in the jumping range is very abrupt, because otherwise it could easily happen that an excessively large volume will be dosed in the region of the equivalence point.

The value should not be less than 1/100 cylinder volume.

Range	0.1 ... 9999.9 μL
Selection	off
Default value	off



Note

It is not advisable to select similar volumes for the minimum and the maximum increment. Monotonic equivalence point titration (MET) is appropriate for these applications.

Dosing rate

This parameter is visible only when **Titration rate = user**.

Rate at which the volume increments are dosed. The maximum dosing rate depends on the cylinder volume (*see Chapter 10.1.1, page 139*).

Range	0.01 ... 166.00 mL/min
Selection	max.
Default value	max.

Signal drift

This parameter is visible only when **Titration rate = user**.

Maximum permissible drift for the measured value acceptance, i.e. maximum change of the measured value per minute. This type of titration is often referred to as equilibrium titration.



Note

A constant measured value is often only achieved after a certain time, as mixing and the reaction itself require a certain time. The response time of an electrode can also increase with time, i.e. reaching a constant measured value takes longer and longer. Drift-controlled measured value acceptance is particularly advisable in such cases, as the measured values are only accepted when equilibrium has almost been achieved.

Measuring mode pH, U and I_{pol}:

Range	0.1 ... 999.0 mV/min
Default value	50.0 mV/min

Selection **off**

off

Measured value acceptance will take place after the maximum waiting time has elapsed. This can be useful when the titration reaction proceeds slowly or the electrode is slow to respond.

Measuring mode Upol:

Range **0.01 ... 99.90 μ A/min**

Default value **50.00 μ A/min**

Selection **off**

off

Measured value acceptance will take place after the maximum waiting time has elapsed. This can be useful when the titration reaction proceeds slowly or the electrode is slow to respond.

Min. waiting time

This parameter is visible only when **Titration rate = user**.

The measured value is not accepted until the minimum waiting time has elapsed, even if the signal drift has already been reached. The minimum waiting time is only important for drift-controlled measurements.

Range **0 ... 999999 s**

Default value **0 s**

Max. waiting time

This parameter is visible only when **Titration rate = user**.

If signal drift has been switched off or has not yet been reached, then the measured value will be accepted when the maximum waiting time has elapsed.

Range **0 ... 999999 s**

Default value **26 s**

Table 1 Standard values for predefined titration rates for DET

	Titration rate		
	slow	optimal	fast
Meas. point density	2	4	6
Min. increment	10.00 μ L	10.00 μ L	30.00 μ L
Max. increment	off	off	off
Dosing rate	maximum	maximum	maximum



	Titration rate		
	slow	optimal	fast
Signal drift			
– pH, U and I _{pol}	20.0 mV/min	50.0 mV/min	80.0 mV/min
– U _{pol}	20.0 μA/min	50.0 μA/min	80.0 μA/min
Min. waiting time	0 s	0 s	0 s
Max. waiting time	38 s	26 s	21 s

Temperature

Manually entered titration temperature. If a temperature sensor is connected then the temperature will be measured continuously. This value is used for temperature correction in pH measurements.

Range	–20.0 ... 150.0 °C
Default value	25.0 °C

Sensor

Selection of the sensor from the sensor list. The selection depends on the measuring mode. Sensors are defined under **System ▶ Sensors**.

Selection	Selection of configured sensors
-----------	--

Solution

Selection of the solution from the solution list. We recommend always selecting the solution. This ensures that the correct data (titer, concentration, etc.) is always used for the calculation. Solutions are defined under **System ▶ Solutions**.

For buret units with integrated data chip, a check is made in the method run to verify whether the correct solution has been attached and whether the type of dosing drive matches. For buret units without integrated data chip, the cylinder volume and the type of dosing drive are checked. The validity of the titer is checked for the selected solution at the start of the determination.

Selection	Selection of configured solution not defined
Default value	not defined

not defined

No check takes place.

I(pol)

The polarization current is the current that is applied at a polarizable electrode during the voltametric measurement. This parameter is available only with I(pol) determinations.

Range	-125 ... 125 μA (Increment: 1)
Default value	5 μA

U(pol)

The polarization voltage is the voltage applied to the polarizable electrode during an amperometric measurement. This parameter is available only with U(pol) determinations.

Range	-1250 ... 1250 mV (Increment: 10)
Default value	400 mV

Electrode test

For polarizable electrodes, an electrode test can be carried out. A check is made that the electrode is properly connected and that no short-circuit is present. The electrode test is carried out when the determination is started. This parameter is available only with I(pol) and U(pol) determinations.

Selection	on off
Default value	off

Stirrer

When this parameter is switched on, then the stirrer is switched on at the start of the determination.

Selection	on off
Default value	on

Stirring rate

Setting the stirring rate. It can be set in steps of -15 to $+15$. The default setting **8** corresponds to 1000 U/min. The formula for calculating the rotational speed is specified in *chapter 10.2, page 140*. The optimum stirring rate can be tested in the manual control.

The algebraic sign of the stirring rate changes the direction in which the stirring is done. When the stirrer is viewed from above, this means:

- "+": counterclockwise rotation
- "-": clockwise rotation

Range	-15 ... 15
Default value	8

7.1.3 Stop conditions

Menu ▶ Parameters ▶ Stop conditions

The conditions for stopping a titration are defined under **Stop conditions**.



Stop volume

The titration is stopped when the specified volume has been dosed since the start of the titration. This volume should be adjusted to the size of the titration vessel in order to prevent the contents from running over.

Range	0.0000 ... 9999.99 mL
Default value	100.000 mL
Selection	off

Stop meas. value

The titration is stopped when the specified measured value has been reached since the start of the titration.

Measuring mode pH:

Range	-20.000 ... 20.000
Selection	off
Default value	off

Measuring mode U, I_{pol}:

Range	-1250.0 ... 1250.0 mV
Selection	off
Default value	off

Measuring mode U_{pol}:

Range	-125.0 ... 125.0 μA
Selection	off
Default value	off

Stop EP

The titration is stopped when the specified number of equivalence points has been found.

Range	1 ... 9
Default value	9
Selection	off

Volume after EP

This volume will be added when the number of equivalence points defined under **Stop EP** has been found. The curve shape after the equivalence point can also be seen this way.

Range	0.01000 ... 9999.99 mL
Selection	off
Default value	off

Stop time

The titration is stopped when the specified time has elapsed since the start of the titration.

Range	0 ... 999999 s
Selection	off
Default value	off

Filling rate

Rate at which the dosing cylinder is filled after the titration. The maximum filling rate depends on the cylinder volume (*see Chapter 10.1.1, page 139*).

Range	0.01 ... 166.00 mL/min
Selection	max.
Default value	max.

7.1.4 Evaluation**Menu ▶ Parameters ▶ Evaluation**

The parameters for the evaluation of the titration curve are defined under **Evaluation**.

Window

Switch this parameter on when equivalence points are to be recognized only in one particular measured value range (window) of the curve. Only one window can be defined.

Selection	on off
Default value	off

Lower limit

This parameter is visible only when **Window = on**.

Measured value for the lower limit.

Measuring mode pH:

Range	-20.000 ... 20.000
Default value	-20.000

Measuring mode U, I_{pol}:

Range	-1250.0 ... 1250.0 mV
Default value	-1250.0 mV

Measuring mode U_{pol}:

Range	-125.00 ... 125.00 μA
Default value	-125.00 μA



Upper limit

This parameter is visible only when **Window = on**.

Measured value for the upper limit.

Measuring mode pH:

Range	-20.000 ... 20.000
Default value	20.000

Measuring mode U, I_{pol}:

Range	-1250.0 ... 1250.0 mV
Default value	1250.0 mV

Measuring mode U_{pol}:

Range	-125.00 ... 125.00 µA
Default value	125.00 µA

EP criterion

The equivalence point criterion found (ERC = Equivalence point Recognition Criterion) is compared with this value. Equivalence points whose ERC is less than the value defined here will be ignored.

Range	0 ... 200
Default value	5

EP recognition

This parameter allows you to filter out only the equivalence points that are being sought.

for Window = off

Selection	all greatest last off
Default value	all

all

All equivalence points will be recognized.

greatest

Only the equivalence point with the greatest ERC value, i.e. the steepest jump, will be recognized.

last

Only the last equivalence point will be recognized.

off

No evaluation takes place.

for Window = on

Selection	first greatest last
Default value	first

first

Only the first equivalence point will be recognized.

greatest

Only the equivalence point with the greatest ERC value, i.e. the steepest jump, will be recognized.

last

Only the last equivalence point will be recognized.

Fixed EP1 at

The associated volume will be interpolated from the measuring point list for the measured value entered. The fixed endpoint must lie between the first and the final entry in the measuring point list.

Measuring mode pH:

Range	-20.000 ... 20.000
Selection	off
Default value	off

Measuring mode U, I_{pol}:

Range	-1250.0 ... 1250.0 mV
Selection	off
Default value	off

Measuring mode U_{pol}:

Range	-125.00 ... 125.00 µA
Selection	off
Default value	off

Fixed EP2 at

See **Fixed EP1 at**.

Evaluation and equivalence point criterion with DET

The equivalence points (EP) are localized in a way similar to the Tubbs method [1][2]. The volume value of the equivalence point (V_E) is shifted from the point of inflection (see arrow) towards the smaller circle of curvature for real asymmetric titration curves.

[1] C. F. Tubbs, *Anal. Chem.* **1954**, 26, 1670–1671.

[2] E. Bartholomé, E. Biekert, H. Hellmann, H. Ley, M. Weigert, E. Weise, *Ullmanns Encyklopädie der technischen Chemie*, Bd. 5, Verlag Chemie, Weinheim, 1980, S. 659.




Figure 30 Tubbs method for determining the equivalence point

The figure shows that the evaluation also requires measuring points from the measuring point list after the equivalence point as well.

For the recognition of the EPs found the set EP criterion is compared with the ERC (Equivalence point Recognition Criterion) found. The ERC is the first derivative of the titration curve combined with a mathematical function which is more sensitive for flat jumps than for steeper ones. EPs whose ERC is smaller than the defined EP criterion will not be recognized. The ERC is displayed in the results dialog for each discovered and recognized EP. If you adjust the EP criterion after the fact in order to recognize more or fewer EPs, then you can initiate the reevaluation in the results dialog with the **[Recalc]** key.

7.1.5 Calculation

7.1.5.1 General

Menu ▶ Parameters ▶ Calculation

A maximum of five calculations can be defined in one method. A series of variables (raw data from the determination, previously calculated results) is available for the calculations. Definition can be made for each calculation as to whether the result is to be saved as a titer or as a common variable.

Calculation	ready
R1: Content	
R2:	
R3:	
R4:	
R5:	
Edit Delete	

The result name is specified in the list for each calculation.

Edit

Edit the data of the selected calculation, see following chapter.

Delete

Delete the selected calculation.

7.1.5.2 Editing a calculation**Result name**

The result name is the text which will be shown in the result view and in the report.

Input	12 characters
Default value	empty

R1=...R5=

Shows the calculation formula. A special editor is opened for the definition (see Chapter 5.3, page 39).

Input	44 characters
Default value	empty

Decimal places

Number of decimal places used to display the result.

Range	0 ... 5
Default value	2

Result unit

The result unit is displayed and saved together with the result.

Selection	% mol/L mmol/L g/L mg/L mg/mL ppm g mg mL mg/piece °C µL mL/min User-defined
Default value	%

User-defined

A user-defined unit can be created. This will be added to the selection list. The previous entry will be overwritten as soon as the new unit has been defined. A blank entry can be generated this way as well.

Save as titer

The result can be saved as titer for this selected solution. If statistics has been switched on, then the current mean value of the determination series will be saved.

Selection	on off
Default value	off

Save as CV

The calculated result can be saved as a method-independent variable, a so-called common variable. The result is then also available in other meth-



ods for calculations. If statistics has been switched on, then the current mean value of the determination series will be saved.

Selection	on off
Default value	off

7.1.6 Statistics

Menu ▶ Parameters ▶ Statistics

The statistics calculation of a multiple determination is switched on under **Statistics** and it is defined how many determinations the series contains.

Statistics

If this function is switched on, then statistics calculations will be carried out for all of the defined results.

Selection	on off
Default value	off

Number of samples

The number of determinations that are carried out for the statistics calculations.

If an additional determination must be added to the determination series, because one determination was incorrect, for example, then this can be accomplished in the statistical overview (*see Chapter 5.9, page 58*).

Range	2 ... 20
Default value	3

7.1.7 Reports

Menu ▶ Parameters ▶ Reports

The reports that will be printed out automatically in connection with a determination are defined under **Reports**.

Results

The result report contains the calculated results, equivalence points and endpoints, sample data, etc.

Selection	on off
Default value	off

Curve

Curve report. The width of the curve is defined in the system settings (*see "Graphics width", page 81*).

Selection	on off
Default value	off

Calculations/Statistics

Output of the calculation formulas for the individual results. The results are printed with the maximum resolution. This makes recalculation with an external program possible. The following data is printed out in addition when statistics is switched on:

- Result and sample size of the individual determinations
- Mean value as well as absolute and relative standard deviation

Selection	on off
Default value	off

Measuring point list

Output of the measuring point list.

Selection	on off
Default value	off

Parameters

All of the parameters of the current method are printed out with the parameter report.

Selection	on off
Default value	off

PC/LIMS

The PC/LIMS report is a machine-readable report with all of the data important for a determination. The PC/LIMS report can be saved as a TXT file on a USB storage medium or sent via an RS-232 interface to a LIMS. The output location is defined in the system settings (*see "PC/LIMS report", page 80*).

The file name of the TXT file is constructed as follows: *PC_LIMS_Report-ID1-YYYYMMDD-hhmmss.txt*.

Selection	on off
Default value	off



7.2 Monotonic equivalence point titrations (MET)

7.2.1 Start conditions

Menu ▶ Parameters ▶ Start conditions

The parameters that are carried out before the start of titration are defined under **Start conditions**.

Activation pulse

Output of an activation pulse on a remote line. This activation pulse starts a connected Dosimat.

Selection	on off
Default value	off

Start delay time

Waiting time after the start of the determination, before titration takes place. An auxiliary solution can, e.g. be added with a Dosimat during this time (parameterization at the Dosimat). The **Activation pulse** parameter must however be switched on for this purpose.

Range	0 ... 999999 s
Default value	0 s

Start volume

Volume that is dosed prior to the start of the titration.

Range	0.00000 ... 9999.99 mL
Default value	0.00000 mL

Dosing rate

Rate at which the start volume is dosed. The maximum dosing rate depends on the cylinder volume (*see Chapter 10.1.1, page 139*).

Range	0.01 ... 166.00 mL/min
Selection	max.
Default value	max.

Pause

Waiting time, e.g. in case of a slow response of the electrode after the start or reaction time after the dosing of a start volume.

Range	0 ... 999999 s
Default value	0 s

Request sample ID

Selection of the sample identification that is queried in the method run.

Selection	off ID1 ID2 ID1&ID2
Default value	off

Request sample size

If this parameter is switched on, then the value for the sample size will be queried.

Selection	on off
Default value	off

Request sample unit

If this parameter is switched on, then the unit for the sample size will be queried.

Selection	on off
Default value	off

Hold at request

If this parameter is switched on, then the method run will be paused during the query. If the parameter is switched off, the titration will be started in the background.

Selection	on off
Default value	on

7.2.2 Titration parameters

Menu ► Parameters ► Titration parameters

The parameters that influence the sequence of the entire titration are defined under **Titration parameters**.

Titration rate

Three predefined sets of parameters can be selected for the titration rate.

Selection	slow optimal fast user
Default value	optimal

slow

For titrations in which the finest details are also to be visible. This could however also lead to an increase in noise, which could result in unwanted equivalence points.

optimal

For all standard titrations. The parameters have been optimized for the most frequent applications.

fast

For fast and less critical titrations.

user

The individual titration parameters can be modified.

**Note**

Select the **optimal** titration rate when you are developing a new titration method. This parameter is suitable for virtually all titrations and only needs adaptation in special cases.

The settings of the individual titration rates are listed in *table 2, page 104*.

Volume increment

This parameter is visible only when **Titration rate = user**.

Volume that is dosed in each dosing step. The choice of the correct volume increment is a basic requirement for achieving high accuracy. A good guideline is 1/20 of the expected EP volume. For steep jumps the volume increment should tend toward 1/100 and for flat jumps toward 1/10 of the EP volume.

Small volume increments are used for determining blank values or with very asymmetrical curves. The accuracy of the evaluation cannot be increased by using smaller increments as the measured value alterations between two measuring points are then of the same order of magnitude as the noise.

Range	0.00005 ... 999.900 mL
Default value	0.10000 mL

Dosing rate

This parameter is visible only when **Titration rate = user**.

Rate at which the volume increments are dosed. The maximum dosing rate depends on the cylinder volume (*see Chapter 10.1.1, page 139*).

Range	0.01 ... 166.00 mL/min
Selection	max.
Default value	max.

Signal drift

This parameter is visible only when **Titration rate = user**.

Maximum permissible drift for the measured value acceptance, i.e. maximum change of the measured value per minute. This type of titration is often referred to as equilibrium titration.



Note

A constant measured value is often only achieved after a certain time, as mixing and the reaction itself require a certain time. The response time of an electrode can also increase with time, i.e. reaching a constant measured value takes longer and longer. Drift-controlled measured value acceptance is particularly advisable in such cases, as the measured values are only accepted when equilibrium has almost been achieved.

Measuring mode pH, U and I_{pol}:

Range	0.1 ... 999.0 mV/min
Default value	50.0 mV/min
Selection	off

off

Measured value acceptance will take place after the maximum waiting time has elapsed. This can be useful when the titration reaction proceeds slowly or the electrode is slow to respond.

Measuring mode U_{pol}:

Range	0.01 ... 99.90 µA/min
Default value	50.00 µA/min
Selection	off

off

Measured value acceptance will take place after the maximum waiting time has elapsed. This can be useful when the titration reaction proceeds slowly or the electrode is slow to respond.

Min. waiting time

This parameter is visible only when **Titration rate = user**.

The measured value is not accepted until the minimum waiting time has elapsed, even if the signal drift has already been reached. The minimum waiting time is only important for drift-controlled measurements.

Range	0 ... 999999 s
Default value	0 s

Max. waiting time

This parameter is visible only when **Titration rate = user**.

If signal drift has been switched off or has not yet been reached, then the measured value will be accepted when the maximum waiting time has elapsed.



Range	0 ... 999999 s
Default value	26 s

Table 2 Standard values for predefined titration rates for MET

	Titration rate		
	slow	optimal	fast
Volume increment	0.05000 mL	0.10000 mL	0.20000 mL
Dosing rate	maximum	maximum	maximum
Signal drift			
– pH, U and I _{pol}	20.0 mV/min	50.0 mV/min	80.0 mV/min
– U _{pol}	20.0 µA/min	50.0 µA/min	80.0 µA/min
Min. waiting time	0 s	0 s	0 s
Max. waiting time	38 s	26 s	21 s

Temperature

Manually entered titration temperature. If a temperature sensor is connected then the temperature will be measured continuously. This value is used for temperature correction in pH measurements.

Range	–20.0 ... 150.0 °C
Default value	25.0 °C

Sensor

Selection of the sensor from the sensor list. The selection depends on the measuring mode. Sensors are defined under **System ► Sensors**.

Selection	Selection of configured sensors
-----------	--

Solution

Selection of the solution from the solution list. We recommend always selecting the solution. This ensures that the correct data (titer, concentration, etc.) is always used for the calculation. Solutions are defined under **System ► Solutions**.

For buret units with integrated data chip, a check is made in the method run to verify whether the correct solution has been attached and whether the type of dosing drive matches. For buret units without integrated data chip, the cylinder volume and the type of dosing drive are checked. The validity of the titer is checked for the selected solution at the start of the determination.

Selection	Selection of configured solution not defined
Default value	not defined

not defined

No check takes place.

I(pol)

The polarization current is the current that is applied at a polarizable electrode during the voltametric measurement. This parameter is available only with I(pol) determinations.

Range	-125 ... 125 μA (Increment: 1)
Default value	5 μA

U(pol)

The polarization voltage is the voltage applied to the polarizable electrode during an amperometric measurement. This parameter is available only with U(pol) determinations.

Range	-1250 ... 1250 mV (Increment: 10)
Default value	400 mV

Electrode test

For polarizable electrodes, an electrode test can be carried out. A check is made that the electrode is properly connected and that no short-circuit is present. The electrode test is carried out when the determination is started. This parameter is available only with I(pol) and U(pol) determinations.

Selection	on off
Default value	off

Stirrer

When this parameter is switched on, then the stirrer is switched on at the start of the determination.

Selection	on off
Default value	on

Stirring rate

Setting the stirring rate. It can be set in steps of -15 to $+15$. The default setting **8** corresponds to 1000 U/min. The formula for calculating the rotational speed is specified in *chapter 10.2, page 140*. The optimum stirring rate can be tested in the manual control.

The algebraic sign of the stirring rate changes the direction in which the stirring is done. When the stirrer is viewed from above, this means:

- "+": counterclockwise rotation
- "-": clockwise rotation

Range	-15 ... 15
Default value	8



7.2.3 Stop conditions

Menu ▶ Parameters ▶ Stop conditions

The conditions for stopping a titration are defined under **Stop conditions**.

Stop volume

The titration is stopped when the specified volume has been dosed since the start of the titration. This volume should be adjusted to the size of the titration vessel in order to prevent the contents from running over.

Range	0.00000 ... 9999.99 mL
Default value	100.000 mL
Selection	off

Stop meas. value

The titration is stopped when the specified measured value has been reached since the start of the titration.

Measuring mode pH:

Range	-20.000 ... 20.000
Selection	off
Default value	off

Measuring mode U, I_{pol}:

Range	-1250.0 ... 1250.0 mV
Selection	off
Default value	off

Measuring mode U_{pol}:

Range	-125.0 ... 125.0 μA
Selection	off
Default value	off

Stop EP

The titration is stopped when the specified number of equivalence points has been found.

Range	1 ... 9
Default value	9
Selection	off

Volume after EP

This volume will be added when the number of equivalence points defined under **Stop EP** has been found. The curve shape after the equivalence point can also be seen this way.

Range	0.01000 ... 9999.99 mL
Selection	off
Default value	off

Stop time

The titration is stopped when the specified time has elapsed since the start of the titration.

Range	0 ... 999999 s
Selection	off
Default value	off

Filling rate

Rate at which the dosing cylinder is filled after the titration. The maximum filling rate depends on the cylinder volume (*see Chapter 10.1.1, page 139*).

Range	0.01 ... 166.00 mL/min
Selection	max.
Default value	max.

7.2.4 Evaluation

Menu ► Parameters ► Evaluation

The parameters for the evaluation of the titration curve are defined under **Evaluation**.

Window

Switch this parameter on when equivalence points are to be recognized only in one particular measured value range (window) of the curve. Only one window can be defined.

Selection	on off
Default value	off

Lower limit

This parameter is visible only when **Window = on**.

Measured value for the lower limit.

Measuring mode pH:

Range	-20.000 ... 20.000
Default value	-20.000

Measuring mode U, I_{pol}:

Range	-1250.0 ... 1250.0 mV
Default value	-1250.0 mV

*Measuring mode Upol:*

Range	-125.00 ... 125.00 μA
Default value	-125.00 μA

Upper limit

This parameter is visible only when **Window = on**.

Measured value for the upper limit.

Measuring mode pH:

Range	-20.000 ... 20.000
Default value	20.000

Measuring mode U, Ipol:

Range	-1250.0 ... 1250.0 mV
Default value	1250.0 mV

Measuring mode Upol:

Range	-125.00 ... 125.00 μA
Default value	125.00 μA

EP criterion

The equivalence point criterion found (ERC = Equivalence point Recognition Criterion) is compared with this value. Equivalence points whose ERC is less than the value defined here will be ignored.

Measuring mode pH:

Range	0.10 ... 9.99
Default value	0.50

Measuring mode U, Ipol:

Range	1 ... 999 mV
Default value	30 mV

Measuring mode Upol:

Range	0.1 ... 99.9 μA
Default value	2.0 μA

EP recognition

This parameter allows you to filter out only the equivalence points that are being sought.

for Window = off

Selection	all greatest last off
Default value	all

all

All equivalence points will be recognized.

greatest

Only the equivalence point with the greatest ERC value, i.e. the steepest jump, will be recognized.

last

Only the last equivalence point will be recognized.

off

No evaluation takes place.

for Window = on

Selection	first greatest last
Default value	first

first

Only the first equivalence point will be recognized.

greatest

Only the equivalence point with the greatest ERC value, i.e. the steepest jump, will be recognized.

last

Only the last equivalence point will be recognized.

Fixed EP1 at

The associated volume will be interpolated from the measuring point list for the measured value entered. The fixed endpoint must lie between the first and the final entry in the measuring point list.

Measuring mode pH:

Range	-20.000 ... 20.000
Selection	off
Default value	off

Measuring mode U, Ipol:

Range	-1250.0 ... 1250.0 mV
Selection	off
Default value	off

Measuring mode Upol:

Range	-125.00 ... 125.00 µA
Selection	off
Default value	off

Fixed EP2 at

See **Fixed EP1 at**.

Edit

Edit the data of the selected calculation, see following chapter.

Delete

Delete the selected calculation.

7.2.5.2 Editing a calculation**Result name**

The result name is the text which will be shown in the result view and in the report.

Input	12 characters
Default value	empty

R1=...R5=

Shows the calculation formula. A special editor is opened for the definition (see Chapter 5.3, page 39).

Input	44 characters
Default value	empty

Decimal places

Number of decimal places used to display the result.

Range	0 ... 5
Default value	2

Result unit

The result unit is displayed and saved together with the result.

Selection	% mol/L mmol/L g/L mg/L mg/mL ppm g mg mL mg/piece °C µL mL/min User-defined
Default value	%

User-defined

A user-defined unit can be created. This will be added to the selection list. The previous entry will be overwritten as soon as the new unit has been defined. A blank entry can be generated this way as well.

Save as titer

The result can be saved as titer for this selected solution. If statistics has been switched on, then the current mean value of the determination series will be saved.

Selection	on off
Default value	off



Save as CV

The calculated result can be saved as a method-independent variable, a so-called common variable. The result is then also available in other methods for calculations. If statistics has been switched on, then the current mean value of the determination series will be saved.

Selection	on off
Default value	off

7.2.6 Statistics

Menu ▶ Parameters ▶ Statistics

The statistics calculation of a multiple determination is switched on under **Statistics** and it is defined how many determinations the series contains.

Statistics

If this function is switched on, then statistics calculations will be carried out for all of the defined results.

Selection	on off
Default value	off

Number of samples

The number of determinations that are carried out for the statistics calculations.

If an additional determination must be added to the determination series, because one determination was incorrect, for example, then this can be accomplished in the statistical overview (*see Chapter 5.9, page 58*).

Range	2 ... 20
Default value	3

7.2.7 Reports

Menu ▶ Parameters ▶ Reports

The reports that will be printed out automatically in connection with a determination are defined under **Reports**.

Results

The result report contains the calculated results, equivalence points and endpoints, sample data, etc.

Selection	on off
Default value	off

Curve

Curve report. The width of the curve is defined in the system settings (see "Graphics width", page 81).

Selection	on off
Default value	off

Calculations/Statistics

Output of the calculation formulas for the individual results. The results are printed with the maximum resolution. This makes recalculation with an external program possible. The following data is printed out in addition when statistics is switched on:

- Result and sample size of the individual determinations
- Mean value as well as absolute and relative standard deviation

Selection	on off
Default value	off

Measuring point list

Output of the measuring point list.

Selection	on off
Default value	off

Parameters

All of the parameters of the current method are printed out with the parameter report.

Selection	on off
Default value	off

PC/LIMS

The PC/LIMS report is a machine-readable report with all of the data important for a determination. The PC/LIMS report can be saved as a TXT file on a USB storage medium or sent via an RS-232 interface to a LIMS. The output location is defined in the system settings (see "PC/LIMS report", page 80).

The file name of the TXT file is constructed as follows: *PC_LIMS_Report-ID1-YYYYMMDD-hhmmss.txt*.

Selection	on off
Default value	off



7.3 Endpoint titrations (SET)

7.3.1 Conditioning

Menu ► Parameters ► Conditioning

The conditions required for conditioning are defined under **Conditioning**.

Conditioning

If this parameter is switched on, then the first time the titration is started the working medium will be titrated to the endpoint with the specified control parameters. The status is kept stable. The actual method run does not begin until **[START]** has been pressed once more. Conditioning will be carried out again automatically after the titration.

Selection	on off
Default value	off

Start drift

Conditioning OK will be displayed as soon as this volume drift has been reached and the titration can be started.

Range	1 ... 999 µL/min
Default value	20 µL/min

Drift correction

The endpoint volume can be corrected by taking the drift value into account. For this, the volume drift is multiplied with the drift correction time and this value is then subtracted from the endpoint volume. The drift correction time is the time interval between the end of conditioning and the end of the determination.

Selection	auto manual off
Default value	off

auto

The value of the current volume drift is applied automatically at the start of the titration.

manual

If the volume drift is known throughout a prolonged period of time, then this can be entered manually.

off

No drift correction takes place.

Drift value

This parameter is visible only when **Drift correction = manual**.

Volume drift for manual drift correction.

Range	0.0 ... 99.9 $\mu\text{L}/\text{min}$
Default value	0.0 $\mu\text{L}/\text{min}$

Cond. stop volume

Maximum permissible volume that can be dosed during conditioning. Conditioning is stopped when the specified volume is dosed. If conditioning is continued by pressing **[START]** once again, then the titrant volume that has already been dosed will not be taken into account; i.e. the dosing starts again at zero. The stop volume should be adjusted to the size of the titration cell in order to prevent any overflow.

Range	0.00000 ... 9999.99 mL
Default value	20.0000 mL
Selection	off

Cond. stop time

Maximum permissible time over which conditioning may take place. Conditioning is stopped when the specified time has elapsed.

Range	0 ... 999999 s
Selection	off
Default value	off

7.3.2 Start conditions

Menu ▶ Parameters ▶ Start conditions

The parameters that are carried out before the start of titration are defined under **Start conditions**.

Activation pulse

Output of an activation pulse on a remote line. This activation pulse starts a connected Dosimat.

Selection	on off
Default value	off

Start delay time

Waiting time after the start of the determination, before titration takes place. An auxiliary solution can, e.g. be added with a Dosimat during this time (parameterization at the Dosimat). The **Activation pulse** parameter must however be switched on for this purpose.

Range	0 ... 999999 s
Default value	0 s

Start volume

Volume that is dosed prior to the start of the titration.



Range	0.00000 ... 9999.99 mL
Default value	0.00000 mL

Dosing rate

Rate at which the start volume is dosed. The maximum dosing rate depends on the cylinder volume (*see Chapter 10.1.1, page 139*).

Range	0.01 ... 166.00 mL/min
Selection	max.
Default value	max.

Pause

Waiting time, e.g. in case of a slow response of the electrode after the start or reaction time after the dosing of a start volume.

Range	0 ... 999999 s
Default value	0 s

Request sample ID

Selection of the sample identification that is queried in the method run.

Selection	off ID1 ID2 ID1&ID2
Default value	off

Request sample size

If this parameter is switched on, then the value for the sample size will be queried.

Selection	on off
Default value	off

Request sample unit

If this parameter is switched on, then the unit for the sample size will be queried.

Selection	on off
Default value	off

Hold at request

If this parameter is switched on, then the method run will be paused during the query. If the parameter is switched off, the titration will be started in the background.

Selection	on off
Default value	on

7.3.3 Titration parameters

Menu ► Parameters ► Titration parameters

The parameters that influence the sequence of the entire titration are defined under **Titration parameters**.

Solution

Selection of the solution from the solution list. We recommend always selecting the solution. This ensures that the correct data (titer, concentration, etc.) is always used for the calculation. Solutions are defined under **System ► Solutions**.

For buret units with integrated data chip, a check is made in the method run to verify whether the correct solution has been attached and whether the type of dosing drive matches. For buret units without integrated data chip, the cylinder volume and the type of dosing drive are checked. The validity of the titer is checked for the selected solution at the start of the determination.

Selection	Selection of configured solution not defined
Default value	not defined

not defined

No check takes place.

Sensor

Selection of the sensor from the sensor list. The selection depends on the measuring mode. Sensors are defined under **System ► Sensors**.

Selection	Selection of configured sensors
-----------	--

I(pol)

The polarization current is the current that is applied at a polarizable electrode during the voltametric measurement. This parameter is available only with I(pol) determinations.

Range	-125 ... 125 μA (Increment: 1)
Default value	5 μA

U(pol)

The polarization voltage is the voltage applied to the polarizable electrode during an amperometric measurement. This parameter is available only with U(pol) determinations.

Range	-1250 ... 1250 mV (Increment: 10)
Default value	400 mV



Electrode test

For polarizable electrodes, an electrode test can be carried out. A check is made that the electrode is properly connected and that no short-circuit is present. The electrode test is carried out when the determination is started. This parameter is available only with I(pol) and U(pol) determinations.

Selection	on off
Default value	off

Stirrer

When this parameter is switched on, then the stirrer is switched on at the start of the determination.

Selection	on off
Default value	on

Stirring rate

Setting the stirring rate. It can be set in steps of -15 to $+15$. The default setting **8** corresponds to 1000 U/min. The formula for calculating the rotational speed is specified in *chapter 10.2, page 140*. The optimum stirring rate can be tested in the manual control.

The algebraic sign of the stirring rate changes the direction in which the stirring is done. When the stirrer is viewed from above, this means:

- "+": counterclockwise rotation
- "-": clockwise rotation

Range	$-15 \dots 15$
Default value	8

Temperature

Manually entered titration temperature. If a temperature sensor is connected then the temperature will be measured continuously. This value is used for temperature correction in pH measurements.

Range	$-20.0 \dots 150.0 \text{ } ^\circ\text{C}$
Default value	$25.0 \text{ } ^\circ\text{C}$

Titration direction

The titration direction is normally determined automatically from the start measured value and the set endpoint. It is recommended that, whenever possible, a positive or negative alteration of the measured value is entered. If two endpoints have been set then the titration direction will be defined automatically. In this case the setting will be ignored.

Selection	+ - auto
Default value	auto

+

Positive measured value alteration, i.e. in the direction of a higher pH value, greater voltage or greater current.

-

Negative measured value alteration, i.e. in the direction of a lower pH value, lesser voltage or lesser current.

auto

The titration direction is determined automatically from the start measured value and the set endpoint.

Extraction time

Minimum duration of the titration. Titration will not be stopped during the extraction time, even if the endpoint has already been reached. Titration is however stopped if a stop condition is fulfilled during this time (*see Chapter 7.3.6, page 122*). The entry of an extraction time may be advisable, e.g. for the titration of sparingly soluble samples.

Range	0 ... 999999 s
Default value	0 s

Time interval MP

Time interval for entering a measuring point in the measuring point list. The measuring point list is limited to 1000 measuring points.

Range	0.1 ... 999999.0 s
Default value	2.0 s

7.3.4 Control parameters EP1**Menu ▶ Parameters ▶ Control parameters EP1**

The control parameters for the first endpoint are defined under **Control parameters EP1**.

Endpoint 1 at

Measured value for the first endpoint.

Measuring mode pH:

Range	-20.000 ... 20.000
Selection	off
Default value	off

Measuring mode U and I_{pol}:

Range	-1250.0 ... 1250.0 mV
Selection	off
Default value	off

*Measuring mode Upol:*

Range	-125.00 ... 125.00 μA
Selection	off
Default value	off

Titration rate

Three predefined sets of parameters can be selected for the titration rate.

Selection	slow optimal fast user
Default value	optimal

slow

For steep titration curves for which dosing must be carried out in small steps at the endpoint.

optimal

For all standard titrations. The parameters have been optimized for the most frequent applications.

fast

For flat titration curves for which the endpoint is reached only slowly.

user

The individual titration parameters can be modified.

The settings of the individual titration rates are listed in *table 3, page 121*.

Dynamics

This parameter is visible only when **Titration rate = user**.

This parameter defines the control range before the specified endpoint. Individual volume steps are dosed in the control range, the dosing is finely controlled. The closer the endpoint, the slower the dosing until the dosing rate defined under **Min. rate** has been reached. The larger the control range, the slower the titration. Outside of the control range, dosing is carried out continuously, and the dosing rate is defined under **Max. rate**.

Measuring mode pH:

Range	0.001 ... 20.000
Default value	2.000
Selection	off

Measuring mode U and Ipol:

Range	0.1 ... 1250.0 mV
Default value	100.0 mV
Selection	off

Measuring mode Upol:

Range	0.01 ... 125.00 μA
Default value	10.00 μA
Selection	off

Max. rate

This parameter is visible only when **Titration rate = user**.

Rate at which dosing is carried out outside of the control range. The maximum dosing rate depends on the cylinder volume (see Chapter 10.1.1, page 139).

Range	0.01 ... 166.00 mL/min
Default value	10.00 mL/min
Selection	max.

Min. rate

This parameter is visible only when **Titration rate = user**.

Rate at which dosing is carried out at the very beginning of the titration and in the control range at the end of the titration. This parameter has a decisive influence on the titration rate and thus also on the accuracy. The smaller the selected minimum rate, the slower the titration.

Range	0.01 ... 9999.00 μL/min
Default value	25.00 μL/min

Table 3 Standard values for predefined titration rates for SET

	Titration rate		
	slow	optimal	fast
Dynamics			
– pH	5.000	2.000	0.500
– U und I _{pol}	300.0 mV	100.0 mV	30.0 mV
– U _{pol}	40.00 μ A	10.00 μ A	5.00 μ A
Max. rate	1.00 mL/min	10.00 mL/min	maximum
Min. rate	5.00 μ L/min	25.00 μ L/min	50.00 μ L/min

Stop criterion

The titration is stopped when the endpoint has been reached and this stop criterion has been fulfilled. If no stop criterion has been selected then the titration will not be stopped. The stop conditions (see Chapter 7.3.6, page 122) always lead to a stop, even if the stop criterion has not been reached.



Selection	drift time off
Default value	drift

drift

The titration is stopped when the stop drift has been reached.

time

The titration is stopped if the endpoint has been exceeded during a certain time period (**Delay time**).

off

The titration will not be stopped until the stop conditions have been fulfilled.

Stop drift

This parameter is visible only when **Stop criterion = drift**.

Titration is stopped when the endpoint and the stop drift have been reached.

Range	1 ... 999 µL/min
Default value	20 µL/min

Delay time

This parameter is visible only when **Stop criterion = time**.

When the endpoint is reached, the specified time is allowed to elapse after the last dosing and the titration is then stopped.

Range	0 ... 999 s
Default value	10 s

7.3.5 Control parameters EP2**Menu ▶ Parameters ▶ Control parameters EP2**

The control parameters for the second endpoint are defined under **Control parameters EP2**. The parameters and the input ranges are the same as for the first endpoint.

7.3.6 Stop conditions**Menu ▶ Parameters ▶ Stop conditions**

The conditions for stopping a titration are defined under **Stop conditions**, if this does not occur automatically. This could be the case when the set endpoint is not reached or if the stop criterion (*see "Stop criterion", page 121*) is not fulfilled.

Stop volume

The titration is stopped when the specified volume has been dosed since the start of the titration. This volume should be adjusted to the size of the titration vessel in order to prevent the contents from running over.

Range	0.00000 ... 9999.99 mL
Default value	100.000 mL
Selection	off

Stop time

The titration is stopped when the specified time has elapsed following the termination of the start conditions.

Range	0 ... 999999 s
Selection	off
Default value	off

Filling rate

Rate at which the dosing cylinder is filled after the titration. The maximum filling rate depends on the cylinder volume (*see Chapter 10.1.1, page 139*).

Range	0.01 ... 166.00 mL/min
Selection	max.
Default value	max.

7.3.7 Calculation

7.3.7.1 General

Menu ▶ Parameters ▶ Calculation

A maximum of five calculations can be defined in one method. A series of variables (raw data from the determination, previously calculated results) is available for the calculations. Definition can be made for each calculation as to whether the result is to be saved as a titer or as a common variable.

Calculation	readw
R1: Content	
R2:	
R3:	
R4:	
R5:	
Edit	Delete

The result name is specified in the list for each calculation.

Edit

Edit the data of the selected calculation, see following chapter.

Delete

Delete the selected calculation.



7.3.7.2 Editing a calculation

Result name

The result name is the text which will be shown in the result view and in the report.

Input	12 characters
Default value	empty

R1=...R5=

Shows the calculation formula. A special editor is opened for the definition (see Chapter 5.3, page 39).

Input	44 characters
Default value	empty

Decimal places

Number of decimal places used to display the result.

Range	0 ... 5
Default value	2

Result unit

The result unit is displayed and saved together with the result.

Selection	% mol/L mmol/L g/L mg/L mg/mL ppm g mg mL mg/piece °C µL mL/min User-defined
Default value	%

User-defined

A user-defined unit can be created. This will be added to the selection list. The previous entry will be overwritten as soon as the new unit has been defined. A blank entry can be generated this way as well.

Save as titer

The result can be saved as titer for this selected solution. If statistics has been switched on, then the current mean value of the determination series will be saved.

Selection	on off
Default value	off

Save as CV

The calculated result can be saved as a method-independent variable, a so-called common variable. The result is then also available in other methods for calculations. If statistics has been switched on, then the current mean value of the determination series will be saved.



Calculations/Statistics

Output of the calculation formulas for the individual results. The results are printed with the maximum resolution. This makes recalculation with an external program possible. The following data is printed out in addition when statistics is switched on:

- Result and sample size of the individual determinations
- Mean value as well as absolute and relative standard deviation

Selection	on off
Default value	off

Measuring point list

Output of the measuring point list.

Selection	on off
Default value	off

Parameters

All of the parameters of the current method are printed out with the parameter report.

Selection	on off
Default value	off

PC/LIMS

The PC/LIMS report is a machine-readable report with all of the data important for a determination. The PC/LIMS report can be saved as a TXT file on a USB storage medium or sent via an RS-232 interface to a LIMS. The output location is defined in the system settings (*see "PC/LIMS report", page 80*).

The file name of the TXT file is constructed as follows: *PC_LIMS_Report-ID1-YYYYMMDD-hhmmss.txt*.

Selection	on off
Default value	off

7.4 pH calibration (CAL)

7.4.1 Calibration parameters

Menu ► Parameters ► Calibration parameters

The parameters that influence the sequence of the entire calibration are defined under **Calibration parameters**.

Signal drift

Maximum permissible drift for the measured value acceptance, i.e. maximum change of the measured value per minute.

Range	0.1 ... 999.0 mV/min
Default value	2.0 mV/min
Selection	off

off

Measured value acceptance will take place after the maximum waiting time has elapsed. This can be useful when the electrode is slow to respond.

Min. waiting time

The measured value is not accepted until the minimum waiting time has elapsed, even if the signal drift has already been reached. The minimum waiting time is only important for drift-controlled measurements.

Range	0 ... 999999 s
Default value	10 s

Max. waiting time

If signal drift has been switched off or has not yet been reached, then the measured value will be accepted when the maximum waiting time has elapsed.

Range	0 ... 999999 s
Default value	110 s

Temperature

Manually entered calibration temperature. If a temperature sensor is connected then the temperature will be measured continuously.

Range	-20.0 ... 150.0 °C
Default value	25.0 °C

Sensor

Selection of the sensor from the sensor list. The selection depends on the measuring mode. Sensors are defined under **System ► Sensors**.



Selection	Selection of configured sensors
-----------	--

Stirrer

When this parameter is switched on, then the stirrer is switched on at the start of the determination.

Selection	on off
Default value	on

Stirring rate

Setting the stirring rate. It can be set in steps of -15 to $+15$. The default setting **8** corresponds to 1000 U/min. The formula for calculating the rotational speed is specified in *chapter 10.2, page 140*. The optimum stirring rate can be tested in the manual control.

The algebraic sign of the stirring rate changes the direction in which the stirring is done. When the stirrer is viewed from above, this means:

- "+": counterclockwise rotation
- "-": clockwise rotation

Range	$-15 \dots 15$
Default value	8

7.4.2 Buffers**Menu ▶ Parameters ▶ Buffers**

The buffer type and the number of buffers is defined under **Buffers**.

Buffer type

Selection of a predefined buffer series or definition of special buffers. In the case of predefined buffer series, the instrument automatically recognizes which buffer is involved.

Selection	Baker Beckman DIN Fisher Fluka Basel Hamilton Merck CertiPUR Merck Titrisol Metrohm Mettler NIST Precisa Radiometer Special
-----------	--

Merck CertiPUR

Reference temperature = 25 °C . The buffer type **Merck Titrisol** must be selected when using Merck CertiPUR buffers (20 °C).

Special

Up to five calibration buffers can be defined in the method. The automatic buffer detection is not activated in this case. The buffers must be measured precisely in the specified sequence.

Number of buffers

Number of buffers that are used for calibration. If calibration is accomplished with more than two buffers, then they can be used repeatedly give them more statistical weight. The first two buffers must however always be different from one another.

Selection	1 2 3 4 5
Default value	2

Buffer 1 pH

This parameter is visible only when **Buffer type = Special**.

Range	-20.000 ... 20.000
Default value	7.000

Buffer 2 pH

This parameter is visible only when **Buffer type = Special**.

Range	-20.000 ... 20.000
Default value	4.000
Selection	off

Buffer 3 pH

This parameter is visible only when **Buffer type = Special**.

Range	-20.000 ... 20.000
Selection	off
Default value	off

Buffer 4 pH

See **Buffer 3 pH**.

Buffer 5 pH

See **Buffer 3 pH**.

7.4.3 Reports

Menu ► Parameters ► Reports

The reports that will be printed out automatically in connection with a calibration are defined under **Reports**.

Results

The result report contains the specifications for the calibration (slope, pH(0), etc.).

Selection	on off
Default value	off



Parameters

All of the parameters of the current method are printed out with the parameter report.

Selection	on off
Default value	off

PC/LIMS

The PC/LIMS report is a machine-readable report with all of the data important for a determination. The PC/LIMS report can be saved as a TXT file on a USB storage medium or sent via an RS-232 interface to a LIMS. The output location is defined in the system settings (see "PC/LIMS report", page 80).

The file name of the TXT file is constructed as follows: *PC_LIMS_Report-ID1-YYYYMMDD-hhmmss.txt*.

Selection	on off
Default value	off

7.5 Automation: Dipping in special

Menu ▶ Parameters ▶ Automation

Automation

Display of the template used for the automation sequence.

Dripping time

Waiting time after the titration head moves out of the sample beaker and out of the rinsing beaker.

Range	0 ... 999 s
Default value	3 s

Rinsing time

Waiting time during which the electrode remains immersed in the rinsing beaker.

Range	0 ... 999 s
Default value	5 s

Stirring rate

Setting the stirring rate. It can be set in steps of -15 to $+15$. The default setting **8** corresponds to 1000 U/min. The formula for calculating the rotational speed is specified in *chapter 10.2, page 140*. The optimum stirring rate can be tested in the manual control.

The algebraic sign of the stirring rate changes the direction in which the stirring is done. When the stirrer is viewed from above, this means:

- "+": counterclockwise rotation
- "-": clockwise rotation

Range	-15 ... 15
Default value	8



Note

The setting of the stirring rate under **Menu ▶ Parameters ▶ Automation** applies only for stirring during the immersion of the electrode in the rinsing beaker. The stirring rate during the determination is set under **Menu ▶ Parameters ▶ Titration parameters**.

7.6 Automation: Dipping in special 2

Menu ▶ Parameters ▶ Automation

Automation

Display of the template used for the automation sequence.

Dripping time

Waiting time after the titration head moves out of the sample beaker and out of the rinsing beaker.

Range	0 ... 999 s
Default value	3 s

Rinsing time

Waiting time during which the electrode remains immersed in the rinsing beaker.

Range	0 ... 999 s
Default value	5 s

Stirring rate

Setting the stirring rate. It can be set in steps of -15 to $+15$. The default setting **8** corresponds to 1000 U/min. The formula for calculating the rotational speed is specified in *chapter 10.2, page 140*. The optimum stirring rate can be tested in the manual control.

The algebraic sign of the stirring rate changes the direction in which the stirring is done. When the stirrer is viewed from above, this means:

- "+": counterclockwise rotation



- "–": clockwise rotation

Range	–15 ... 15
Default value	8



Note

The setting of the stirring rate under **Menu ▶ Parameters ▶ Automation** applies only for stirring during the immersion of the electrode in the rinsing beaker. The stirring rate during the determination is set under **Menu ▶ Parameters ▶ Titration parameters**.

7.7 Automation: Double dipping

Menu ▶ Parameters ▶ Automation

Automation

Display of the template used for the automation sequence.

Dripping time

Waiting time after the titration head moves out of the sample beaker and out of the rinsing beaker.

Range	0 ... 999 s
Default value	3 s

Rinsing time

Waiting time during which the electrode remains immersed in the rinsing beaker.

Range	0 ... 999 s
Default value	5 s

Stirring rate

Setting the stirring rate. It can be set in steps of –15 to +15. The default setting **8** corresponds to 1000 U/min. The formula for calculating the rotational speed is specified in *chapter 10.2, page 140*. The optimum stirring rate can be tested in the manual control.

The algebraic sign of the stirring rate changes the direction in which the stirring is done. When the stirrer is viewed from above, this means:

- "+": counterclockwise rotation
- "–": clockwise rotation

Range	–15 ... 15
Default value	8



Note

The setting of the stirring rate under **Menu ▶ Parameters ▶ Automation** applies only for stirring during the immersion of the electrode in the rinsing beaker. The stirring rate during the determination is set under **Menu ▶ Parameters ▶ Titration parameters**.

7.8 Automation: Rinsing in sample

Menu ▶ Parameters ▶ Automation

Automation

Display of the template used for the automation sequence.

Dripping time

Waiting time after the titration head moves out of the sample beaker.

Range	0 ... 999 s
Default value	3 s

Aspiration time

Aspiration time of Pump 2, in case an 843 Pump Station is connected. It runs before the rinsing and after the rinsing.

Range	0 ... 999 s
Default value	10 s

Rinsing time

Rinsing time of Pump 1, if an 843 Pump Station is connected. Rinsing pump **and** aspiration pump run during the rinsing time.

Range	0 ... 999 s
Default value	5 s



7.9 Automation: Rinsing in special

Menu ▶ Parameters ▶ Automation

Automation

Display of the template used for the automation sequence.

Dripping time

Waiting time after the titration head moves out of the sample beaker and out of the rinsing beaker.

Range	0 ... 999 s
Default value	3 s

Aspiration time

Aspiration time of Pump 2, in case an 843 Pump Station is connected. It runs after the rinsing time.

Range	0 ... 999 s
Default value	10 s

Rinsing time

Rinsing time of Pump 1, in case an 843 Pump Station is connected. It runs before the aspiration time. Rinsing pump **and** aspiration pump run during the rinsing time.

At the end of a sample series, the rinsing time determines how long rinsing solution will be filled into the rinsing beaker.

Range	0 ... 999 s
Default value	5 s

8 Handling and maintenance

8.1 General

The 862 Compact Titrosampler requires appropriate care. Excess contamination of the instrument may result in malfunctions and a reduction in the service life of the sturdy mechanics and electronics of the instrument.

Severe contamination can also have an influence on the measured results. Regular cleaning of exposed parts can prevent this to a large extent.

Spilled chemicals and solvents must be removed immediately. Above all, the plug connections (in particular the mains connection socket) should be protected from contamination.

Check all tubing connections regularly for leakage.

8.2 Quality Management and validation with Metrohm

Quality Management

Metrohm offers you comprehensive support in implementing quality management measures for instruments and software. Further information on this can be found in the brochure «**Quality Management with Metrohm**» available from your local Metrohm agent.

Validation

Please contact your local Metrohm agent for support in validating instruments and software. Here you can also obtain validation documentation to provide help for carrying out the **Installation Qualification** (IQ) and the **Operational Qualification** (OQ). IQ and OQ are also offered as a service by the Metrohm agents. In addition, various application bulletins are also available on the subject, which also contain **Standard Operating Procedures** (SOP) for testing analytical measuring instruments for reproducibility and correctness.

Maintenance

Electronic and mechanical functional groups in Metrohm instruments can and should be checked as part of regular maintenance by specialist personnel from Metrohm. Please ask your local Metrohm agent regarding the precise terms and conditions involved in concluding a corresponding maintenance agreement.



Note

You can find information on the subjects of quality management, validation and maintenance as well as an overview of the documents currently available at www.metrohm.com/com/ under **Support**.

9 Troubleshooting

9.1 SET titration

Problem	Cause	Remedy
The titration will not be finished.	<i>The minimum dosing rate is too low.</i>	Define Titration rate = user and increase minimum rate (Min. rate) (see Chapter 7.3.4, page 119).
	<i>The stop criterion is unsuitable.</i>	Adjust the control parameters (see Chapter 7.3.4, page 119): <ul style="list-style-type: none"> ▪ Increase the stop drift. ▪ Select a short delay time.
The sample is over-titrated.	<i>The control parameters are unsuitable.</i>	Adjust the control parameters (see Chapter 7.3.4, page 119): <ul style="list-style-type: none"> ▪ Select Titration rate = slow. ▪ Define Titration rate = user and increase the control range. ▪ Define Titration rate = user and reduce maximum rate (Max. rate). ▪ Define Titration rate = user and reduce minimum rate (Min. rate). ▪ Stir faster. ▪ Arrange the electrode and buret tip to an optimum.
	<i>The electrode responds too slowly.</i>	Replace the electrode.
The titration time is too long.	<i>The control parameters are unsuitable.</i>	Adjust the control parameters (see Chapter 7.3.4, page 119): <ul style="list-style-type: none"> ▪ Select Titration rate = optimal or fast. ▪ Define Titration rate = user and reduce the control range. ▪ Define Titration rate = user and increase maximum rate (Max. rate). ▪ Define Titration rate = user and increase minimum rate (Min. rate).
	<i>The minimum dosing rate is too high.</i>	Define Titration rate = user and reduce minimum rate (Min. rate) (see Chapter 7.3.4, page 119).
The results are spread widely.	<i>The minimum dosing rate is too high.</i>	Define Titration rate = user and reduce minimum rate (Min. rate) (see Chapter 7.3.4, page 119).



Problem	Cause	Remedy
	<i>The electrode responds too slowly.</i>	Replace the electrode.

9.2 Miscellaneous

Problem	Cause	Remedy
No report is printed.	<i>The printer is not recognized by the instrument.</i>	<ul style="list-style-type: none"> ▪ Switch the 862 Compact Titrosampler off and on again. ▪ Use a hub on the USB connector and connect the printer to the hub.
	<i>The printer model is not compatible.</i>	Use a printer which fulfills the required specifications (<i>see Chapter 10.4.4, page 143</i>).
The USB keyboard or the PC mouse does not function.	<i>The keyboard or mouse is not recognized by the instrument.</i>	<ul style="list-style-type: none"> ▪ Switch the 862 Compact Titrosampler off and on again. ▪ Use a hub on the USB connector and connect the keyboard or mouse to the hub.
	<i>The keyboard or the mouse is not compatible.</i>	Use a model which fulfills the required specifications (<i>see Chapter 10.4, page 141</i>).
The display is no longer readable.	<i>The contrast is set incorrectly.</i>	Adjust the contrast correctly (<i>see Chapter 6.1, page 69</i>).
Mettler XP balances send "R" or "O ----" as ID1.	<i>The automatic calibration of the balance is switched on.</i>	Switch off the automatic calibration.
Message 020-511 "Action not possible" is displayed.	<i>The USB stick is no longer connected.</i>	<ol style="list-style-type: none"> 1. Connect the USB stick. 2. Switch the instrument off and on again.
	<i>The USB stick is full.</i>	<ul style="list-style-type: none"> ▪ Use a different USB stick. ▪ Delete files with the aid of a PC.

10 Appendix

10.1 Dosing unit

10.1.1 Maximum dosing and filling rate

The maximum dosing rate and maximum filling rate for the dosing unit depend on the cylinder volume:

Cylinder volume	Maximum rate
2 mL	6.67 mL/min
5 mL	16.67 mL/min
10 mL	33.33 mL/min
20 mL	66.67 mL/min
50 mL	166.00 mL/min

Independent of the cylinder volume, values ranging from 0.01 to 166.00 mL/min can always be entered. When the function is carried out the rate will be, if necessary, decreased automatically to the highest possible value.

10.1.2 Parameters for preparing (PREP) and emptying (EMPTY)

The **PREP** function (Preparing) is used to rinse the cylinder and tubings of the dosing unit and fill them air bubble-free. You should carry out this function before the first determination or once per day.

The **EMPTY** function empties the cylinder and the tubings of the dosing unit.

With these functions you can easily change the reagent in the dosing unit without coming into contact with chemicals.

Preparing and emptying are carried out with the following, non-alterable settings:

- The entire cylinder volume is dosed via port 1.
- For filling/emptying the tubings the following configuration is assumed:
 - Tubing on port 1: length = 70 cm, diameter = 2 mm
 - Tubing on port 2: length = 25 cm, diameter = 2 mm
- Dosing and filling are both carried out with the same maximum rate.



10.2 Stirring rate

The stirring rate can be adjusted in steps of -15 to $+15$.

The approximate rotational speed can be calculated with the following formula:

$$\text{Rotational speed/min (r/min)} = 125 \cdot \text{stirring rate}$$

Example:

Stirring rate set: 8

Rotational speed in rpm = $125 \cdot 8 = 1000$




Figure 31 Rotational speed depending on stirring rate

10.3 Balance

The sample size and the associated unit can be sent from a connected balance. The sample size is transmitted as a number with up to ten characters (including algebraic sign and decimal point).

Sample size and unit are sent as a single character string. They are separated by a space character. The string is terminated with the ASCII characters **CR** and **LF**.

If the balance sends a negative sample size (e.g. when reweighing), then the algebraic sign is adopted. The algebraic sign is however ignored for the calculations.



Note

With some balances, the sample identification and the method can be sent in addition to the sample size.

Make sure that the balance does not send the sample size until the end.

Mettler AX

For the Mettler AX balance, the fields that contain the sample identification or the method must be designated as follows:

- Designation for the field with the method name: **METHOD**
- Designation for the field with sample identification 1: **ID1**
- Designation for the field with sample identification 2: **ID2**

10.4 USB devices



Note

USB peripheral devices that are to be connected directly must support the USB 1.0/1.1 (Full Speed) or USB 2.0 (High Speed) standard. The maximum data transfer rate is however in any case 12 MBit/s.

Keyboards, PC mice and barcode readers are so-called HID devices (**H**uman **I**nterface **D**evice) and can be connected via a USB hub **only**.

Printers should also be connected via a USB hub. Depending on the manufacturer or printer type a direct connection is however possible.

10.4.1 Numerical USB keypad 6.2147.000

The key **[Num Lock]** must be pressed for navigating in the dialog. Then the arrow keys are effective.

The respective editing dialog must be opened for the number input.

Table 4 Keyboard assignment

Key of the 862 Compact Titro-sampler or function in the editing dialog	Key on the numerical USB keypad
[BACK]	[Home]
[↑] [↓]	[↑] [↓]
[←] [→]	[←] [→]
[OK]	[Enter]
[+-]	[BS] (backspace)
Delete	[Del]
Accept	[Home]



10.4.2 Key assignment of a USB keyboard

A commercially available USB keyboard can be connected to make it easier to enter text and numbers.

The respective editing dialog must be opened for text and number inputs.

Table 5 Keyboard assignment

Key of the 862 Compact Titrosampler or function in the editing dialog	Key on the USB keyboard
[BACK]	[Esc]
[↑] [↕]	[↑] [↓]
[⇐] [⇒]	[←] [→]
[OK]	[↵] (enter key) or [Enter] on the numerical keypad
[STOP]	[Ctrl] + [S]
[START]	[Ctrl] + [G]
[+-]	[←] (backspace)
Delete	[Delete]
Cancel	[Strg] + [Q]
Accept	[Esc]



Note

The lettering of the USB keyboard can differ from above lettering, depending on the country-specific keyboard used.

10.4.3 PC mouse

In order to make navigating in the dialog of the 862 Compact Titrosampler easier, a PC mouse can be connected.

Table 6 Mouse functions

Key of the 862 Compact Titrosampler	Mouse function
[OK]	Left mouse button

Key of the 862 Compact Titro-sampler	Mouse function
[BACK]	Right mouse button
[↑] [↓] [←] [→]	Mouse movement vertical/horizontal
[↑] [↓]	Scroll wheel vertical

10.4.4 Printer

The range of USB printers available is extremely varied and constantly changing. The following points must be taken into account when selecting a printer:

- USB interface necessary
- Printer language: HP-PCL, Canon BJI Commands, Epson ESC P/2 or ESC/POS



Note

Inexpensive printers are often designed solely for use with a PC and may not be equipped with one of the printer languages listed above. Such models are not suitable for this reason.

10.5 System initialization

In very rare instances, it could happen that a faulty file system (e.g. because of a program crash) will lead to an impairment of program functioning. The internal file system needs to be initialized in such cases.



Caution

All user data (methods, solutions, etc.) are deleted if a system initialization is carried out. Afterwards, the instrument will have the factory settings again.

We recommend creating a backup copy of the system at regular intervals in order to avoid data losses.

After a system initialization the program versions and language files do not have to be reloaded. Only the selection of the dialog language may have to be reset in the system settings.



Proceed as follows for system initialization:

1 Switch off the instrument

- Keep the red **[STOP]** key pressed down for at least 3 s.

A progress bar is displayed. If the key is released during this time, then the instrument will not be switched off.

2 Switch on the instrument

- Keep the red **[STOP]** key pressed down for approx. 10 s.

The dialog for confirmation of the initialization is displayed for 8 s. The initialization must be confirmed during this time.

```
System reset request detected.  
>> Press [BACK] key twice  
to confirm !  
>> Time remaining: 8 sec
```

3 Confirm the initialization



Note

If the query is not confirmed within 8 s, then the procedure will be interrupted.

- Press **[BACK]** twice.

Initialization is started. The process takes approximately 80 s. The instrument will be automatically restarted after successful initialization.

10.6 Remote interface

10.6.1 Pin assignment of the remote interface





Figure 32 Pin assignment of remote socket and plug

The above figure of the pin assignment applies for all Metrohm instruments with 9-pin D-Sub remote connectors.

Table 7 Inputs and outputs of the remote interface

Pin No.	Assignment	Function
1	Output 0	Sample ready
2	Output 1	Dosimat
3	Output 2	Pump 1
4	Output 3	Pump 2
5	Output 4	Error
6	0 Volt (GND)	
7	+5 volts	
8	Input 0	Start
9	Input 1	Stop

Inputs




approx. 5 k Ω Pull-up

$t_p > 100$ ms

active = low, inactive = high



Outputs



Open Collector

$$t_p \text{ } \overline{t_p} > 200 \text{ ms}$$

active = low, inactive = high

$$I_C = 20 \text{ mA}, V_{CE0} = 40 \text{ V}$$

+5 V: maximum load = 20 mA

10.6.2 Status diagram of the remote interface




Figure 33 Remote status diagram

11 Technical specifications

11.1 Measuring inputs

11.1.1 Potentiometry

A high-ohm measuring input (**Ind.**) for pH and redox electrodes and a measuring input for separate reference electrodes (**Ref.**).

Input resistance > $1 \cdot 10^{12} \Omega$

Offset current < $1 \cdot 10^{-12} \text{ A}$

Measuring mode
pH

Measuring range -13...20

Resolution 0.001

Measuring accuracy ± 0.003
(± 1 digit, without sensor error, under reference conditions)

Measuring mode
U

Measuring range -1200...1200 mV

Resolution 0.1 mV

Measuring accuracy $\pm 0.2 \text{ mV}$
(± 1 digit, without sensor error, under reference conditions)

11.1.2 Polarizer

A measuring input (**Pol.**) for polarizable electrodes.

Measuring mode
I_{pol} Determination with adjustable polarization current.

Polarization current -120...120 μA (increment: 1 μA)
-125...-121 μA / +121...+125 μA : non-guaranteed values, dependent on reference voltage +2.5 V

Measuring range -1200...1200 mV

Resolution 0.1 mV

Measuring accuracy $\pm 0.2 \text{ mV}$
(± 1 digit, without sensor error, under reference conditions)



<i>Measuring mode</i>	Determination with adjustable polarization voltage.
<i>Upol</i>	
<i>Polarization voltage</i>	–1200...1200 mV (increment: 10 mV) –1250...–1210 mV / +1210...+1250 mV: non-guaranteed values, dependent on reference voltage +2.5 V
<i>Measuring range</i>	–120...120 μ A
<i>Resolution</i>	0.01 μ A
<i>Measuring accuracy</i>	–

11.1.3 Temperature

A measuring input (**Temp.**) for temperature sensors of the Pt1000 or NTC type with automatic temperature compensation.

R (25 °C) and B value can be configured for NTC sensors.

<i>Measuring range</i>	
<i>Pt1000</i>	–150...250 °C
<i>NTC</i>	–5...250 °C (For an NTC sensor with R (25 °C) = 30000 Ω and B (25/50) = 4100 K)
<i>Resolution</i>	
<i>Pt1000</i>	0.1 °C
<i>NTC</i>	0.1 °C
<i>Measuring accuracy</i>	
<i>Pt1000</i>	± 0.2 °C (applies for measuring range –20...150 °C)
<i>NTC</i>	± 0.6 °C (applies for measuring range –10...40 °C)

11.2 Dosing drive

<i>Resolution</i>	10000 steps per cylinder volume
<i>Dosing unit</i>	
<i>Cylinder volume</i>	<ul style="list-style-type: none"> ▪ 2 mL ▪ 5 mL ▪ 10 mL 10 mL ▪ 20 mL ▪ 50 mL
<i>Accuracy</i>	Fulfills ISO/DIN standard 8655-3

11.3 Lift

<i>Stroke path</i>	132 mm
<i>Maximum load</i>	5 N
<i>Lift rate</i>	15 mm/s (typical)

11.4 Turntable

<i>Rack positions</i>	12
<i>Maximum load</i>	17 N
<i>Shift rate</i>	13 degrees/s (typical)

11.5 Interfaces and connectors

<i>Stirrer connector</i>	DIN socket
<i>Stirring rate</i>	Rod Stirrer 722/802: 180...3000 rpm Adjustable in 15 steps each in both directions of rotation.
<i>USB (OTG) connector</i>	For connecting USB devices.
<i>Remote connector</i>	For connecting instruments with a remote interface.



11.6 Mains connection

<i>Voltage</i>	100...240 V
<i>Frequency</i>	50...60 Hz
<i>Power consumption</i>	45 W
<i>Fuse</i>	2.0 ATH

11.7 Safety specifications

<i>Design and testing</i>	According to EN/IEC/UL 61010-1, CSA-C22.2 No. 61010-1, EN/IEC 61010-2-081, protection class I
<i>Safety instructions</i>	This document contains safety instructions which have to be followed by the user in order to ensure safe operation of the instrument.

11.8 Electromagnetic compatibility (EMC)

<i>Emission</i>	Standards fulfilled: <ul style="list-style-type: none">▪ EN/IEC 61326-1▪ EN/IEC 61000-6-3▪ EN 55022 / CISPR 22
<i>Immunity</i>	Standards fulfilled: <ul style="list-style-type: none">▪ EN/IEC 61326-1▪ EN/IEC 61000-6-1▪ EN/IEC 61000-4-2▪ EN/IEC 61000-4-3▪ EN/IEC 61000-4-4▪ EN/IEC 61000-4-5▪ EN/IEC 61000-4-6▪ EN/IEC 61000-4-11▪ EN/IEC 61000-4-14▪ NAMUR

11.9 Ambient temperature

<i>Nominal function range</i>	5...45 °C Humidity < 85 %
<i>Storage</i>	-20...60 °C
<i>Transport</i>	-40...60 °C

11.10 Reference conditions

<i>Ambient temperature</i>	25 °C (±3 °C)
<i>Relative humidity</i>	≤ 60 %

11.11 Dimensions

<i>Width</i>	0.26 m
<i>Height</i>	0.47 m
<i>Depth</i>	0.43 m
<i>Weight</i>	9.09 kg (without accessories)
<i>Material</i>	
<i>Housing</i>	Lower part: Crastin PBTP Lift: Metal, surface-treated
<i>Rack</i>	PVC



12 Conformity and warranty

12.1 Declaration of Conformity

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

Name of commodity

862 Compact Titrosampler

Combined sample changer and titrator for automated processing of small sample series in analytical laboratories.

This instrument has been built and has undergone final type testing according to the standards:

Electromagnetic compatibility

Emission: EN/IEC 61326-1: 2006, EN/IEC 61000-6-3: 2004, EN 55022 / CISPR 22: 2006

Immunity: EN/IEC 61326-1: 2006, EN/IEC 61000-6-1: 2007, EN/IEC 61000-4-2: 2001, EN/IEC 61000-4-3: 2002, EN/IEC 61000-4-4: 2004, EN/IEC 61000-4-5: 2001, EN/IEC 61000-4-6: 2001, EN/IEC 61000-4-11: 2004, EN/IEC 61000-4-14: 2004, NAMUR: 2004

Safety specifications

EN/IEC 61010-1: 2001, UL 61010-1: 2004, CSA-C22.2 No. 61010-1: 2004, EN/IEC 61010-2-081: 2003, protection class I



This instrument meets the requirements of the CE mark as contained in the EU directives 2006/95/EC (LVD), 2004/108/EC (EMC). It fulfils the following specifications:

EN 61326-1	Electrical equipment for measurement, control and laboratory use – EMC requirements
EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use
EN 61010-2-081	Particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes



Manufacturer

This instrument meets the requirements of the ETL Listed Mark for the North American market. It conforms to the electrical safety standards UL 61010-1 and CSA-C22.2 No. 61010-1. This product is listed in Intertek's Directory of Listed Products.

Metrohm Ltd., CH-9101 Herisau/Switzerland

Metrohm Ltd. is holder of the SQS certificate ISO 9001:2000 Quality management system for development, production and sales of instruments and accessories for ion analysis.

Herisau, 28 September 2009

D. Strohm

Vice President, Head of R & D

A. Dellenbach

Head of Quality Management

12.2 Warranty (guarantee)

Metrohm guarantees that the deliveries and services it provides are free from material, design or manufacturing errors. The warranty period is 36 months from the day of delivery; for day and night operation it is 18 months. The warranty remains valid on condition that the service is provided by an authorized Metrohm service organization.

Glass breakage is excluded from the warranty for electrodes and other glassware. The warranty for the accuracy corresponds to the technical specifications given in this manual. For components from third parties that make up a considerable part of our instrument, the manufacturer's warranty provisions apply. Warranty claims cannot be pursued if the Customer has not complied with the obligations to make payment on time.

During the warranty period Metrohm undertakes, at its own choice, to either repair at its own premises, free of charge, any instruments that can be shown to be faulty or to replace them. Transport costs are to the Customer's account.

Faults arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc. are expressly excluded from the warranty.



12.3 Quality Management Principles

Metrohm Ltd. holds the ISO 9001:2000 Certificate, registration number 10872-02, issued by SQS (Swiss Association for Quality and Management Systems). Internal and external audits are carried out periodically to assure that the standards defined by Metrohm's QM Manual are maintained.

The steps involved in the design, manufacture and servicing of instruments are fully documented and the resulting reports are archived for ten years. The development of software for PCs and instruments is also duly documented and the documents and source codes are archived. Both remain the possession of Metrohm. A non-disclosure agreement may be asked to be provided by those requiring access to them.

The implementation of the ISO 9001:2000 quality management system is described in Metrohm's QM Manual, which comprises detailed instructions on the following fields of activity:

Instrument development

The organization of the instrument design, its planning and the intermediate controls are fully documented and traceable. Laboratory testing accompanies all phases of instrument development.

Software development

Software development occurs in terms of the software life cycle. Tests are performed to detect programming errors and to assess the program's functionality in a laboratory environment.

Components

All components used in the Metrohm instruments have to satisfy the quality standards that are defined and implemented for our products. Suppliers of components are audited by Metrohm as the need arises.

Manufacture

The measures put into practice in the production of our instruments guarantee a constant quality standard. Production planning and manufacturing procedures, maintenance of production means and testing of components, intermediate and finished products are prescribed.

Customer support and service

Customer support involves all phases of instrument acquisition and use by the customer, i.e. consulting to define the adequate equipment for the analytical problem at hand, delivery of the equipment, user manuals, training, after-sales service and processing of customer complaints. The Metrohm service organization is equipped to support customers in implementing standards such as GLP, GMP, ISO 900X, in performing Opera-



tional Qualification and Performance Verification of the system components or in carrying out the System Validation for the quantitative determination of a substance in a given matrix.





13 Accessories







Note

Subject to change without notice.



13.1 Scope of delivery 2.862.0010

Qty.	Order no.	Description
1	1.802.0010	802 Stirrer Rod stirrer with 6.1909.010 Propeller stirrer.
		
1	1.862.0010	862 Compact Titrator Compact and economic instrument for the automatic determination of dynamic (DET) and monotonic (MET) titrations with automatic equivalence point detection as well as SET endpoint titration.
		

Qty.	Order no.	Description
2	6.1446.000	SGJ stopper / B-14/(15)
	Material:	PP
	Height (mm):	30.5
	SGJ size:	B-14/(15)
		
1	6.1459.300	Sample beaker / 120 mL / 100 pieces
	Sample beaker with screw cap for sampling, 100 pieces.	
	Material:	PP
	Height (mm):	113
	Outer diameter (mm):	40
	Volume (mL):	120
		
1	6.1815.010	Spiral band 0.5 m
	For holding together different cables or tubing.	
	Length (m):	0.5
		
1	6.1909.050	Propeller for 120 mL beaker (6.1459.300)
	Propeller for 802 Stirrer when this is used with sample racks for 120 mL beakers (6.1459.300)	
	Material:	ETFE
	Length (mm):	111
		







Qty.	Order no.	Description	
1	6.2103.130	Adapter red 2 mm plug / 4 mm socket For connecting plug B (4 mm) to 2 mm socket.	
1	6.2103.140	Adapter black 2 mm plug / B socket 4 mm For connecting plug B (4 mm) to 2 mm socket.	
1	6.2147.000	Numerical USB keypad Numerical USB keypad for Titrino plus, Dosimat plus, 862 Compact Titrosampler and 869 Compact Sample Changer. Incl. USB Hub with two connectors.	
1	6.2151.100	Adapter USB MINI (OTG) - USB A For connecting USB instruments.	

Qty.	Order no.	Description
1	6.2621.130	Hexagon key 2 mm 2 mm.
		
1	6.2621.140	Hexagon key 2.5 mm
		
1	A.725.0003	Metrohm USB Stick, 256 MB
1	6.2122.0x0	Mains cable with C13 line socket IEC-60320-C13
		Cable plug according to customer requirements.
		Switzerland: Type SEV 12 6.2122.020
		Germany, ...: Type CEE(7), VII 6.2122.040
		USA, ...: Type NEMA/ASA 6.2122.070
1	8.862.8002DE	862 Compact Titrosampler Manual







13.2 Optional accessories

Order no.	Description	
2.141.0100	USB Thermal printer NEO's Compact printer with USB interface, paper width 60 mm (40 characters). Including 6.2151.120 USB cable.	
2.800.0010	800 Dosino Drive with write/read hardware for intelligent dosing units. With fixed cable.	
2.843.0030	843 Pump Station (membrane) - rinse/aspirate The 843 Pump Station (membrane) has two built-in membrane pumps. These can be controlled directly from the 869 Compact Sample Changer via remote signals. The Rinse / Aspirate version is provided with all the accessories needed for automatically emptying the titration beaker and rinsing the titration equipment.	
2.843.0130	843 Pump Station (peristaltic) - rinse/aspirate The 843 Pump Station (peristaltic) has two built-in peristaltic pumps. These can be controlled directly from the 869 Compact Sample Changer via remote signals. The Rinse / Aspirate version is provided with all the accessories needed for automatically emptying the titration beaker and rinsing the titration equipment.	

Order no.	Description	
2.865.0010	865 Dosimat plus	
<p>Dispensing Unit for manual titration and dosing applications. Including manual push-button cable for manual dispensing control and 6.3026.220 Exchange Unit 20 mL.</p>		
2.876.0010	876 Dosimat plus	
<p>876 Dosimat plus for manual titration and dosing applications. Including push-button cable for manual dispensing control and 6.3026.220 Exchange Unit 20 mL.</p>		
6.2104.020	Electrode cable / 1 m / F	
<p>For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).</p>		
Length (m):	1	
6.2136.010	Remote cable (9p) - Dosimat (2xSt B)	
<p>For 870 Titrino Plus and 862 Compact Titrosampler</p>		
6.2141.230	Remote cable Compact Sample Changer/Compact Autosampler - Titrino plus/Pump Station	
<p>Cable for the remote connection between Compact Sample Changer/Compact Autosampler and Titrino plus and/or Pump Station.</p>		



Order no.	Description	
6.2141.240	Remote cable Dosimat plus - Compact Sample Changer/Titrino plus	
	Cable for the remote connection between Dosimat plus and Compact Sample Changer and/or Titrino plus using 6.2141.260 Cable.	
	6.3032.150	Dosing Unit 5 mL
	Dosing unit with integrated data chip with 5 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, antidiffusion buret tip.	Volume (mL): 5
6.3032.210	Dosing Unit 10 mL	
Dosing unit with integrated data chip with 10 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, antidiffusion buret tip.	Volume (mL): 10	
6.3032.220	Dosing Unit 20 mL	
Dosing unit with integrated data chip with 20 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, antidiffusion buret tip.	Volume (mL): 20	
6.3032.250	Dosing unit 50 mL	
Dosing unit with integrated data chip with 50 mL glass cylinder and light protection, mountable on reagent bottle with ISO/DIN GL45 thread. FEP tubing connection, antidiffusion buret tip.	Volume (mL): 50	

Index

Numbers/Symbols

802 Stirrer 17

A

Adapter

Connecting 18

Aspiration time 133, 134

Aspiration tip 15

Aspiration tubing 15, 16

Automation 29, 130, 131, 132,

133, 134

Auxiliary solution 24, 25

Auxiliary solutions 34

B

Backup 79

Balance 47, 140

Connect 20

Balance type 81

Buffers

CAL 128

Buret unit

Preparing (PREP) 62

C

Cable cover

Mounting 26

Removing 11

CAL

Parameters 127

Calculation 96, 110, 123

Delete 97, 111, 123

Edit 96, 111, 123

Calibration parameters

CAL 127

Canister 16

Common variable 77

Delete 78

Edit 78

Conditioning

SET 114

Contrast 71

Control parameters

SET 119

Curve 57

D

DET 2, 28

Parameters 85

Determination

Result 56

Diagnosis 84

Dialog language

Load 83

Dialog type 70

Dipping 130, 131, 132

Dipping in special 29, 30

Directory structure 79

Distributor

Mount 16

Dosimat 24, 25, 34

Dosing 34

Continuous (DOS) 63

Fixed volume (ADD) 65

Dosing rate

Maximum 139

Dosino

Connecting 13

Mounting 11

Double dipping 31

Dripping time 130, 131, 132, 133,

134

E

Electrode 14

Connect 21

Electrode cable 14

Electrode connector

pH electrode 9

Polarizable electrode 9

Redox electrode 9

Reference electrode 9

Electrostatic charge 5

EMPTY

Parameters 139

Emptying

Parameters 139

ERC 96, 110

Evaluation

DET 93, 95

MET 107, 110

Expert dialog 70

F

FEP tubing 14

FEP tubings 15

Filling rate

Maximum 139

Formula editor 39

G

GLP 135

Graphics width 81

Guarantee 153

I

Initialization 143

Initial position 73

Instrument

Switch off 36

Switch on 35

Instrument diagnosis 83

K

Keyboard

Connecting 18

Key assignment 142

Keyboard layout 81

Keypad

Key assignment 141

L

Language 69

Language file

Load 83

Lift

Move 62

Settings 72

Live modifications 53

Live parameters 55

M

Mains voltage 5

Maintenance 135

Manual control

Dosing 62

Measuring 66

Stir 67

MET 2, 28

Parameters 100

Method 41

Create 41

Delete 79

Export 44

Import 79

Load 43

Save 42

Index

MSB
Connector 9

N

Navigating 37
Next sample pos. 49
Number of samples 49, 51
Numerical input 38

O

Operation
General 36

P

Parameters
CAL 127
DET 85
MET 100
Modify live 55
SET 114
PC/LIMS report 80
Pin assignment 145
PREP
Parameters 139
Preparing
Parameters 139
Print 59
Printer 143
Connecting 18
Printer type 80
Program version
Update 83
Propeller stirrer 17
PTFE tubings 16
Pumps 33

Q

Quality Management 135

R

Recalculate 57
Remote
Cables 23
Connection 24
Connector 9
Interface 145
Pin assignment 145
Status diagram 146
Report
Print manually 59
Reports 98, 112, 125, 129
Result 56
Save as common variable 78,
97, 112, 124
Save as titer 97, 111, 124

Rinsing 133, 134
Rinsing beaker 8
Rinsing in sample 32
Rinsing in special 32
Rinsing time 130, 131, 132, 133,
134
Rod stirrer 14, 17
Routine dialog 70

S

Safety instructions 4
Safety shield 8, 17
Mounting 26
Removing 11
Sample data 45
Modify live 53
Sample rack
Rotating 61
Sample series
Cancel 50
Interrupt 50
Perform 49
Prolong 51
Resume 51
Shorten 51
Start 49
Stop 50
Sample size
Send from a balance 140
Sending from the balance ... 47
Sample table 45
Sensor 73
Add 74
Connect 21
Delete 74
Edit 73
Sensor connector
pH electrode 9
Polarizable electrode 9
Redox electrode 9
Reference electrode 9
Temperature sensor 9
Service 4
Service Agreement 135
SET 2, 29
Parameters 114
Solution 75
Add 76
Delete 76
Edit 76
Special position 8
Spiral band 16
Splash protection 8
Spray nozzles 15

Start conditions
DET 85
MET 100
SET 115
Statistics 58, 98, 112, 125
Stirrer
Connecting 17
Mount 14
Stirring propeller 14
Stirring rate 140
Stop conditions
DET 91
MET 106
SET 122
Switch off 36
Switch on 35
System initialization 143

T

Temperature sensor
Connect 21
Text input 38
Titrant 11
Titration head
Setting up 14
With rinsing equipment 15
Titration mode
DET 2
MET 2
SET 2
Titration parameters
DET 86
MET 101
SET 117

U

Update 83
USB 141
Keyboard 142
Keypad 141
USB (OTG)
Connector 9
USB device
Adapter 18
Connecting 18
User name 69

V

Validation 135

W

Warranty 153
Working height 72
Work position 72