

875 KF Gas Analyzer



Manual
8.875.8001EN



Metrohm AG
CH-9100 Herisau
Switzerland
Phone +41 71 353 85 85
Fax +41 71 353 89 01
info@metrohm.com
www.metrohm.com

875 KF Gas Analyzer

Manual

Teachware
Metrohm AG
CH-9100 Herisau
teachware@metrohm.com

This documentation is protected by copyright. All rights reserved.

Although all the information given in this documentation has been checked with great care, errors cannot be entirely excluded. Should you notice any mistakes please send us your comments using the address given above.

Documentation in additional languages can be found on
<http://documents.metrohm.com>.

[illegible]

■■■■■■■ |||

4.2.2	Working steps for carrying out a measurement	20
4.2.3	Explanations regarding the shape of the gas flow and titra- tion curves	22
4.2.4	Method "Reference measurement"	23
4.2.5	Changing the gas type	23
4.2.6	Calibrating a new gas type	24
4.2.7	Automatic addition of methanol, automatic reagent replace- ment (optional accessories)	31
4.2.8	Rinsing with solvent (optional accessories)	3
4.3	QUICKSTOP module	32
5	Operation and maintenance	33
5.1	General notes	33
5.1.1	Care	33
5.1.2	Maintenance by Metrohm Service	34
5.2	Quality Management and validation with Metrohm	34
6	Troubleshooting	36
7	Technical specifications	38
7.1	Temperature ranges	38
7.2	Pressure ranges	38
7.3	Supply voltage	38
7.4	Safety specifications	38
7.5	Electromagnetic compatibility (EMC)	39
7.6	Dimensions	39
7.7	Weight	39
8	Conformity and warranty	40
8.1	Quality Management Principles	40
8.2	Warranty (Guarantee)	41
9	Accessories	43
9.1	Scope of delivery	43
9.2	Spare parts for the basic unit	44
9.3	Spare parts for the base plate	44
9.4	Spare parts for integrating the 851 Titrandos	45
9.5	Spare parts for the 851 Titrandos	45
9.6	Optional accessories	45
Index		47

1 Introduction

1.1 Instrument description

The 875 KF Gas Analyzer is a robust, modularly designed analysis system based on **tiamo**™ for routine analysis at a site.

The system described on the following pages has been devised for the coulometric water content determination according to Karl Fischer in gases and allows the analysis of both liquefied gases and permanent gases. This method is also suitable for very low water contents.

The system comprises an operating unit and an analysis module. The analysis module is equipped with a base plate to convey the gas and with a water content determination cell as well as internally with an 851 Titrand in order to carry out all required analysis steps fully automatically. For this process, an amount of gas defined by the user is precisely measured with the flow meter and fed to the connected coulometer cell. Sample residue and water that might be present in the piping system are rinsed with dry nitrogen. The water is absorbed by the coulometric reagent and determined there by way of Karl Fischer titration. In coulometry, the iodine required for titration is produced by anodic oxidation, and the water content is subsequently determined. For the determination of liquefied gases, the samples are first vaporized in a controlled manner and then conveyed to sample determination.

Please also refer to the manuals and the documentation regarding the individual components (851 Titrand, mass flow controller, individual components) in addition to this documentation of the KF Gas Analyzer.

1.2 System description

- Robust analysis system with high-quality components for routine analysis tailored to the requirements of users.
- Gas-carrying system separate from the electronics and the power supply.
- The base plate with the system components is mounted behind a hood.
- The base plate comprises all components of gas conveyance and preparation as well as the coulometer cell.
- The base plate's gas system is pressure-tested.
- Nitrogen feed line with drying cartridge for predrying and check valve.
- Sample input filter preventing particles from entering the gas system.
- Deaeration bypass for pressure release during gas change.



- Integrated, adjustable vaporizer for liquefied gases.
- Heated oil filter with stainless steel filter element for analyzing used refrigerants with chiller oil contents.
- Rinsing connector for removing oil residue.
- Precise gas measurement with mass flow controller (MFC).
- Automated analysis process thanks to the use of solenoid valves.
- Predefined analysis method with a prerinsing, gas feed and postrinsing phase.
- Coulometric procedure for direct water content determination.
- Industrial PC and TFT panel (available as an option).
- All components except for the TFT panel are contained in one housing.
- Flexible control, user-friendly method creation and management and extensive data management using the **tiamo™** software. The operation of tiamo™ is described in the online help. Complete integration and control of all system components via the software.

1.3 System specification

- The system must be operated in a fume cupboard.
- Maximum sample input pressure: 40 bar.
- Maximum vaporization temperature: 80 °C.
- Nitrogen is required as auxiliary gas. The molecular sieve is used for predrying in the 875 KF Gas Analyzer. The input pressure must correspond to the vapor pressure of the samples.
- Gas connectors for nitrogen, rinsing medium, high-pressure waste gas: 6 mm Swagelok ferrule screw connector.
- Sample gas connector: 1/16" or 6 mm Swagelok ferrule screw connector.
- Gas type: The system is suitable for the liquefied gases and permanent gases listed below. The gas system must be rinsed with nitrogen after each measurement. Additional gases may be approved on request and after testing.
 - Propane, propene, butane, butene, butadiene, LPG, natural gas
 - Dimethyl ether, ethylene oxide
 - Chlorinated hydrocarbons: methyl chloride, ethyl chloride, vinyl chloride
 - Refrigerants: various chlorofluorocarbon (CFC), hydrofluorocarbon (HFC) and chlorinated hydrocarbon (CHC) compounds.
Fresh and used refrigerants with chiller oil contents.
- Safety specification: degree of protection IP54.





Note

The materials of the components used have been carefully selected in accordance with the aforementioned gases. According to the current state of technology and the material manufacturers' resistance lists, these materials are resistant to the aforementioned gases.

However, a general guarantee is impossible to give, as we cannot predict how the gas mixtures will behave in the system and we do not know the concentration, degree of purity and aggregate state of the various gases that flow through the system.

1.4 About the documentation




Caution

Please read through this documentation carefully before putting the instrument into operation. The documentation contains information and warnings which the user must follow in order to ensure safe operation of the instrument.

1.4.1 Symbols and conventions

The following symbols and formatting may appear in this documentation:

(5-12)	Cross-reference to figure legend The first number refers to the figure number, the second to the instrument part in the figure.
1	Instruction step Carry out these steps in the sequence shown.
Method	Dialog text, parameter in the software
File ► New	Menu or menu item
[Next]	Button or key
	Warning This symbol draws attention to a possible life hazard or risk of injury.



Warning

This symbol draws attention to a possible hazard due to electrical current.



Warning

This symbol draws attention to a possible hazard due to heat or hot instrument parts.



Warning

This symbol draws attention to a possible biological hazard.



Caution

This symbol draws attention to a possible damage of instruments or instrument parts.



Note

This symbol marks additional information and tips.

1.5 Safety instructions

1.5.1 General notes on safety



Warning

This instrument may only be operated in accordance with the specifications in this documentation.

The present system is suitable for processing gases and liquefied gases. In addition, hazardous substances are used in the wet end. Usage therefore requires the user to have basic knowledge and experience in handling liquefied gases, gases and pressurized media. Knowledge with respect to the application of the fire prevention measures prescribed for laboratories is also mandatory. The system may be operated only by trained staff. The operator must be trained both with regard to these operating instructions and the customer's laboratory rules and regulations.

This instrument has left the factory in a flawless state in terms of technical safety. To maintain this state and ensure non-hazardous operation of the instrument, the following instructions must be observed carefully.



The gas system is under pressure. It contains both pressurized gases and liquefied gases.

Before the sample vessel can be changed, the pressure must be released in the piping system and the latter may need to be rinsed with nitrogen.

Observe the applicable regulations.



The oven used for vaporizing the liquefied gases and the oil filter downstream of the oven may exhibit a temperature of up to 70 °C. Avoid direct skin contact. Wear heat-insulating gloves, if necessary.

Clean the oil filter and rinse the piping carrying gas through the oven only with the instrument switched off and while it is cold.

The electrical safety when working with the instrument is ensured as part of the international standard IEC 61010.



Only personnel qualified by Metrohm are authorized to carry out service work on electronic components.



Never open the housing of the instrument. The instrument could be damaged by this. There is also a risk of serious injury if live components are touched.

There are no parts inside the housing which can be serviced or replaced by the user.



Mains voltage



Warning

An incorrect mains voltage can damage the instrument.

Only operate this instrument with a mains voltage specified for it (see rear panel of the instrument).

Protection against electrostatic charges



Warning

Electronic components are sensitive to electrostatic charges and can be destroyed by discharges.

Do not fail to pull the mains cable out of the mains connection socket before you set up or disconnect electrical plug connections at the rear of the instrument.

1.5.3 Flammable solvents and chemicals




Warning

All relevant safety measures are to be observed when working with flammable solvents and chemicals.

- The instrument must be set up in a fume cupboard.
- Keep all sources of flame far from the workplace.
- Clean up spilled liquids and solids immediately.
- Follow the safety instructions of the chemical manufacturer.

1.5.4 Recycling and disposal



This product is covered by European Directive 2002/96/EC, WEEE – Waste from Electrical and Electronic Equipment.

The correct disposal of your old equipment will help to prevent negative effects on the environment and public health.

More details about the disposal of your old equipment can be obtained from your local authorities, from waste disposal companies or from your local dealer.

2 Overview of the instrument

2.1 Instruments




Figure 1 Operating unit and analysis module

2.2 Piping diagram




Figure 2 Schematic arrangement of the system

A Nitrogen	B Rinsing with solvent
-------------------	-------------------------------

Terminal	Function	Port	Port description
KL2424-2-5	A2	DigOut_1_2_2	Valve 2 - sample
KL2424-2-6	GND		
KL2424-2-7	GND		
KL2424-2-8	A4	DigOut_1_2_4	Valve 4 - measurement
Protective ground conductor terminal, 4-pin	Earth	Terminals 1 - 4	Earth for each of the 4 valves
KL2424-3-1	A1	DigOut_1_3_1	-
KL2424-3-2	GND		
KL2424-3-3	GND		
KL2424-3-4	A3	DigOut_1_3_3	MFC
KL2424-3-5	A2	DigOut_1_3_2	Heater
KL2424-3-6	GND		
KL2424-3-7	GND		
KL2424-3-8	A4	DigOut_1_3_4	-

Analog inputs

Table 3 Analog inputs

Terminal	Function	Port	Port description
KL3204-4-1	+I1	AnIn_1_4_1	Oven temperature
KL2424-4-2			
KL2424-4-3	+I3	AnIn_1_4_3	-
KL2424-4-4	GND		
KL2424-4-5	+I2	AnIn_1_4_2	-
KL2424-4-6	GND		
KL2424-4-7	+I4	AnIn_1_4_4	-
KL2424-1-8	GND		

3.2 General

The 875 KF Gas Analyzer is delivered in a largely preconfigured state.

As a rule, the installation steps described in the individual manuals have been carried out prior to delivery.

Additional notes are described in the subchapters below.

Fill the nitrogen drying cartridge with molecular sieve.

Establish the gas connections for nitrogen and, if required, for rinsing medium with 6 mm Swagelok ferrule screw connectors.

Establish the gas connection for the sample with 1/16" Swagelok ferrule screw connector.

Connect the high-pressure waste gas and the waste gas of the coulometer cell to the extraction system.

3.3 Power connection



Warning

The on-site supply voltage has to match the voltage specified on the 875 KF Gas Analyzer's housing.

The instrument is set to either 110 V or 230 V at the power supply unit.

The 875 KF Gas Analyzer and the operating unit are both connected to a socket with the preinstalled power feeder.



Warning

Electrical connections may only be made by authorized specialist personnel.

3.4 Connecting control lines



Warning

Always disconnect the instrument from the supply voltage.

Only shielded cables may be used for digital outputs, digital inputs, analog outputs and analog inputs.

The cable shielding must be connected to the grounding terminal.

The lines are connected directly to the I/O controller (see *Chapter 2.3, page 8*).

In order to open the contact springs, insert a 2.5 x 0.4 mm screw driver vertically into the rectangular actuation opening and press towards the LED.

A prefabricated cable has to be connected to the computer's network card directly if the 875 KF Gas Analyzer is being integrated into a LAN.




Figure 3 I/O controller

- | | | | |
|----------|--------------------------------|----------|---------------------------------|
| 1 | Analog output terminals | 2 | Digital output terminals |
| 3 | Digital input terminals | | |

3.5 Connecting the PC and the operating unit

The operating unit is connected directly to the industrial PC at the labeled locations.

The cable entry plate is screwed onto the 875 KF Gas Analyzer's housing.

3.6 Windows passwords

User	Password	Group
ProcessLab		User
Administrator	ADMINISTRATOR	Administrator
Metrohm	*****	Administrator

3.7 Gas connections



Warning

Lines must be laid in such a way that they cannot be pulled out.

Nitrogen, sample, high-pressure waste gas

The connection between sample vessel and the Gas Analyzer's sample input must be as short as possible, have as little dead volume as possible, be absolutely tight and consist of suitable material. Observe the notes in the enclosed assembly instructions from Swagelok on connecting the Swagelok ferrule screw connectors.

Waste gas lines

The waste gas lines must be routed to the exhaust air system with no counterpressure.



3.8 Drying cartridge for nitrogen

Depending on the residual water content of the nitrogen, the cartridge is to be filled with dried molecular sieve.

Secure the filling in place with a glass wool plug on both sides. In addition, insert a sieve disk (enclosed in the delivery as an accessory) on the output side (right).

3.9 851 Titrandu



Note

For installation and preparation, refer to the manual of the 851 Titrando.

Both electrodes (indicator electrode and generator electrode) are protected from being pushed out with an SGJ clip.

Given the gas flow, only the adsorber tube with enlarged bore supplied is to be used (see Chapter 9.5, page 45).

The adsorber tube and the stopper of the gas infeed tip are not secured in order to prevent an uncontrolled pressure rise.

3.10 Shutting down

If the system is shut down for an extended period of time, then the entire gas system (gas flow to the coulometer, waste gas, rinsing and bypass piping) has to be rinsed with nitrogen ("Shut down system" method) and the coulometer cell has to be cleared of reagent and rinsed with dry methanol or ethanol. The cell can then be stored in a dry place.



4 Operation

4.1 Arrangement of the gas-carrying system

The valve arrangement mounted on the front plate of the 875 KF Gas Analyzer permits a safe and complete transfer of the sample and the water contained in it into the coulometric titration cell. The diagram (*see Figure 4, page 16*) shows the schematic arrangement of the gas-carrying system.

The sample is introduced into the apparatus via valve 2 (**4-7**) and vaporized at the precision control valve (regulator). The heating block (**4-11**) compensates the heat that is lost in the system due to the enthalpy of vaporization and thus prevents the water to be analyzed from condensing or cooling.

The gas-carrying components are automatically rinsed with nitrogen that is predried in a drying cartridge (**4-1**) via valve 1 (**4-3**) before and after the sample is introduced. This nitrogen rinsing completely removes sample gas from the piping, so that no errors resulting from dead volumes can occur. Furthermore, rinsing with inert gas ensures that the water load on seals and internal metal surfaces in the apparatus is equal before and after sample introduction. Memory effects can be ruled out in this way.

The sample amount is metered with a mass flow controller (**4-5**), which records the amount of gas flowing in and regulates the volumetric flow. During the introduction of liquefied gases, no pressure may build up downstream of the precision control valve, as this would entail the risk of sample condensing upstream of and within the mass flow controller and possibly interfere with the flow control and damage the instrument. For this reason, the precision control valve should be adjusted in such a way that the setpoint value for the mass flow controller is not achieved. As an additional safety, the system is equipped with a control that closes the sample input valve if the gas flow exceeds a threshold value defined as common variable.

When a new sample is connected, the feed line is first prerinsed with sample via valve 3 (**4-12**). This is necessary because, initially, the connection fittings of gas bottles generally release water into the passing sample and the results of the first measurement without sample rinsing are generally higher. At the end of the measurement, the user can release the pressure from the sample infeed via stopcock 1 (**4-2**) in a controlled manner. The infeed line is then no longer under pressure when the gas container is disconnected.

If samples contain nonvolatile parts, such as oil contaminations, then these parts are held back by the filter element (4-**13**). Contamination of the mass flow controller is thus excluded.

A thermally conductive connection exists between the oil filter and the heating block, which significantly increases the filter temperature. The retarding effect of oils on water is reduced in this way. The filters and the vaporizer are cleaned by rinsing the lines with a suitable solvent via stopcock 2 (4-6). The corresponding dosing device forms part of the optional scope of delivery of the 875 KF Gas Analyzer.




Figure 4 Schematic arrangement of the system

A	Nitrogen	B	Rinsing with solvent
C	Sample	D	Waste gas
E	To the coulometer cell		
1	Drying cartridge (nitrogen)	2	Stopcock 1 (deaeration)
3	Valve 1 (nitrogen)	4	Check valve
5	Mass flow controller	6	Stopcock 2 (rinsing with solvent)
7	Valve 2 (sample)	8	Valve 4 (measurement)
9	Sample input filter	10	Precision control valve (vaporizer regulator)

12 Valve 3 (waste gas)

The diagram shows two stirred tank reactors, labeled A and B. Reactor A is a vertical cylindrical tank with a rounded bottom, containing a liquid with blue wavy lines representing the stirrer. It has a single outlet pipe at the bottom leading to a control valve with two ports labeled 1 and 2 in red. Reactor B is a vertical cylindrical tank with a rounded top, also containing a liquid with blue wavy lines. It has a single outlet pipe at the top leading to a control valve with two ports labeled 1 and 2 in red.

Figure 5 Sample vessel connector

B Gas phase

2 Sample input to 875 KF Gas Analyzer

4.2 Methods



Warning

The gas system is under pressure. It contains both pressurized gas and liquefied gas.

The prescribed analysis procedure may not be modified. Users must have detailed knowledge of the gas conveyance in order to use the manual operation. Uncontrolled operation of the valves may result in a sudden vaporization of the liquefied gas or in pressure surges.



Note

The correct position of the precision control valve has a decisive effect on the precision of the analysis. The exact position has to be determined for each gas type.

As standard, the 875 KF Gas Analyzer is delivered with the following methods (control programs of the **tiamo™** software):

- Sample measurement
- Reference measurement
- Precision control valve setting
- Gas calibration_liquefied gas

- Gas calibration_gas
- Shut down system
- Drift diagnosis
- System preparation

The following methods form part of the optional scope of delivery:

- Rinsing with solvent
- Reagent replacement
- Addition of methanol



Note

Please note:

The **tiamo**[™] method can only be run if the **Flow** program has been started.

4.2.1 Sequence of the "Sample measurement" method

The water content determination of the samples is controlled by the **Sample measurement** method, which basically consists of three steps:

- Prerinsing the line route with nitrogen
- Feeding in the sample
- Postrinsing with nitrogen

The method is designed in such a way that the pressure prevailing in the area before the regulator (line volume between precision control, nitrogen and sample valve) is released during the change from prerinsing to sample introduction and from sample introduction to postrinsing. In this way, a mixing of nitrogen and sample that could result in faulty measurements is prevented. The entire sequence is shown in (*see Table 4, page 18*).

The flow diagrams of the analysis are visualized in figure 6. Some partial steps are only run through if the corresponding scans are set to "yes" in the sample table. The dosing device for methanol addition and reagent replacement as well as for the automated rinsing with solvent is an optional equipment of the 875 KF Gas Analyzer.

Table 4 Gas conveyance and valve control during the analysis

Partial step	Condition	Opened valves	Stop condition
Prerinsing with sample	Method variable "first sample measurement?" is set to "yes"	Sample valve Waste gas valve	90 seconds expired

Partial step	Condition	Opened valves	Stop condition
Draining of the sample that flowed into the area upstream of the regulator	Method variable "first sample measurement?" is set to "yes"	Waste gas valve	60 seconds expired
Rinsing out the waste gas line with nitrogen	Method variable "first sample measurement?" is set to "yes"	Nitrogen valve Waste gas valve	45 seconds expired
Prerinsing with nitrogen	None	Nitrogen valve Measurement valve	Status message from the coulometer "Conditioning OK", but at least 60 seconds
Pressure release nitrogen	None	Measurement valve	20 seconds expired
Sample introduction	None	Sample input valve Measurement valve	Value entered for minimum sample amount (mg) in the method variable is achieved
Pressure release sample	None	Measurement valve	Gas flow falls below 30 mL/min for more than 6 seconds
Postrinsing with nitrogen	None	Nitrogen valve Measurement valve	Stop criteria of the coulometric KF titration are met (extraction time and relative drift)
Relieving the sample infeed	Method variable "disconnect gas container after measurement?" is set to "yes"	Sample valve Stopcock 1	




Figure 6 Schematic representation of the gas flows during an analysis

Red marking = sample flow		Green marking = nitrogen	
1	Prerinsing with sample	2	Rinsing out the waste gas line with nitrogen
3	Prerinsing and postrinsing with nitrogen	4	Sample introduction
5	Relieving the sample infeed	6	Rinsing the feed line with nitrogen

4.2.2 Working steps for carrying out a measurement

Load the sample table **Standard sample table gas measurement** in the run window of your **tiamo™** workplace under **Determination series ▶ Sample table ▶ Load**. This sample table is preset in such a way that you can make the entries that are relevant for you. The input window opens by double-clicking in the first line of the table template.

Method: sample measurement

gas type: LPG

sample number: 19649

minimum sample (mg): 1000

first sample measurement?: yes

disconnect sample gas after measurement?: yes

Line: 1 of 15

Apply Close

☐ Display application note

Gas type

Designation of a sample (substance or substance mixture), such as butadiene or propane, selected from the drop-down bar. The gas type is linked to the calibration factor that is stored under the same name as common variable.

Sample number

Sample ID used to identify your sample. The designation may be changed. It is also possible to assign further sample identifications. These must be created in the method and in the sample table.

Minimum sample amount

Valve 2 closes after the amount of sample entered in this field has been fed in.

Sample infeed is only completed after the sample contained in the area upstream of the regulator has flowed out.

Recommended range: approx. 1,000 to 2,500 mg, depending on the water content.

First sample measurement?

(yes/no)

Enter **yes** here in the case of the first measurement after a gas bottle has been connected. In this case, the feed line is rinsed with sample first.

Disconnect gas container after measurement?

(yes/no)

Enter **yes** here if you would like to disconnect the gas bottle after the measurement. The pressure is then released from the feed line via valve 1 in a controlled manner after the analysis and the feed line is subsequently rinsed with nitrogen.

4.2.3 Explanations regarding the shape of the gas flow and titration curves

The analysis procedure described in (*see Chapter 4.2.1, page 18*) results in a characteristic shape of the gas flow and titration curves. The sample infeed phase concludes with the gas flow dropping to a value close to zero. The titration rate (drift) follows this drop with a delay of approx. 10 seconds. If the gas flow is below a threshold value defined as common variable for 6 seconds, then the nitrogen valve opens and post-rinsing commences.

The amount of water detected in the postrinsing phase increases if the samples contain nonvolatile components that remain in the vaporizer and the oil filter. The distribution of the liquid and the gas phase balances out during the infeed phase, so that, at the end of the infeed phase, a part of the water contained in the sample is still present in the instrument's piping. Postrinsing serves to remove the retained water. Hydrophilic, nonvolatile sample components, such as glycol ester oils used in the refrigerant industry, for instance, therefore lead to a flattening of the drift curve during the infeed phase and as a result to an extension of the analysis time. As a general rule, the minimum titration time (extraction time) has to extend beyond the beginning of the postrinsing phase, as the titration would otherwise be finished in the "trough" between infeed and postrinsing. The control program uses the following formula to calculate the extraction time:

$$t_e = 60 \cdot \left(\frac{m \cdot 1000 \text{ mg}}{v} \right) \cdot t_n \cdot m / 6000$$

Figure 7 Formula for calculating the extraction time

te	Extraction time	m	Minimum sample amount in mg
v	value in mg/min saved under CV.mean mass flow	tn	Value in sec entered under CV.time for postrinsing

The default value of the **time for postrinsing** common variable is 3 minutes. If a sample requires a longer postrinsing phase, then the value must be increased accordingly.




Figure 8 Typical shape of the gas flow curve and drift curve

Gas type Propene	Sample amount 1.25 g
Minimum sample amount 0.5 g	Vaporization temperature 70 °C

4.2.4 Method "Reference measurement"

The trueness of the analysis can be checked by measuring water-spiked reference gases using the **Reference measurement** method.

Control gases with certified water contents are commercially available.

The **Reference measurement** method relies on the nitrogen calibration of the mass flow controller integrated in the instrument; i.e., it only delivers correct values if nitrogen is used as reference gas. The procedure for reference measurement is the same as the one applied for **sample measurement**. The result is indicated as a recovery rate in percent.

4.2.5 Changing the gas type

If the measurement of a new sample coincides with a change of the gas type, then the flow rate of the precision control valve has to be adjusted to the current sample using the **Precision control valve setting** method. This method sets the setpoint value at the MFC to the maximum value of 5 L/min and graphically displays the current flow by utilizing the internal nitrogen calibration. In order to prevent a pressure rise in the area after the regulator, the precision control valve has to be set in such a way that its vaporization rate is lower than 5 L/min and the setpoint value is not reached at the MFC. After the start of the method, follow the instructions of the text messages and adjust the precision control valve so that the gas flow is within the required limits (definition by common variable).



Note

Please note:

This method does not use the calibration factor that is assigned to this gas type. The mass flow displayed during the subsequent analysis may therefore considerably deviate from the value that was set when the precision control valve was adjusted.

4.2.6 Calibrating a new gas type

At the factory, the mass flow controller is calibrated to nitrogen. If the instrument is to be operated with a different gas, then the flow value has to be corrected by an appropriate factor. These correction factors are determined gravimetrically by letting larger amounts of gas flow through the MFC and monitoring the weight reduction of the gas container. The quotient of the gas volume indicated and the weight difference is the correction factor. This factor is in the range between 0.5 and 1.5 mL for most liquefied gases. The correction factors have to be individually determined for each flow controller using the **Gas calibration** method. This method saves the correction factor in the **tiamo™** configuration as common variable. In order to achieve a sufficient level of accuracy, the sample weight difference should have at least three significant places. The balance used therefore has to offer a corresponding resolution and maximum weight in accordance with the gas bottle size. For the determination of the calibration factor, the gas container has to be connected to the 875 KF Gas Analyzer with the flexible plastic capillary (OD 1/16") enclosed in the scope of delivery, as steel capillaries transmit vibrations to the balance.

Samples should be taken from the gas phase **Gas calibration_gas** rather than the liquid phase of the gas container **Gas calibration_liquefied gas** for calibrations, because the flow pattern is much more uniform if vaporization does not take place in the 875 KF Gas Analyzer. The **Gas calibration_liquefied gas** method is only to be used if a water content determination is to be done for the same gas container after calibration.

The procedure to determine the calibration factor is described below step by step using butadiene as an example:

- 1** You can find the correction factors for the gases you have used so far in the Common Variable subwindow in the **tiamo**™ configuration. Templates with the designation "additional gas type x" (x = 1 to 9) are stored for adding further gases. The common variables can be rendered editable via **Edit ► Properties**. Replace the blank variable **additional gas type x** with the lowest number x by the term butadiene.

Common Variable - additional gas type 4

Common Variable History

Name: additional gas type 4

Type: Number

Value: 1,0 mL/mg

Comment:

Assignment date: 2012-08-23 14:40:22 UTC+2

Assignment method: manual

User: Metrohm

☐ **Common Variable monitoring**

Validity: 999 days

Next assignment: 2015-05-19 ...

Message:

☐ Message by e-mail E-mail...

☐ Acoustic signal

Action:

☐ Record message

☒ Display message

☐ Cancel determination

OK Cancel

Common Variable - additional gas type 4

Common Variable History

Name: Butadien

Type: Number

Value: 1,0 mL/mg

Comment:

Assignment date: 2012-08-23 14:40:22 UTC+2

Assignment method: manual

User: Metrohm

☐ **Common Variable monitoring**

Validity: 999 days

Next assignment: 2015-05-19 ...

Message:

☐ Message by e-mail E-mail...

☐ Acoustic signal

Action:



☐ Record message

☒ Display message


☐ Cancel determination



OK Cancel

- 2 Enter the name of your gas type also in the **additional gas type x** text template under **Tools ► Text templates ► Text templates for ID** in the workplace of *tiamo*TM.



- 3 Open the **Gas calibration_liquefied gas** method under **File ► Open** in the Methods part of *tiamo*TM. The method consists of tracks that run from the top to the bottom. Each track is labeled with a letter. The individual commands are numbered consecutively from the top to the bottom. Search the **R4 call** command in the exit track. Double-click on the command to edit it. Overwrite the first line saying **additional gas type** by editing the line via the properties. Click on the ÷ symbol to open the formula editor. Replace the term **additional gas type x** in inverted commas with butadiene.





- 4 Edit the CALC command of the track to which the previously modified call command refers (in the example above, the track name was K additional gas type 4). Double-clicking in the calculation line opens a subwindow for the result properties. Click on the **Options** tab, select butadiene as common variable and then save the method with **File ► Save**.



CALC - K1 calculation

Command nameK1 calculation

	Result name	Formula	Unit	Decimal places	Assignment	Statistics	Result monitoring
1	K2 calibration factor ...	= 'R5.R3 volume sample gas (mL/mg)'	mL/mg	2	R520	<input type="checkbox"/>	<input type="checkbox"/>

New

Properties

Delete

Templates

OK

Cancel

Result K2 calibration factor (mL/mg) - R520

DefinitionMonitoringOptions

Result nameK2 calibration factor (mL/mg)

Formula

= 'R5.R3 volume sample gas (mL/mg)'

UnitmL/mg

Decimal places2

AssignmentR520

☐ Statistics

Description


R5.'Resultatname'[.VAL]
Resultatwert.

Save as template


OK

Cancel


28




- 5 Now load the **Gas calibration_liquefied gas** method in the sample table of your *tiamo*TM workplace. Select the designation of the gas type that was newly added and enter a target value for the sample gas volume (recommended range: approx. 20 L). This is the value that is displayed with the internal nitrogen calibration and not the actual gas volume of your sample. This value should be approx. 1.5 times the gas amount (in grams) which you want to convey through the instrument.




- 6 Tare the balance and start the method. After the target value has been reached, a prompt appears in which you have to enter the weight difference after gas infeed. The prefix does not matter for this.



- 7 Check whether a valid value is entered under the corresponding common variable in the configuration.
- 8 Load the **Sample measurement** method in the **tiamo™** Methods part with **File ► Open** and double-click on the **A6 calculation** command to open it.
- 9 Edit the line **A7 checking gas volume** via the properties and open the formula editor by clicking on the \div symbol. An if-then query (nested CASE function) then opens; in this query the additional gas type x that you have replaced with butadiene is listed twice in a row. Replace the term additional gas type x with butadiene also here and save the method using **File ► Save**. You can now select your new gas type for the subsequent analyses, and the method automatically uses the appropriate correction factor for the calculations.





4.2.7 Automatic addition of methanol, automatic reagent replacement (optional accessories)

The analyte in the coulometer cell consists mainly of methanol, which is removed to a considerable extent by the sample gas and the rinsing gas. The fill level of the measuring solution therefore decreases by approx. 8 mL per hour under normal operating conditions. In order to avoid malfunctions and faulty measurements, the missing methanol must be added regularly. This can be done manually with a syringe. Alternatively, the KF Gas Analyzer can be equipped with a dosing device to add methanol cyclically that is part of the optional scope of delivery. The rate at which the fill level decreases depends on the composition and temperature of the analyte. The fill level can be increased if necessary using the **Addition of methanol** method. The **Reagent replacement** method is used for a complete exchange of analyte.

4.2.8 Rinsing with solvent (optional accessories)

If liquefied gases contain nonvolatile components, these components precipitate in the piping of the KF Gas Analyzer. This is particularly the case for used refrigerants, which are usually contaminated with compressor oils. To prevent the sensitivity of the mass flow controller's sensors being compromised by such substances, an oil filter made of sintered stainless steel is located beneath the vaporizer. However, an inflood of larger amounts of oil results in a measurable retardation of the water in the piping and additionally increases the flow resistance of the oil filter, as its pores are covered by the oil. If samples contaminated with oil are to be measured, the system has to be rinsed with a suitable solvent from time to time.



The rinsing medium has to fulfill the following requirements:

- It has to be a suitable solvent for the nonvolatile residues.
- It has to exhibit a low boiling point, as it can be removed from the piping only by nitrogen rinsing.

Petroleum ether with a boiling range between 40 °C and 60 °C is recommended for oil contaminations. The rinsing medium is dosed with a dosing device that is optionally available. The system can be cleaned with the **Rinsing with solvent** method. The precision control valve must be entirely open during rinsing. For the subsequent sample measurements, the precision control valve has to be adjusted to the corresponding sample again using the **Precision control valve setting** method.

4.3 QUICKSTOP module

The red button on the left side of the housing resets all modules that are connected to the I/O controller to their default state (this usually means switched off), e.g., heater, valves and potential-free signal contacts.

The button locks in place and has to be pushed again to unlock.

Dosinos, stirrers and other devices that are connected directly to the 851 Titrando are not affected. They must be stopped directly in the software.

If an automatic analysis is running, then the quickstop module input can be queried in this **tiamo**™ method. Thus, the devices connected to the 851 Titrandro can also be stopped in this method.



5 Operation and maintenance

5.1 General notes

5.1.1 Care



Warning

Appropriate personal safety measures must be taken for any work during which hazardous substances may be released (e.g., removing connection piping, disassembling or modifying the gas-carrying system). Examples of these safety measures include wearing personal protective equipment in accordance with the laboratory regulations: protective glasses, gloves and clothing.

Rinse with nitrogen and release the pressure from the system prior to carrying out work on the gas system.

The 875 KF Gas Analyzer requires appropriate care. Excess contamination of the instrument may result in functional disruptions and a reduction in the lifetime of the otherwise sturdy mechanics and electronics.

Spilled chemicals and solvents should be removed immediately. Above all, the plug connectors on the rear of the instrument (in particular the power socket) should be protected from contamination.



Caution

Although this is largely prevented by design measures, the power plug should be unplugged immediately if aggressive media have penetrated the inside of the instrument, so as to avoid serious damage to the instrument electronics. In such cases, Metrohm Service must be informed.

The molecular sieve of the predrying cartridge must be exchanged at regular intervals (in accordance with the residual water content of the nitrogen used).

Please refer to the 851 Titrando manual for information on maintenance and care of the coulometer cell.

A careful visual inspection of the gas-carrying system and the wet end has to be performed before an analysis series is started (e.g., status of the coulometer cell, gas connections and exhaust lines, leak-tightness). Check all connections of the system for leakage at regular intervals and particu-

larly after having made any modifications. If leakage is detected, this has to be eliminated immediately so as to prevent instrument damages.

If the necessity to clean the oil filter should arise periodically as a result of analyzing liquefied gases with nonvolatile components, the **rinsing with solvent** (see Chapter 4.2.8, page 31) option is particularly recommended. Given the automated rinsing, no mechanical work is required on the gas-carrying system. The risk of leakage is thus eliminated. If the filter is cleaned manually, the system's tightness should be checked again after the filter is built in, like after any changes to the gas system.



Note

The nitrogen inlet's check valve, which is a safety feature in case of an operating error, must be subjected to a functional check at least once a year. It has to be checked whether an additional check valve is required for the nitrogen supply.

5.1.2 Maintenance by Metrohm Service

Maintenance of the 875 KF Gas Analyzer is best carried out as part of an annual service, which is performed by specialist personnel of the Metrohm company. If you are frequently working with caustic and corrosive chemicals, we recommend a shorter maintenance interval.

Metrohm Service offers every form of technical advice for maintenance and service of all Metrohm instruments.

5.2 Quality Management and validation with Metrohm

Quality Management

Metrohm offers you comprehensive support in implementing quality management measures for instruments and software. Further information on this can be found in the brochure **«Quality Management with Metrohm»** available from your local Metrohm agent.

Validation

Please contact your local Metrohm agent for support in validating instruments and software. Here you can also obtain validation documentation to provide help for carrying out the **Installation Qualification** (IQ) and the **Operational Qualification** (OQ). IQ and OQ are also offered as a service by the Metrohm agents. In addition, various application bulletins are also available on the subject, which also contain **Standard Operating Procedures** (SOP) for testing analytical measuring instruments for reproducibility and correctness.

Maintenance

Electronic and mechanical functional groups in Metrohm instruments can and should be checked as part of regular maintenance by specialist personnel from Metrohm. Please ask your local Metrohm agent regarding the precise terms and conditions involved in concluding a corresponding maintenance agreement.



Note

You can find information on the subjects of quality management, validation and maintenance as well as an overview of the documents currently available at www.metrohm.com/com/ under **Support**.

In an equilibrated state, the cell drift lies in a range between 1 and 3 $\mu\text{g}/\text{min}$. If a volumetric stream of 1 L/min of nitrogen that has been dried through the molecular sieve is set for the titration cell, the cell drift increases to approx. 2 to 4 $\mu\text{g}/\text{min}$.

The carrier gas' share in the total drift can be determined with the **Drift diagnosis** method. This share should not exceed 2 µg/min.

Cause	Remedy
Cell drift rise due to the infeed of reactive matrix components	Exchange the anolyte
Cell drift rise due to the accumulation of water and H ₂ S in the catholyte	Exchange the catholyte
Water concentration rise in the rinsing gas due to exhaustion of the molecular sieve	Check the nitrogen quality, exchange the molecular sieve cartridge
Retardation of the water due to accumulation in the vaporizer and the oil filter	Rinse the gas-carrying system with solvent




Figure 9 Systematic procedure for identifying the cause of drift rises

8 Conformity and warranty

8.1 Quality Management Principles

Metrohm Ltd. holds the ISO 9001:2000 Certificate, registration number 10872-02, issued by SQS (Swiss Association for Quality and Management Systems). Internal and external audits are carried out periodically to assure that the standards defined by Metrohm's QM Manual are maintained.

The steps involved in the design, manufacture and servicing of instruments are fully documented and the resulting reports are archived for ten years. The development of software for PCs and instruments is also duly documented and the documents and source codes are archived. Both remain the possession of Metrohm. A non-disclosure agreement may be asked to be provided by those requiring access to them.

The implementation of the ISO 9001:2000 quality management system is described in Metrohm's QM Manual, which comprises detailed instructions on the following fields of activity:

Instrument development

The organization of the instrument design, its planning and the intermediate controls are fully documented and traceable. Laboratory testing accompanies all phases of instrument development.

Software development

Software development occurs in terms of the software life cycle. Tests are performed to detect programming errors and to assess the program's functionality in a laboratory environment.

Components

All components used in the Metrohm instruments have to satisfy the quality standards that are defined and implemented for our products. Suppliers of components are audited by Metrohm as the need arises.

Manufacture

The measures put into practice in the production of our instruments guarantee a constant quality standard. Production planning and manufacturing procedures, maintenance of production means and testing of components, intermediate and finished products are prescribed.

Customer support and service

Customer support involves all phases of instrument acquisition and use by the customer, i.e. consulting to define the adequate equipment for the analytical problem at hand, delivery of the equipment, user manuals, train-

ing, after-sales service and processing of customer complaints. The Metrohm service organization is equipped to support customers in implementing standards such as GLP, GMP, ISO 900X, in performing Operational Qualification and Performance Verification of the system components or in carrying out the System Validation for the quantitative determination of a substance in a given matrix.

8.2 Warranty (Guarantee)

Metrohm guarantees that the deliveries and services it provides are free of errors in materials, design or manufacturing.

The general warranty period is 36 months (exclusions below) from the date of delivery or 18 months in the event of continuous operation. The warranty remains valid on the condition that the servicing is provided by a Service Organization authorized by Metrohm at defined intervals and with a defined scope.

The warranty period for anion suppressors is 120 months from the date of delivery or 60 months in the event of continuous operation.

The warranty period for IC separation columns is 90 days after start-up.

For third-party components that are recognizable as such, the manufacturer's warranty regulations apply.

Consumables and materials with limited storage life and glass breakage in the case of electrodes or other glass parts are excluded from the warranty.

Warranty claims cannot be asserted if the customer has failed to meet his payment obligations according to schedule.

During the warranty period, Metrohm undertakes either to replace free of charge or to credit the purchaser for any assemblies or components that can be shown to be faulty. Any transport or customs fees that may apply are the ordering party's responsibility.

The precondition for this is that the ordering party must use the Return Material Authorization (RMA) to report the faulty part, along with specification of the article number, the article designation, an adequate error description, the delivery date and (if applicable) the serial number or the chip data, respectively. In addition, the ordering party undertakes to store the faulty part for at least 24 months in accordance with current storage directives (in compliance with ESD guidelines) and to hold it in readiness for onsite inspection or for return shipment to Metrohm. Metrohm reserves the right to invoice the ordering party for these articles, including retroactively, in the event of noncompliance with these pre-conditions.

The original warranty periods for the original part apply to parts that are replaced or repaired under the above-referenced warranties (no extension of the warranty period).



Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty.

Metrohm also offers a 120-month spare parts availability guarantee and a 60-month PC software support warranty, calculated from the date on which the product is withdrawn from the market. The content of this warranty is the ability of the customer to obtain functioning spare parts or appropriate software support at market prices during the time of the warranty period.

If Metrohm AG is unable to meet this obligation due to circumstances beyond the control of Metrohm AG, then the ordering party shall be offered alternative solutions at preferential conditions.

9 Accessories



Accessories and spare parts for the 875 KF Gas Analyzer are exclusively available from Metrohm Germany.



Subject to change without notice.

9.1 Scope of delivery



After receiving the instrument, check the shipment to ensure that it is complete.

Qty.	Order no.	Description
		851 Titrando and coulometer cell
		Please refer to the 851 Titrando manual for standard accessories.

Accessories 875 KF Gas Analyzer analysis module

Qty.	Order no.	Description
1 m	6.1803.040	PTFE capillary tubing 1/16", 0.5 mm ID
2	ZPLGA01010	Sieve disks for drying cartridge
1	ZPLGA01000	Connection set, consisting of: The following articles form part of this set.
1 m	ZPLGA01020	Tube 1/16" * 0.0147"
1	ZPLGA01030	Reduction nozzle from 6 mm to 1/16"
3	ZPLGA01040	Ferrule set 1/16"

Qty.	Order no.	Description
	ZPLGA10210	MFC (Mass Flow Controller)
	ZPLGA10420	Heating cartridge
	ZPLGA10440	Bimetal switch
	ZPLGA10430	Resistance thermometer

9.4 Spare parts for integrating the 851 Titrand

Qty.	Order no.	Description
	6.1820.020	M6-M10 screw connector
	6.1808.020	Tubing olive with M6
	6.1805.090	M6 FEP tubing connection, 31 cm
	6.1805.120	M6 FEP tubing connection, 100 cm

9.5 Spare parts for the 851 Titrand

Refer to the 851 Titrand manual.

Qty.	Order no.	Description
	ZPLGA10700	Adsorber tube coulometer cell with large bore

9.6 Optional accessories

For automatically replacing the Coulomat reagent and adding methanol for continuous operation.

Qty.	Order no.	Description
	2.800.0100	Dosino 800
	8.5617.000	Reagent replacement and methanol dosing

For rinsing with solvent in the presence of nonvolatile components.



Qty.	Order no.	Description
	2.800.0100	Dosino 800
	8.5617.001	Rinsing with solvent

Index

A

Arrangement 16

C

Calibration 24

Care 33

Control lines

 Connect 12

D

Design 15

Documentation 3

E

Electrical connection 11

Electrodes 14

Electrostatic charge 6

Emission 39

Extraction time 22

F

Fume cupboard 6, 10

G

Gas flow 20

GLP 34

Guarantee 41

H

High-pressure waste gas 13

I

Immunity 39

Input pressure 2

L

Laboratory location 10

Liquefied gases 1, 2

M

Mains voltage 6

Mass flow controller 2

Metrohm Service 34

Molecular sieve 2, 11, 14, 33

N

Nitrogen 13

P

Password 13

Permanent gases 1, 2

Power connection 11

Pressure range 38

Protective equipment 33

Q

Quality Management 34

S

Safety instructions 4, 6

Safety specification 38

Sample 13, 18

Sample measurement 18

Sample table 20

Sequence 18

Service 5

Service Agreement 34

Shutting down 14

Supply voltage 11, 38

T

Temperature range 38

V

Validation 34

Vaporization temperature 2

W

Warranty 41

Waste gas lines 13

Water content determination

 Karl Fischer 1

Working steps 20