890 Titrando



Manual 8.890.8002EN





Metrohm AG CH-9100 Herisau Switzerland Phone +41 71 353 85 85 Fax +41 71 353 89 01 info@metrohm.com www.metrohm.com

890 Titrando

Manual

8.890.8002EN

11.2011 ek

Teachware Metrohm AG CH-9100 Herisau teachware@metrohm.com

This documentation is protected by copyright. All rights reserved.

Although all the information given in this documentation has been checked with great care, errors cannot be entirely excluded. Should you notice any mistakes please send us your comments using the address given above.

Documentation in additional languages can be found on *http://products.metrohm.com* under **Literature/Technical documenta-tion**.

Table of contents

1	Introductio	n	1
	1.1	The Titrando system	1
	1.2	Instrument description	2
	1.3	Titration modes — Measuring modes — Dosing com- mands	2
	1.4	Instrument versions	3
	1.5 1.5.1	About the documentation	3 3
	1.6 1.6.1 1.6.2 1.6.3 1.6.4	Safety instructions General notes on safety Electrical safety Working with liquids Recycling and disposal	4 4 5 5
2	Overview o	of the instrument	7
3	Installatior	1	9
	3.1 3.1.1 3.1.2 3.1.3	Setting up the instrument Packaging Checks Location	9 9 9
	3.2 3.2.1	Connecting controller Operation	9 9
	3.3 3.3.1 3.3.2 3.3.3	Connecting MSB devices Connecting dosing devices Connecting a stirrer or titration stand Connecting a remote box	13 14 15 16
	3.4 3.4.1 3.4.2 3.4.3 3.4.4 3.4.5 3.4.6 3.4.7 3.5	Connecting USB devices	17 17 18 18 20 21 22 23 25
	3.5.1 3.5.2	General Connecting a polarizable electrode	25 25

4	4 Working with the exchange unit		26	
5	Handli	ng ai	nd maintenance	28
	!	5.1	General information	28
	1	5.1.1 5.1.2	Care Maintenance by Metrohm Service	28 28
	Į	5.2	Quality Management and validation with Metrohm	29
e	Troubl	ocho	ating	20
0	ITOUDI	esno	Drohleme	30
	,	D. I	Problems	30
7	Appen	dix		32
	-	7.1	Remote interface	32
		/.1.1	Pin assignment of the remote interface	32
8	Techni	cal s	pecifications	36
	8	8.1	Measuring interface	36
	3	8.1.1	Temperature	36
	۲	8.7.	Internal dosing device	50 37
	•	Q.2	Mains connection	37
	•	0.5	Safety specifications	57
	•	0.4 0 E	Electromagnetic compatibility (EMC)	
	•	0.5 0.6	Ambient town evolution	
	•	8.0	Amblent temperature	38
	č	8.7	Reference conditions	38
	ł	8.8	Dimensions	39
	2	8.9	Interfaces	39
9	Conformity	rmity	and warranty	40
	9	9.1	Declaration of Conformity	40
	9	9.2	Quality Management Principles	41
	9	9.3	Warranty (guarantee)	42
10) Access	ories		43
		10.1	Scope of delivery	43
		10.1.1	890 Titrando with Touch Control 28900110	43
		10.1.2	890 litrando with tiamo light 28900210	54
	Index			66

Table of figures

Figure 1	The Titrando system	1
Figure 2	Front 890 Titrando	7
Figure 3	Rear 890 Titrando	8
Figure 4	Connecting a Touch Control	10
Figure 5	Connecting a computer	12
Figure 6	MSB connectors	13
Figure 7	Connecting a dosing device	15
Figure 8	Connecting an MSB stirrer	16
Figure 9	Rod Stirrer and titration stand	16
Figure 10	Connecting a remote box	17
Figure 11	Connecting a printer	18
Figure 12	Connecting a balance	20
Figure 13	Connecting a polarizable electrode	25
Figure 14	Attaching the exchange unit	26
Figure 15	Connectors of the remote box	32
Figure 16	Pin assignment of remote socket and plug	32

1 Introduction

1.1 The Titrando system

The Titrando is the core of the modular Titrando system. Operation is carried out by a Touch Control with touch-sensitive screen ("Stand alone" titrator) or by a computer with a corresponding software.

A Titrando system can contain numerous, various kinds of devices. The following figure provides an overview of the peripheral devices you can connect to the 890 Titrando.



Figure 1 The Titrando system

Up to three control devices (Titrandos, Dosing Interfaces, USB Sample Processors etc.) can be controlled via USB connection by PC Control/Touch

Control. With the tiamo software the system can arbitrarily be extended with control devices.

You can request information on special applications in the "Application Bulletins" and "Application Notes"; available free of charge via the Metrohm agent responsible. Various monographs on the subjects of titration techniques and electrodes are also available.

Updating the device software is described in the manual for PC Control/ Touch Control or in the tiamo help, respectively.

1.2 Instrument description

The 890 Titrando has the following characteristics:

Operation

Operation is carried out by means of a touch-sensitive Touch Control or with a high-performance PC software, e.g. tiamo.

Dosing

An internal dosing drive for exchange units.

MSB connectors

Four MSB connectors (Metrohm Serial Bus) to control dosing devices (Dosimat with exchange unit or Dosino with dosing unit), stirrer or titration stand and remote boxes.

USB connectors

Two USB connectors, through which devices such as printers, PC keyboards, barcode readers or additional control devices (Titrando, USB Sample Processor, Dosing Interface, etc.) can be connected.

• **Measuring interface** A measuring interface with a connector for polarizable electrodes.

1.3 Titration modes — Measuring modes — Dosing commands

KFT

Karl Fischer titration for determining water with automatic pre- and post-conditioning. As measuring mode can be selected:

- Ipol (voltametric measurement with adjustable polarization current)
- Upol (amperometric measurement with adjustable polarization voltage)
- MEAS

As measuring mode can be selected:

- Ipol (voltametric measurement with adjustable polarization current)
- Upol (amperometric measurement with adjustable polarization voltage)

Dosing commands

The following commands for dosing can be selected:

- **ADD** (adding a predefined volume)
- **PREP** (rinsing cylinder and tubings of an exchange or dosing unit)
- **EMPTY** (for emptying cylinder and tubings)

1.4 Instrument versions

The 890 Titrando is available in two different instrument versions:

2.890.0110	Titrando with Touch Control
2.890.0210	Titrando with tiamo light

1.5 About the documentation



Caution

Please read through this documentation carefully before putting the instrument into operation. The documentation contains information and warnings which have to be followed by the user in order to ensure safe operation of the instrument.

1.5.1 Symbols and conventions

The following symbols and styles are used in this documentation:

(5- 12)	Cross-reference to figure legend	
	The first number refers to the figure number, the second to the instrument part in the figure.	
1	Instruction step	
_	Carry out these steps in the sequence shown.	
	Warning	
	This symbol draws attention to a possible life hazard or risk of injury.	
	Warning	
	This symbol draws attention to a possible hazard due to electrical current.	

Warning
This symbol draws attention to a possible hazard due to heat or hot instrument parts.
Warning
This symbol draws attention to a possible biological hazard.
Caution
This symbol draws attention to a possible damage of instruments or instrument parts.
Note
This symbol marks additional information and tips.

1.6 Safety instructions

1.6.1 General notes on safety

Warning

This instrument may only be operated in accordance with the specifications in this documentation.

This instrument has left the factory in a flawless state in terms of technical safety. To maintain this state and ensure non-hazardous operation of the instrument, the following instructions must be observed carefully.

1.6.2 Electrical safety

The electrical safety when working with the instrument is ensured as part of the international standard IEC 61010.



Warning

Only personnel qualified by Metrohm are authorized to carry out service work on electronic components.



Warning

Never open the housing of the instrument. The instrument could be damaged by this. There is also a risk of serious injury if live components are touched.

There are no parts inside the housing which can be serviced or replaced by the user.

Mains voltage



Warning

An incorrect mains voltage can damage the instrument.

Only operate this instrument with a mains voltage specified for it (see rear panel of the instrument).

Protection against electrostatic charges



Warning

Electronic components are sensitive to electrostatic charges and can be destroyed by discharges.

Always pull the mains cable out of the mains connection socket before connecting or disconnecting electrical appliances on the rear panel of the instrument.

1.6.3 Working with liquids



Caution

Periodically check all system connections for leaks. Observe the relevant regulations in respect to working with flammable and/or toxic fluids and their disposal.

1.6.4 Recycling and disposal



This product is covered by European Directive 2002/96/EC, WEEE – Waste from Electrical and Electronic Equipment.

The correct disposal of your old equipment will help to prevent negative effects on the environment and public health.

More details about the disposal of your old equipment can be obtained from your local authorities, from waste disposal companies or from your local dealer.

2 Overview of the instrument



Figure 2 Front 890 Titrando

1 **Guide openings** 2 Spindle Moves the piston of the exchange unit up For centering the exchange unit. and down. 3 **Contact pins** 4 Coupling For switching the flat cock. For the data chip. LED "On" LED "Status" 5 6 Lights up when the Titrando and a controller Indicates the current status of the internal (Touch Control or computer) are connected dosing drive. to the mains and switched on.



Figure 3 Rear 890 Titrando

 Type plate Contains specifications concerning mains voltage, instrument type and serial number.
 Controller connector (Controller) Connection of a Touch Control or PC with installed PC software. Mini DIN, 9-pin.
 Mains connection socket
 Electrode connector (Pol.)

For connecting polarizable electrodes, e.g.

double Pt electrodes. Socket F.

- 2 USB connector (USB 1 and USB 2) USB ports (type A) for connecting printer, keyboard, barcode reader, additional Titrandos, USB Sample Processor, etc.
- 4 MSB connector (MSB 1 to MSB 4) Metrohm Serial Bus. For connecting external dosing devices, stirrers to remote boxes. Mini DIN, 9-pin.
- **6** Measuring interface 1 (Input 1)

3 Installation

3.1 Setting up the instrument

3.1.1 Packaging

The instrument is supplied in highly protective special packaging together with the separately packed accessories. Keep this packaging, as only this ensures safe transportation of the instrument.

3.1.2 Checks

Immediately after receipt, check whether the shipment has arrived complete and without damage by comparing it with the delivery note.

3.1.3 Location

The instrument has been developed for operation indoors and may not be used in explosive environments.

Place the instrument in a location of the laboratory which is suitable for operation, free of vibrations, protected from corrosive atmosphere, and contamination by chemicals.

The instrument should be protected against excessive temperature fluctuations and direct sunlight.

3.2 Connecting controller

3.2.1 Operation

Two different versions are available for operating the 890 Titrando:

- With a Touch Control with a touch-sensitive screen. It forms a "standalone" instrument together with the 890 Titrando.
- A computer enables handling of the 890 Titrando with the aid of a PC software, e.g. PC Control.



Caution

Take care to ensure that the mains cable is pulled out of the mains connection socket before either setting up or disconnecting connections between the instruments.

3.2.1.1 Connecting a Touch Control

1

Connect the Touch Control as follows:

Note



The plug is protected against accidental disconnection of the cable by means of a "pull-out protection" feature. If you wish to pull out the plug, you will first need to pull back the outer plug sleeve marked with arrows.

Insert the plug of the Touch Control connection cable into the controller socket.



Figure 4 Connecting a Touch Control

2 Connect all peripheral devices.

- Connect MSB devices (see Chapter 3.3, page 13)
- Connect USB devices (see Chapter 3.4, page 17)
- **3** Connect the Titrando to the mains.
- **4** Switch on the Touch Control.

The Touch Control power supply is supplied through the Titrando. At the time of switch on system tests are run automatically on both the Titrando and the Touch Control. The LED **On** on the front of the Titrando lights up when the system test has been completed and the instrument is ready for operation.

Caution

The Touch Control must be shut down properly by switching off with the ON/OFF switch on the rear of the instrument before the current supply is interrupted. If this is not done, then there is a danger of data loss. Because of the fact that the power supply for the Touch Control is provided through the Titrando, you must never separate the Titrando from the mains (e.g. by switching off with a connector strip) before you have switched off the Touch Control.

If you would prefer not to position the Touch Control directly next to the Titrando, then you can lengthen the connection between Titrando and Touch Control using the 6.2151.010 cable. The maximum connection length permitted is 5 m.

3.2.1.2 Connecting a computer

The 890 Titrando requires a USB connection to a computer in order to be able to be controlled by a PC software. When a controller cable 6.2151.000 is used, the instrument can be connected directly, either to a USB socket on a computer, to a connected USB hub or to a different Metrohm control instrument.

Cable connection and driver installation

A driver installation is required in order to ensure that the 890 Titrando is recognized by the PC software. To accomplish this, you must comply with the procedures specified. The following steps are necessary:

1 Install software

- Insert the PC software installation CD and carry out the installation program directions.
- Exit the program if you have started it after the installation.

2 Establish cable connections

- Connect all peripheral devices to the instrument.
 - Connect MSB devices (see Chapter 3.3, page 13)
 - Connect USB devices (see Chapter 3.4, page 17)
- Connect the instrument to the mains supply if you have not already done this.

The LED **On** on the 890 Titrando is not yet illuminated!

Connect the instrument to a USB connector (Type A) of your computer (see manual of your computer). The 6.2151.000 cable is used for this purpose.



Figure 5 Connecting a computer

For Windows 2000: the instrument is recognized and the driver is installed automatically.

For Windows XP: the instrument is recognized and the installation assistant for the driver is started automatically. Select the option "Install software automatically" and click on **[Next >]**. Exit the assistant with **[Finish]**.

For Windows Vista: the instrument is recognized and the installation assistant for the driver is started automatically. Select the option "Find and install driver software". Agree to all of the requests that follow. The installation assistant will be exited automatically.

The **On** LED on the 890 Titrando lights up when the driver installation has been exited and the device is ready for operation.



Note

The plug on the instrument end of the 6.2151.000 controller cable is protected against accidental disconnection by means of a pull-out protection feature. If you wish to pull out the plug, you will first need to pull back the outer plug sleeve marked with arrows.

Registering and configuring the instrument in the PC software

The instrument must be registered in the configuration of your PC software. Once that has been done, you can then configure it according to your requirements. Proceed as follows:

1 Set up the instrument

- Start up PC software.
 The instrument is recognized automatically. The configuration dialog for the instrument is displayed.
- Make configuration settings for the instrument and its connectors.

More detailed information concerning the configuration of the instrument can be found in the documentation for the respective PC software.

3.3 Connecting MSB devices

In order to connect MSB devices, e.g. stirrers or dosing devices, Metrohm instruments are equipped with up a maximum of four connectors at what is referred to as the *Metrohm Serial Bus* (MSB). Various kinds of peripheral devices can be connected in sequence (in series, as a "daisy chain") at a single MSB connector (8-pin Mini DIN socket) and controlled simultaneously by the respective control instrument. In addition to the connection cable, stirrers and the remote box are each equipped with their own MSB socket for this purpose.

The following illustration provides an overview of the devices that can be connected to an MSB socket, along with a number of different cabling variations.



Dosimat / Dosino



The question of which peripheral devices are supported depends on the control instrument.

i

Note

When connecting MSB devices together, the following must be observed:

- Only one device of the same type can be used at a single MSB connector at one time.
- Type 700 Dosino and 685 Dosimat dosing devices cannot be connected together with other MSB instruments on a shared connector. These dosing devices must be connected separately.



Caution

Exit the control software before you plug MSB instruments in. The control instrument recognizes when it is switched on which instrument is connected at which MSB connector. The operating unit or the control software enters the connected MSB devices into the system configuration (Device manager).

MSB connections can be extended with the 6.2151.010 cable. The length of the connection must not exceed a maximum of 15 m.

3.3.1 Connecting dosing devices

Three dosing devices can be connected to the instrument (**MSB 2 to MSB 4**).

The types of dosing devices that are supported are:

- 800 Dosino
- 700 Dosino
- 805 Dosimat
- 685 Dosimat

Proceed as follows:

1 Connect a dosing device

- Exit the control software.
- Connect the connection cable to one of the sockets marked with **MSB** on the rear of the control instrument.
- Start the control software.



Figure 7 Connecting a dosing device

If you connect a stirrer and a dosing device of the type 800 Dosino (or 805 Dosimat), then it is to be recommended to connect the stirrer to **MSB 1** and the dosing device to the MSB connector of the stirrer, as these correspond to the default settings in the method templates.

3.3.2 Connecting a stirrer or titration stand

You can use a magnetic stirrer 801 Stirrer or 803 Ti Stand (stirring "from below") or the 804 Ti Stand with a rod stirrer 802 Stirrer (stirring "from above").

Connect a stirrer or a titration stand as follows:

1 Connect a stirrer or titration stand

- Exit the control software.
- Connect the connection cable of the magnetic stirrer or of the titration stand to one of the sockets marked with **MSB** on the rear of the control instrument.
- If desired, connect the rod stirrer to the stirrer socket (with stirrer symbol) of the titration stand.
- Start the control software.



Figure 8 Connecting an MSB stirrer



Figure 9 Rod Stirrer and titration stand

3.3.3 Connecting a remote box

Instruments that are controlled via remote lines and/or which send control signals via remote lines can be connected using the 6.2148.010 remote box. In addition to Metrohm, other instrument manufacturers also use similar connectors that make it possible to connect different instruments together. These interfaces are also frequently given the designations "TTL Logic", "I/O Control" or "Relay Control" and generally have a signal level of 5 volts.

Control signals are understood to be electrical line statuses or brief (> 200 ms) electrical pulses which display the operational state of an instrument or which trigger or report an event. Sequences on a variety of instruments can thus be coordinated in a single complex automation system. No exchange of data is possible, however.

Proceed as follows:

1 Connect a remote box

• Exit the control software.

- Connect the remote box connection cable to one of the sockets marked with **MSB** on the rear of the control instrument.
- Start the control software.



Figure 10 Connecting a remote box

You can, for example, connect an 849 Level Control Box (fill level monitor in a waste canister) or a 731 Relay Box (switch box for 230/110 volt alternating current sockets and low-voltage direct current outlets). The remote box also has an MSB socket at which a further MSB instrument, e.g. a dosing device or a stirrer, can be connected.

You will find precise information concerning the pin assignment of the interface on the remote box in Appendix (see Chapter 7.1, page 32).

3.4 Connecting USB devices

3.4.1 General

The 890 Titrando has two USB connectors (type A sockets) for peripheral devices with USB interfaces. The Titrando functions as a USB hub (distributor) no matter how it is operated. If you wish to connect more than two devices to the USB, you can also use an additional commercially available USB hub .



Caution

If you operate the Titrando with the aid of the Touch Control, take care to ensure that the Touch Control is switched off when you set up or disconnect connections between the various instruments. If you use a PC software to control the Titrando, you should exit the program before you set up or disconnect the USB connections.

3.4.2 Connecting a printer

Printers that are connected to the 890 Titrando with Touch Control must meet the following requirements:

- Printer languages: HP-PCL, Canon BJL Commands or Epson ESC P/2
- Printer resolution: 300 dots/inch or 360 dots/inch (Epson)
- Paper format: A4 or Letter, single-sheet feed.

Currently connectable printer models can be found on the Internet under *www.titrando.com*.

Connect the printer as follows:

- **1** Switch off the Touch Control.
- **2** With the aid of the 6.2151.020 cable, connect the USB connector of the Titrando (type A) with the USB connector of the printer (type B, see manual for the printer).
- **3** Switch on the printer first, then the Touch Control.
- **4** Configure the printer in the device manager of the Touch Control (see manual for PC Control/Touch Control).



Figure 11 Connecting a printer

3.4.3 Connecting a balance

If you control the 890 Titrando with a PC software, connect the balance directly to the serial port (COM) on the computer. This is usually 9-pin and marked with the symbol **IOIOI**. If you operate the Titrando with the help of the Touch Control, you will need the 6.2148.020 RS-232/USB Box for connecting a balance.

The following table offers an overview of the balances that you can use together with the Titrando system and of which cable you will need for connection to the RS-232 interface:

Balance	Cable	
AND ER-60, 120, 180, 182FR-200, 300FX 200, 300, 320 with RS-232 interface (OP-03)	6.2125.020 + 6.2125.010	
Mettler AB, AG, PR (LC-RS9)	In the scope of delivery for the balance	
Mettler AM, PM, PE with interface	6.2146.020 + 6.2125.010	
option 016	also from Mettler: ME 47473	
or	adapter and either ME 42500	
Mettler AJ, PJ with interface option 018	switch	
Mettler AT	6.2146.020 + 6.2125.010	
	also from Mettler: ME 42500 hand switch or ME 46278 foot switch	
Mettler AX, MX, UMX, PG, AB-S, PB-S	6.2134.120	
Mettler AE with interface option	6.2125.020 + 6.2125.010	
011 or 012	also from Mettler: ME 42500 hand switch or ME 46278 foot switch	
Ohaus Voyager, Explorer, Analyti- cal Plus	Cable AS017-09 from Ohaus	
Precisa balances with RS-232-C interface	6.2125.080 + 6.2125.010	
Sartorius MP8, MC1	6.2134.060	
Shimadzu BX, BW	6.2125.080 + 6.2125.010	

Operation with Touch Control

Connect the balance as follows:

- **1** With the aid of the 6.2151.030 cable, connect the USB connector of the Titrando (type A) with the USB connector of the RS-232/USB Box (type B).
- **2** Connect one of the RS-interfaces of the RS-232/USB box with the RS-232 interface of the balance (see table for cable).



Figure 12 Connecting a balance

- **3** Switch on the Touch Control.
- **4** Switch on the balance.
- **5** Switch on the RS-232 interface of the balance.
- **6** Configure the RS-232 interface of the RS-232/USB Box in the device manager of the Touch Control as described in the manual for the PC Control / Touch Control.

Operation with PC software

Connect the balance as follows:

- **1** Connect the RS-232 interface of the computer with the RS-232 interface of the balance (see table for cable).
- **2** Switch on the balance.
- **3** Configure the RS-232 interface of the balance.
- **4** Configure the RS-232 interface of the computer in the device manager of the PC Control software as described in the manual for the PC Control / Touch Control.

3.4.4 Connecting a PC keyboard (only Titrando with Touch Control)

The PC keyboard is used as an aid for text and numerical input. If you operate the 890 Titrando with the aid of the Touch Control, then you can connect a PC keyboard with USB interface to the Titrando. Currently connectable keyboard models can be found on the Internet under *www.titrando.com*.

Connect the PC keyboard as follows:

- **1** Insert the USB plug of the keyboard (type A) into one of the USB sockets of the Titrando.
- **2** Switch on the Touch Control.

The keyboard will be recognized automatically and entered in the device manager.

3 Configure the keyboard in the device manager of the Touch Control as described in the manual for the PC Control / Touch Control.

3.4.5 Connecting a barcode reader

The barcode reader is used as an aid for text and numerical input. You can connect a barcode reader with USB interface. Currently connectable barcode reader models can be found on the Internet under *www.titrando.com.*

Operation with Touch Control

Connect the barcode reader as follows:

- **1** Insert the USB plug of the barcode reader (type A) into one of the USB sockets of the Titrando.
- **2** Switch on the Touch Control.

The barcode reader will be recognized automatically and entered in the device manager of the Touch Control.

3 Configure the barcode reader in the device manager as described in the manual for the PC Control / Touch Control.

Operation with a PC software

Connect the barcode reader as follows:

- **1** Insert the USB plug of the barcode reader (type A) into one of the USB sockets of the Titrando or the computer.
- 2 Start the PC software.
- **3** Configure the barcode reader in the device manager as described in the manual for the PC Control / Touch Control.

Settings on the barcode reader:

Program the barcode reader as follows (also see manual for the barcode reader).

- **1** Switch the barcode reader over to programming mode.
- **2** Specify the desired layout for the keyboard (USA, Germany, France, Spain, German-speaking Switzerland).

This setting must match the setting in the device manager (see manual for PC Control / Touch Control).

- **3** Make sure that the barcode reader is set in such a way that **[Ctrl]** characters (ASCII 00 to 31) can be sent.
- **4** Program the barcode reader in such a way that the ASCII character 02 (STX or **[Ctrl-B]**) is sent as the first character. This first character is normally referred to as the "Preamble" or "Prefix Code".
- 5 Program the barcode reader in such a way that the ASCII character 04 (EOT or [Ctrl-D]) is sent as the last character. This last character is normally referred to as the "Postamble", "Record Suffix" or "Postfix Code".

6 Exit the programming mode.

3.4.6 Connecting a USB hub

If you wish to connect more than two devices to the USB connector of the Titrando, you can also use an additional commercially available USB hub (distributor). If you operate the Titrando with the aid of the Touch Control, then you should use a USB hub with its own power supply.

Connect the USB hub as follows:

- **1** Switch off the Touch Control and/or close the PC software.
- **2** With the aid of the 6.2151.020 cable, connect the USB connector of the Titrando (type A) with the USB connector of the hub (type B, see manual for the hub).
- **3** Switch on the Touch Control, or start PC software respectively.

The USB hub is recognized automatically.

3.4.7 Connecting a Bluetooth[®] adapter

Printers and balances (or other devices with RS-232 connection) can be optionally connected using a wireless Bluetooth[®] connection. Printers and balance models with integrated Bluetooth[®] functionality are to be recommended for this purpose. Bluetooth[®] printer adapters for the USB interface can be obtained commercially, as can also Bluetooth[®] serial adapters for RS-232 connections (e.g. for balances).



Bluetooth[®] is a registered and copyrighted name of the Bluetooth[®]

Special Interest Group (Bluetooth® SIG, Inc.).

PC Control

If the Titrando system is operated with the PC Control software, then a Bluetooth USB adapter can be connected to one of the USB sockets on the computer (or of the Titrando / USB Sample Processor). The driver software (for MS Windows 2000/XP) included in the scope of delivery by the manufacturer of the Bluetooth adapter must be installed in accordance with the directions contained in the associated manual. A Bluetooth USB adapter must support the Bluetooth specifications **HCRP (Hardcopy Cable Replacement Profile)** for printers and/or **SPP (Serial Port Profile)** for balances or RS-232 connections. Printer drivers must be set up prior to the installation of the Bluetooth adapter.

Touch Control

If the Titrando system is operated as a stand-alone system with a Touch Control, then the **Metrohm Bluetooth USB adapter for the 840** (6.2162.000, Bluetooth[®] V1.1 qualified Class 2 device) is required for a Bluetooth connection.



Note

The **Metrohm Bluetooth USB adapter** cannot be operated on a computer. The adapter is intended exclusively for use in a Titrando system in stand-alone operation, i.e. with a Touch Control operating unit as controller. The prerequisite is the Touch Control Software version 5.840.0140 or higher.

The Metrohm Bluetooth USB adapter permits wireless data transfer for distances up to 10 m.

Install the Metrohm Bluetooth USB adapter as follows:

- **1** Plug the Bluetooth adapter into an open USB interface on the rear of the Titrando.
- **2** Switch on the Touch Control.

The Bluetooth USB adapter is recognized automatically.

3 Configure the adapter in the device manager of the Touch Control as described in the manual for the PC Control / Touch Control.

Printer and Bluetooth

A Bluetooth-compatible printer or a Bluetooth printer adapter must support the **HCRP (Hardcopy Cable Replacement Profile)**.

The necessary settings for a Bluetooth-capable printer can be found in the user manual for it. Bluetooth printer adapters can usually be connected without configuration to the USB connector of the respective printer. Comply with the instructions contained in the documentation for the printer adapter. Enter the definition of the printer type in the device manager of the Touch Control.

Balances and Bluetooth

A Bluetooth-capable balance or a Bluetooth serial adapter must comply with the **SPP (Serial Port Profile)** in accordance to Bluetooth specifications. If the balance manufacturer offers a specific Bluetooth serial adapter, then this is to be preferred to a commercially available adapter.

Bluetooth serial adapters must be configured on the PC with the help of an auxiliary program supplied by the manufacturer. The data transfer parameters of the instrument and the adapter must match one another. The Bluetooth serial adapter must be defined as the acceptor and not as the initiator of a serial connection. Authentication with a PIN code is not supported.

3.5 Connecting sensors

3.5.1 General

The measuring interface includes one measuring input (**Pol.**) for a polarizable electrode (*see Figure 3, page 8*).

3.5.2 Connecting a polarizable electrode

Connect a polarizable electrode as follows:

1 Plug the electrode plug into the **Pol.** socket of the 890 Titrando.





Note



The electrode cable is protected against accidental disconnection of the cable by means of a pull-out protection feature. If you wish to remove the plug, then you must first retract the outer plug sleeve.

4 Working with the exchange unit

The 806 exchange unit has an integrated data chip which allows to save data for the exchange unit and for the reagent. The data is being edited in the Touch Control or in the PC software. The assembly of the exchange unit is described in the manual for the 806 exchange unit.



Figure 14 Attaching the exchange unit

To attach the exchange unit, proceed as follows:

1 Slide the exchange unit from the front onto the 890 Titrando and push all the way to the rear, so that it snaps in and the LED **Status** flashes slowly.

It must snap in audibly.

If the exchange unit is attached correctly, its initialization is being activated by a micro switch which is triggered by the guide bolts of the exchange unit. The exchange unit is being recognized and the data is being read automatically from the data chip. After this, the LED **Status** lights up constantly.

The following table provides an overview showing all operating states of the internal dosing device which are indicated with the LED **Status**:

LED Status	Operating status of the dosing device
off	No exchange unit is attached.
flashes constantly	The Titrando is ready to dose or titrate. The exchange unit has been attached correctly and rec- ognized and is now in the exchange position, i.e. the exchange unit can be removed.
flashes slowly	The Titrando is dosing or filling or the exchange unit is not in the exchange position, respectively.
	An 806 intelligent exchange unit has been attached and the data on the integrated data chip is being read or written.
flashes fast	Error at the dosing drive (<i>see</i> <i>Chap. 5</i>)

5 Handling and maintenance

5.1 General information

5.1.1 Care

The 890 Titrando requires appropriate care. Excess contamination of the instrument may result in functional disruptions and a reduction in the service life of the sturdy mechanics and electronics.

Spilled chemicals and solvents must be removed immediately. Above all, the plug connections on the rear of the instrument (in particular the mains connection socket) should be protected from contamination.



Caution

Although this is extensively prevented by design measures, the mains plug should be unplugged immediately if aggressive media has penetrated the inside of the instrument, so as to avoid serious damage to the instrument electronics. In such cases, the Metrohm Service must be informed.

5.1.2 Maintenance by Metrohm Service

Maintenance of the 890 Titrando is best carried out as part of an annual service, which is performed by technicians of the Metrohm company. If working frequently with caustic and corrosive chemicals, a shorter maintenance interval could be necessary.

The Metrohm service department offers every form of technical advice for maintenance and service of all Metrohm instruments.
5.2 Quality Management and validation with Metrohm

Quality Management

Metrohm offers you comprehensive support in implementing quality management measures for instruments and software. Further information on this can be found in the brochure **«Quality Management with Metrohm»** available from your local Metrohm agent.

Validation

Please contact your local Metrohm agent for support in validating instruments and software. Here you can also obtain validation documentation to provide help for carrying out the **Installation Qualification** (IQ) and the **Operational Qualification** (OQ). IQ and OQ are also offered as a service by the Metrohm agents. In addition, various application bulletins are also available on the subject, which also contain **Standard Operating Procedures** (SOP) for testing analytical measuring instruments for reproducibility and correctness.

Maintenance

Electronic and mechanical functional groups in Metrohm instruments can and should be checked as part of regular maintenance by specialist personnel from Metrohm. Please ask your local Metrohm agent regarding the precise terms and conditions involved in concluding a corresponding maintenance agreement.



Note

You can find information on the subjects of quality management, validation and maintenance as well as an overview of the documents currently available at <u>www.metrohm.com/com/</u> under **Support**.

6 Troubleshooting

6.1 **Problems**

Problem	Cause	Remedy
LED "On" is not illu- minated although the Titrando is con- nected to the mains.	The Touch Control or the computer has not been switched on yet or the plugs are not correctly connected.	Check the plug connections and switch on the Touch Control or the computer.
LED "Status" is flashing fast.	The dosing drive is overloa- ded because the stopcock is jammed.	Switch off the Touch Control or exit PC Con- trol/tiamo [™] , respectively. Check whether the exchange unit can be removed. If the exchange unit cannot be removed, check whether the flat stopcock still can be turned. Switch it manually to the exchange position by turning it to the right (see manual for the exchange unit). Remove the exchange unit and proceed as described in the manual for the exchange unit.
	The dosing drive is overloa- ded because the piston is jammed. This error is being displayed in the software (Touch Control or PC Con- trol / tiamo ™).	Switch the control device off and back on again. The dosing device is being initialized during switching on. Remove the exchange unit and clean it as described in the manual for the exchange unit in chapter "Care and main- tenance". Contact the Metrohm service if removing the exchange unit is not possible.
	The data of the exchange unit cannot be read because the data chip has been damaged mechani- cally or by chemicals.	Have the data chip replaced by the Metrohm service. Until the data chip is being replaced you can remove the data chip yourself in order to be able to still use the exchange unit. The cylinder volume is automatically recognized nevertheless, but no data can be read from the exchange unit or be saved on it anymore.
LED "Status" is not illuminated although an exchange unit is attached.	The exchange unit has not correctly been attached.	Remove the exchange unit and attach it again until it snaps in. The LED flashes during data is read out from an intelligent exchange unit (806) and lights up constantly if the exchange unit has correctly been recognized.

Problem	Cause	Remedy
The exchange unit cannot be attached.	The flat stopcock of the exchange unit is not in the exchange position.	Switch the flat stopcock manually to the exchange position (switching lever directed to the right).
	The piston rod in the exchange unit is not in the right position.	Switch the piston rod to the right position (see manual for the exchange unit).
The exchange unit cannot be removed and the LED "Sta- tus" is flashing slowly.	It is currently being dosed or filled and/or the Titrando is not in the exchange position.	Stop the run or carry out a "Filling".

7 Appendix

7.1 Remote interface

The 6.2148.010 remote box allows devices to be controlled which cannot be connected directly to the MSB interface of the Titrando.



Figure 15 Connectors of the remote box

1 Cable For connecting the Titrando. 2 MSB connector Metrohm Serial Bus. For connecting external dosing devices or stirrers.

3 Remote connector

For connecting devices with a remote interface.

7.1.1 Pin assignment of the remote interface



Figure 16 Pin assignment of remote socket and plug

The above presentation of the pin assignment of a Metrohm remote interface applies not only for the remote box, but also for all Metrohm devices with 25-pin D-Sub remote connection.

Inputs



approx. 50 k Ω Pull-up ___tp___t_p >20 ms active = low, inactive = high

Outputs



active = low, inactive = high $I_{C} = 20$ mA, $V_{CEO} = 40$ V

+5 V: maximum load = 20 mA

The following tables offer information concerning the allocation of the individual pins and their function.

Table 1Inputs and outputs of the remote interface

Assignment	Pin No.	Function
Input 0	21	Start
Input 1	9	Stop
Input 2	22	Quit
Input 3	10	
Input 4	23	-
Input 5	11	
Input 6	24	
Input 7	12	
Output 0	5	Ready
Output 1	18	Conditioning OK
Output 2	4	Determination
Output 3	17	EOD
Output 4	3	
Output 5	16	Error
Output 6	1	
Output 7	2	Warning

Assignment	Pin No.	Function
Output 8	6	
Output 9	7	
Output 10	8	
Output 11	13	
Output 12	19	
Output 13	20	
0 volts / GND	14	
+5 volts	15	
0 volts / GND	25	

Table 2Explanation of the individual functions

Function	Explanation
Start	The current method is started at the time of activation.
	t _{pulse} > 100 ms
Stop	The current method is canceled (Stop) at the time of activation.
	t _{pulse} > 100 ms
Quit	The current command in the determination run will be canceled at the time of activation.
	t _{pulse} > 100 ms
Ready	The instrument is ready to receive a start sig- nal.
Conditioning OK	The line is set when Conditioning with SET and KFT titration is at OK. The line remains set until the determination is started with [START].
Determination	The instrument performs a data-generating determination.
EOD	End of Determination.
	Pulse (t _{pulse} = 200 ms) after a determination or after a buffer/standard during calibration using a Sample Processor.
Error	The line is set for error message display.

Function	Explanation
Warning	The line is set for warning message display.

Technical specifications 8

Measuring interface 8.1

The 890 Titrando has one measuring input for polarizable electrodes.

The measuring cycle is 100 ms for all measuring modes.

8.1.1 **Temperature**

One measuring input (Temp.) for a Pt1000 or NTC temperature sensor. Automatic temperature compensation for NTC sensors; R (25 °C) and B value can be configured.

Measuring range	
Pt1000	−150…+250 °C
NTC	−5…+250 °C
	(R (25 °C) = 30000 Ω and B (25/50) = 4100 K)
Resolution	
Pt1000	0.1 °C
NTC	0.1 °C
Measuring accu-	
racy	
Pt1000	±0.2 °C
	(Applies for measuring range $-20+150 \text{ °C}$; $\pm 1 \text{ digit}$; without sensor error, under reference conditions)
NTC	±0.6 °C
	(Applies for measuring range +10+40 °C; ±1 digit; without sensor error, under reference conditions)

8.1.2 Polarizer

	A measuring input (Pol.) for a polarizable electrode.
Measuring mode Ipol	Determination with adjustable polarization current
Polarization current	-122.5+122.5 μ A (in steps of 0.5 μ A) -125.0+125.0 μ A: non-guaranteed values, dependent on reference voltage +2.5 V
Measuring range	-1200+1200 mV
Resolution	0.1 mV

Measuring accuracy	±0.2 mV (±1 digit, without sensor error, under reference conditions)
Measuring mode Upol	Determination with adjustable polarization voltage
Polarization voltage	-1225+1225 mV (in steps of 25 mV) -1250+1250 mV: non-guaranteed values, dependent on reference voltage +2.5 V
Measuring range	-120+120 μA
Resolution	0.1 µA

8.2 Internal dosing device

Cylinder volume of	1 mL, 5 mL, 10 mL, 20 mL or 50 mL
the exchange unit	
Resolution	20000 steps per cylinder volume

8.3 Mains connection

Voltage	100240 V
Frequency	5060 Hz
Power consump- tion	maximum 45 W
Fuse	Electronic overload protection

8.4 Safety specifications

Design and testing

- EN/IEC 61010-1: 2001
- UL 61010-1: 2004
- CSA-C22.2 No. 61010-1: 2004
- Protection class I

Safety instructions

The documentation contains safety instructions which have to be followed by the user in order to ensure safe operation of the instrument.

8.5 Electromagnetic compatibility (EMC)

Emission

Standards fulfilled:

- EN/IEC 61326-1: 2006
- EN/IEC 61000-6-3: 2005
- EN/IEC 61000-6-4: 2005
- EN 55022 / CISPR 22: 2006

Immunity

Standards fulfilled:

- EN/IEC 61326-1: 2006
- EN/IEC 61000-6-2: 2006
- EN/IEC 61000-4-2: 2001
- EN/IEC 61000-4-3: 2002
- EN/IEC 61000-4-4: 2004
- EN/IEC 61000-4-5: 2001
- EN/IEC 61000-4-6: 2001
- EN/IEC 61000-4-11: 2004
- EN/IEC 61000-4-14: 2004
- NAMUR: 2004

8.6 Ambient temperature

Nominal function	+5+45 °C
range	
Storage	−20+60 °C
Transport	-40…+60 °C

8.7 Reference conditions

Ambient tempera- ture	+25 °C (±3 °C)
Relative humidity	≤ 60 %
Instrument status	Instrument in operation at least 30 min
Validity of the data	After adjustment

8.8 Dimensions

Width	142 mm
Height	227 mm
Depth	231 mm
Weight	2817 g
Material (Housing)	Polybutylene terephthalate (PBT)

8.9 Interfaces

USB connectors USB ports 2 USB downstream ports (model A sockets), each 500 mA, for connecting peripheral devices such as printers, keyboards, barcode readers or RS-232/USB boxes (Metrohm order no. 6.2148.020). **Controller connector** USB upstream port with auxiliary power supply (Mini DIN socket) for Controller port connecting Touch Control or computer for controlling the 890 Titrando. With integrated Touch Control cable Touch Control with cable 6.2151.000 Computer **MSB connectors (Metrohm Serial Bus)** Dosing device Connection for a maximum of 3 external dosing devices, models Dosimat or Dosino (MSB 2 to MSB 4). Stirrer Connection for a maximum of 4 stirrers. Stirrer control: switching on/off manually or coordinated with the titration sequence. Rate in 15 steps and shift direction can be selected. Remote Box Connection for a maximum of four Remote Boxes. Remote Boxes can be used to actuate and monitor external devices.

9 Conformity and warranty

9.1 Declaration of Conformity

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

Name of commodity	890 Titrando		
	Universal titrator, controlled by a Touch Control or a comp Software.		
	This instrument has been built and has undergone final type testing according to the standards:		
Electromagnetic compatibility	Emission:	EN/IEC 61326-1: 2006, EN/IEC 61000-6-3: 2004, EN 55022 / CISPR 22: 2006	
	Immunity:	EN/IEC 61326-1: 2006, EN/IEC 61000-6-2: 2005, EN/IEC 61000-4-2: 2001, EN/IEC 61000-4-3: 2002, EN/IEC 61000-4-4: 2004, EN/IEC 61000-4-5: 2001, EN/IEC 61000-4-6: 2001, EN/IEC 61000-4-11: 2004, EN/IEC 61000-4-14: 2004, NAMUR: 2004	
Safety specifications	EN/IEC 61010-1: 2001 CSA-C22.2 No. 61010	, UL 61010-1: 2004, -1: 2004, protection class I	
CE	This instrument meets the requirements of the CE mark as contain the EU directives 2006/95/EC (LVD), 2004/108/EC (EMC). It fulfils lowing specifications:		
	EN 61326-1	Electrical equipment for measurement, control and laboratory use – EMC requirements	
	EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use	
Manufacturer	Metrohm Ltd., CH-910	1 Herisau/Switzerland	
	Metrohm Ltd. is holder of the SQS certificate ISO 9001:2000 Quali agement system for development, production and sales of instrum and accessories for ion analysis.		

Herisau, 10 January, 2008

D. Strohm Vice President, Head of R&D

Pacom ann

Ch. Buchmann Vice President, Head of Production Responsible for Quality Assurance

9.2 Quality Management Principles

Metrohm Ltd. holds the ISO 9001:2000 Certificate, registration number 10872-02, issued by SQS (Swiss Association for Quality and Management Systems). Internal and external audits are carried out periodically to assure that the standards defined by Metrohm's QM Manual are maintained.

The steps involved in the design, manufacture and servicing of instruments are fully documented and the resulting reports are archived for ten years. The development of software for PCs and instruments is also duly documented and the documents and source codes are archived. Both remain the possession of Metrohm. A non-disclosure agreement may be asked to be provided by those requiring access to them.

The implementation of the ISO 9001:2000 quality management system is described in Metrohm's QM Manual, which comprises detailed instructions on the following fields of activity:

Instrument development

The organization of the instrument design, its planning and the intermediate controls are fully documented and traceable. Laboratory testing accompanies all phases of instrument development.

Software development

Software development occurs in terms of the software life cycle. Tests are performed to detect programming errors and to assess the program's functionality in a laboratory environment.

Components

All components used in the Metrohm instruments have to satisfy the quality standards that are defined and implemented for our products. Suppliers of components are audited by Metrohm as the need arises.

Manufacture

The measures put into practice in the production of our instruments guarantee a constant quality standard. Production planning and manufacturing procedures, maintenance of production means and testing of components, intermediate and finished products are prescribed.

Customer support and service

Customer support involves all phases of instrument acquisition and use by the customer, i.e. consulting to define the adequate equipment for the analytical problem at hand, delivery of the equipment, user manuals, training, after-sales service and processing of customer complaints. The Metrohm service organization is equipped to support customers in implementing standards such as GLP, GMP, ISO 900X, in performing Operational Qualification and Performance Verification of the system components or in carrying out the System Validation for the quantitative determination of a substance in a given matrix.

9.3 Warranty (guarantee)

Metrohm guarantees that the deliveries and services it provides are free from material, design or manufacturing errors. The warranty period is 36 months from the day of delivery; for day and night operation it is 18 months. The warranty remains valid on condition that the service is provided by an authorized Metrohm service organization.

Glass breakage is excluded from the warranty for electrodes and other glassware. The warranty for the accuracy corresponds to the technical specifications given in this manual. For components from third parties that make up a considerable part of our instrument, the manufacturer's warranty provisions apply. Warranty claims cannot be pursued if the Customer has not complied with the obligations to make payment on time.

During the warranty period Metrohm undertakes, at its own choice, to either repair at its own premises, free of charge, any instruments that can be shown to be faulty or to replace them. Transport costs are to the Customer's account.

Faults arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc. are expressly excluded from the warranty.

10 Accessories

10.1 Scope of delivery



10.1.1 890 Titrando with Touch Control 28900110

Qty.	Order no.	Description	n	
1	1.890.0010	890 Titrand	lo	
1	1.840.0100	Touch Cont	rol	
	Operating eleme Standard dialog guage can be ret	nt with touch-s languages Engli trofitted.	ensitive color display for Titrandos. ish and German. One further lan-	
1	1.803.0010	803 Ti-Stan	d	
	Titration stand w addition and asp	ith magnetic st iration of the ti	irrer and built-in pump for solvent tration vessel content.	
1	6.0338.100	Double Pt-v	vire electrode	
	Double Pt-wire e Karl Fischer titrat Measuring rar Measuring un Temperature e Shaft length te Indicator elect Indicator elect Plug for filling Pressure tolera	lectrode, Metro ion and 684/73 nge: it: range (°C): o head (mm): trode type : trode shape: opening: ance:	ohm plug-in head G, for volumetric 87 Coulometer (cell with diaphragm) -2000 2000 mV -20 70 96 Pt Wire None > 0	Contraction of the second

Qty.	Order no.	Description	
1	6.1244.040	Set of seals	000
1	6.1403.040	KF adsorber tube	
	Used with KF equ Height (mm): Outer diamete	ipment. Including cover and O-ring. 113 r (mm): 8	
1	6.1414.030	KF titration vessel lid	
	For volumetric KF Material: Material remar	titrations. PTFE k: Insert	
1	6.1415.220	Titration vessel / 20-90 mL	
	Material: Height (mm): Outer diamete Volume (mL):	Clear glass 80 r (mm): 78 20 90	6.1415.220

Qty.	Order no.	Description	
1	6.1415.250	Titration vessel / 50-150 mL	
	Material:	Clear glass	
	Height (mm):	80	8.1475 B
	Outer diameter	(mm): 78	-450
	Volume (mL):	50 150	
2	6.1446.040	Thread stopper / M6	
_			
	Material:	PVDF	
	Height (mm):	21.5	
	Outer diameter	(mm): 4.9	
1	6 1446 060	Stopper / B-14/15 / M10	
•	0.1440.000		
	Material:	PP	
	Height (mm):	22	
	SGJ size:	B-14/15	
1	6.1446.090	Stopper / B-14/15 / M8	
	Material:	PP	
	Height (mm):	23	
	SGJ size:	B-14/15	

Qty.	Order no.	Description	
2	6.1448.010	Septum 12 mm, 5 pieces	
	Set of 5 items. Material: Height (mm): Outer diameter	Silicone rubber 2 (mm): 12	
1	6.1543.110	Aspiration tip / M8 thread	
	Used with 676 San and 703 Ti Stand	nple Changer together with 6.1805.200 tubing	
	Material:	ETFE/PTFE	
1	6.1543.120	Aspiration tip / M8 thread	
	Aspiration tip if 6. Material: Length (mm):	1805.200 siphoning tubing is used ETFE/PTFE 151	
2	6.1602.105	Bottle neck attachment / GL 45 / green	
	For direct dosing f Material: Opening groun	rom reagent bottles with GL 45 thread. PFA/PP d joint: A-14/15	

10 Accessories

6.1602.115 For direct dosing	Bottle neck attachment / S 40 / green	
For direct dosing		
Material:	from reagent bottles with S 40 thread (Merck). PFA/PP	
6.1608.023	Amber glass bottle / 1000 mL / GL 45	
For Exchange Uni	ts. Bottle for auxiliary solutions.	
Material:	Amber alass	
Width (mm):	96	
Height (mm):	223	
Volume (mL):	1000	
6.1608.030	Round glass bottle / 1000 mL / GL 45	
To Exchange Unit	S.	
Material:	Clear glass	
Height (mm):	223	
Volume (mL):	1000	5100 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 100 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1000 1
6.1609.010	Adsorber tube with tubing connection	
Cover with tubing	g connection.	
Height (mm):	155	
Inner diameter	(mm): 32	
SGJ size:	B-14/15	Toologily of the second s
	6.1608.023 For Exchange Unit Material: Width (mm): Height (mm): Volume (mL): 6.1608.030 To Exchange Unit Material: Height (mm): Volume (mL): 6.1609.010 Cover with tubing Height (mm): Inner diameter SGJ size:	6.1608.023 Amber glass For Exchange Units. Bottle for auxiliary solutions. Material: Amber glass Width (mm): 96 Height (mm): 223 Volume (mL): 1000 6.1608.030 Round glass buttle / 1000 mL / GL 45 To Exchange Units. Clear glass Material: Clear glass Height (mm): 223 Volume (mL): 1000 Material: Clear glass Height (mm): 223 Volume (mL): 1000 Foloson of the second s

Otv	Order no	Description		
	6 1622 000	Overflow	ataction for Ti Stand	
I	6.1623.000	Overnow pr	otection for 11 Stand	
	For 803 Ti Stand			
1	6.1801.120	PVC tubing	/ 4 mm / 6 mm / 2 m	
	Material:		PVC (transparent)	
	Outer diamete	r (mm):	6	
	Inner diameter	· (mm):	4	
	Length (m):		2	
2	6.1805.200	PTFE tubing	/ M8 / 0.5 m	
	Material:		PTFE	
	Outer diamete	r (mm):	4	
	Inner diameter	· (mm):	3	
	Length (mm):		500	
				1
				щ
1	6.1808.050	Connector t	ubing nozzle M8	
	1 M8 outer threa	d and 1 tubing	olive. E.g. for thermostat jacket of	
	Exchange Units a	nd stability mea	asuring instruments.	
	Material:		PVDF	
	Length (mm):		31.5	

Qty.	Order no.	Descripti	on	
1	6.1819.030	Tube		
	For M8 threade 6.1608.0XX bot	d holes. For 6. tles	1602.1XX bottle attachments and	
	Material:		PTFE	
	Outer diamet	ter (mm):	4	
	Inner diamete	er (mm):	3	
	Length (mm)	:	250	
1	6 1819 050	Tube		
•	Eor M8 throadou	d bolos of 6 10	602 10X bottle attachments with	
	6.1608.030 asp	iration bottle.	For 703 and 803 Ti Stand, 681 Pump	
	Material:		PTEE	
	Outer diamet	ter (mm):	4	
	Inner diamete	er (mm):	3	
	Length (mm)		88	
2	6.1903.020	Stirring ba	ar / 16 mm	
	Stirring Bar with	magnetic cor	e, PTFE coated, length 16 mm	
	Material:		PTFE	
	Outer diamet	ter (mm):	4	
	Length (mm)		16	

Qty.	Order no.	Description	
2	6.1903.030	Stirring bar / 25 mm	740 500
	Stirring Bar with i Material: Outer diamete Length (mm):	magnetic core, PTFE coated, length 25 mm PTFE er (mm): 5 25	
1	6.2001.060	Base-plate without support rod	
-	For mounting stir	rers.	
1	6.2013.010	Clamping ring	
	For stand rods of Material: Width (mm): Height (mm):	10 mm diameter Metal 20 16	
1	6.2016.070	Support rod / 400 mm	ß
	Material: Outer diamete Length (mm):	Stainless steel 18/9 er (mm): 10 400	

10	Διισί	OLIDC
10	ACCC33	ULICS

Qty.	Order no.	Description	
1	6.2023.020	Clip for SGJ 14/15	
	Clip for SGJ 14/1	6	
	Material:	POM	CAR A
2	6.2043.005	Holding clip for bottles	
	Holding clip for r	eagent bottle of Exchange Units	
1	6.2104.020	Electrode cable / 1 m / F	
	For connecting e	lectrodes with Metrohm plug-in head G to Metrohm	
	instruments (sock	ret F).	
	Length (m).	I	
1	6.2412.000	Glass weighing spoon	
	Used with 6.141	4.030 Titration vessel lid or 6.1455.31X,	
	6.1464.32X, 6.14	165.320 Titration vessels. For Karl Fischer water	
	Material	Glass	
	Length (mm):	120	
	2		15/
			(C).
1	6.2621.070	Hexagon key 5 mm	
	5 mm.		
	Length (mm):	80	

---- 51

Qty.	Order no.	Description
1	6.2621.130	Hexagon key 2 mm
	2 mm	
2	6 2730 010	Screw ninnle
-	Made of PP Toge	other with 6 1448 010 Septum or for mounting
	6.1403.040 Adsc	rber tube or 6.0338.100 Electrode in 6.1414.030
	and 6.1414.040	KF Titration vessel lids. Used for KF titrations.
	Inner diameter	(mm): 11.5
	Length (mm):	23.5
	6 2720 020	Cambum abannan
Ĩ	6.2730.020	Septum stopper
	Titration vessel lic	I. Used for KF applications.
	Outer diamete	r (mm): 18
	Length (mm):	30
3	6.2730.030	Stopper
	Nipple and E.301	.0043 O-ring included. For 6.1414.030 titration ves-
	Advanced IC Dilu	tion Sample Processor
	Length (mm):	38

10 Accessories

Qty.	Order no.	Description	
1	6.2739.000	Wrench	<u> </u>
	For tightening co	nnectors	
	Length (mm):	68	
1	6.2811.000	Molecular sieve	
	Molecular sieve. I	Bottle containing 250 g. Pore size: 0.3 nm. Without	
	moisture indicato	r. For Rancimats and Karl Fischer instruments.	Same
			De Contraction of the second
			And the state of t
			Assert du Constant Constant
1	6.3026.210	Exchange unit / 10 mL	
	Exchange unit wi	th integrated data chip with 10 mL glass cylinder	
	tion, antidiffusior	on. PCTFE/PTFE flat stopcock, FEP tubing connec- buret tip and standard amber glass reagent bottle.	
	Volume (mL):	10	
			-
			5
1	6.6051.000	CompactFlash card with Titrando methods	
	Memory card wit	h application methods	
1	6 2122 0v0	Mains cable with C12 line socket	
•	0.2122.000	IEC-60320-C13	
	Cable plug accord	ding to customer requirements.	
	Switzerland:	Type SEV 12 6.2122.020	
	Germany,:	Type CEE(7), VII	
	USA,:	6.2122.040 Type NEMA/ASA	
		6.2122.070	

---- 53

Qty.	Order no.	Description
1	8.890.8002EN	890 Titrando Manual

10.1.2 890 Titrando with tiamo light 28900210

Qty.	Order no.	Description	ı	
1	1.890.0010	890 Titrand	lo	
1	1.803.0010	803 Ti-Stan	d	
	Titration stand w addition and asp	vith magnetic st piration of the ti	irrer and built-in pump for solvent tration vessel content.	
1	6.0338.100	Double Pt-v	vire electrode	
	Double Pt-wire e Karl Fischer titrat	electrode, Metro tion and 684/73	ohm plug-in head G, for volumetric 7 Coulometer (cell with diaphragm)	
	Measuring rai	nge:	-2000 2000	a a t
	Measuring un	it:	mV	
	Temperature	rang (°C):	-20 70	
	Shaft length t	o head (mm):	96	
	Indicator elec	trode type :	Pt	
	Indicator elec	trode shape:	Wire	
	Plug for filling	opening:	None	
	Pressure toler	ance:	> 0	1. Alexandre and the second se
1	6.1244.040	Set of seals		

10 Accessories

QLY.	Order no. Descript	tion	
1	6.1403.040 KF adso	rber tube	
	Used with KF equipment. Ind Height (mm): Outer diameter (mm):	cluding cover and O-ring. 113 8	
1	6.1414.030 KF titrat	ion vessel lid	13
	For volumetric KF titrations. Material: Material remark:	PTFE Insert	
1	6.1415.220 Titration	vessel / 20-90 mL	
	Material: Height (mm): Outer diameter (mm): Volume (mL):	Clear glass 80 78 20 90	6.1415.220
1	6.1415.250 Titration	vessel / 50-150 mL	
	Material: Height (mm): Outer diameter (mm): Volume (mL):	Clear glass 80 78 50 150	61115.250

Qty.	Order no.	Description	
2	6.1446.040	Thread stopper / M6	and the second s
	Material:	PVDF	
	Height (mm):	21.5	
	Outer diameter	(mm): 4.9	
1	6.1446.060	Stopper / B-14/15 / M10	
	Mataviali	••	
	Naterial: Height (mm):	77 22	
	SGL size:	B-14/15	
	505 5120.		
1	6.1446.090	Stopper / B-14/15 / M8	
	Material:	PP	
	Height (mm):	23	
	SGJ size:	B-14/15	
2	6.1448.010	Septum 12 mm, 5 pieces	
	Set of 5 items.		
	Material:	Silicone rubber	
	Height (mm):	2	
	Outer diameter	(mm): 12	

Qty.	Order no.	Description	
1	6.1543.110	Aspiration tip / M8 thread	17
	Used with 676 Sa	mple Changer together with 6.1805.200 tubing	
	Material [.]	FTFF/PTFF	
	Length (mm):	151	
1	6.1543.120	Aspiration tip / M8 thread	
	Aspiration tip if 6	.1805.200 siphoning tubing is used	
	Material:	ETFE/PTFE	
	Length (mm):	151	
2	6.1602.105	Bottle neck attachment / GL 45 / green	
	For direct dosing	from reagent bottles with GL 45 thread.	
	Material:	PFA/PP	
	Opening grour	nd joint: A-14/15	
1	6 1602 115	Bottle neck attachment / S 40 / green	
•	For direct dosing	from reagent bottles with S 40 thread (Merck $)$	
	Material:	PFA/PP	
			C.S.

890 Titrando

---- 57

Qty	y. Order no.	Description	
1	6.1608.023	Amber glass bottle / 1000 mL / GL 45	
	For Exchange Un	its. Bottle for auxiliary solutions.	
	Material:	Amber glass	
	Width (mm):	96	
	Height (mm):	223	
	Volume (mL):	1000	
			Comments of the

1	6.1608.030	Round glass bottle / 1000 mL / GL 45	
	To Exchange Unit	S.	e
	Material:	Clear glass	1
	Height (mm):	223	
	Volume (mL):	1000	
			10 HOT SOU
			100 mi - 700 600
			500

1	6.1609.010 Ad	Isorber tube with tubing connection	6
	Cover with tubing co	nnection.	
	Height (mm):	155	
	Inner diameter (mr	n): 32	
	SGJ size:	B-14/15	010°6091.9

I

10 Accessories

Qty.	Order no.	Description	
1	6.1623.000	Overflow protection for Ti Stan	d
	For 803 Ti Stand		
1	6.1801.120	PVC tubing / 4 mm / 6 mm / 2 n	n
	Material [.]	PVC (transparent)	
	Outer diamete	(mm): 6	
	Inner diameter	(mm): 4	
	Length (m):	2	
2	6.1805.200	PTFE tubing / M8 / 0.5 m	1
	Material:	PTFE	
	Outer diamete	(mm): 4	
	Inner diameter	(mm): 3	
	Length (mm):	500	
			4
			II.
1	6.1808.050	Connector tubing nozzle M8	
	1 M8 outer threa	and 1 tubing olive. E.g. for thermos	tat jacket of
	Exchange Units a	d stability measuring instruments.	
	Material:	PVDF	
	Length (mm):	31.5	

Qty.	Order no.	Description	
1	6.1819.030	Tube	
	For M8 threaded 6.1608.0XX bott	holes. For 6.1602.1XX bottle attachments and es	
	Material:	PTFE	
	Outer diamete	r (mm): 4	
	Inner diamete	(mm): 3	
	Length (mm):	250	

1 6.1819.050 Tube

For M8 threaded holes of 6.1602.10X bottle attachments with 6.1608.030 aspiration bottle. For 703 and 803 Ti Stand, 681 Pump Unit and 679 Rancimat

Material:	PTFE
Outer diameter (mm):	4
Inner diameter (mm):	3
Length (mm):	88

2	6.1903.020	Stirring ba	r / 16 mm	
	Stirring Bar with	n magnetic core	, PTFE coated, length 16 mm	
	Material:		PTFE	
	Outer diamet	ter (inches):	4	
	Length (mm)	:	16	

Qty.	Order no.	Description
2	6.1903.030	Stirring bar / 25 mm
	Stirring Bar with	magnetic core, PTFE coated, length 25 mm
	Material:	PTFE
	Outer diamete	er (mm): 5
	Length (mm):	25
1	6.2001.060	Base-plate without support rod
	For mounting sti	rrers.
1	6.2013.010	Clamping ring
	For stand rods of	f 10 mm diameter
	Material:	Metal
	Width (mm):	20
	Height (mm):	16
1	6.2016.070	Support rod / 400 mm
	Material:	Stainless steel 18/9
	Outer diamete	er (mm): 10
	Length (mm):	400

Qty.	Order no.	Description	
1	6.2023.020	Clip for SGJ 14/15	
	Clip for SGJ 14/16	6	
	Material:	POM	San And Stranger
			T JA
2	6.2043.005	Holding clip for bottles	
	Holding clip for re	eagent bottle of Exchange Units	
1	6.2104.020	Electrode cable / 1 m / F	
	For connecting el	ectrodes with Metrohm plug-in head G to Metrohm	
	Instruments (sock	(et F). 1	
	Length (m).		
1	6.2151.000	Cable USB A - mini-DIN 8P	
	Cable connecting	a Titrando, USB Sample Processors, Dosing Inter-	
	face, Robotic Titre	osampler to a PC (USB connection, type A) and for	
	Length (m):	1.8	
	2		
			<u> </u>
1	6.2412.000	Glass weighing spoon	a
	0.1464.32X, 6.14	4.030 Titration vessel lid or 6.1455.31X, 165.320 Titration vessels. For Karl Fischer water	
	determinations. In	ncluding protective tube.	
	Material:	Glass	
	Length (mm).	120	15 7
			111
			10/1

Qty.	Order no.	Description	
1	6.2621.070	Hexagon key 5 mm	
	5 mm.		
	Length (mm):	80	
1	6.2621.130	Hexagon key 2 mm	
	2 mm		
2	6.2730.010	Screw nipple	
	Made of PP. Toge	ether with 6.1448.010 Septum or for mounting	
	6.1403.040 Dryir	ng tube or 6.0338.100 Electrode in 6.1414.030 and	
	6.1414.040 KF II	tration vessel lids. Used for KF titrations.	
	Materiai:	PP 11.5	
	Length (mm):	23.5	Ju.
	_eg ().		
1	6.2730.020	Septum stopper	
	With E.301.0041	O-ring. For 6.1448.010 Septum and 6.1414.030	
	Titration vessel lic	. Used for KF applications.	
	Outer diamete	r (mm): 18	
	Length (mm):	30	P

Qty.	Order no.	Description	
3	6.2730.030	Stopper	
	Nipple and E.301 sel lid. Used with Advanced IC Dilu Length (mm):	.0043 O-ring included. For 6.1414.030 titration ves- KF Titrators, VA Stands, VA Computrace and 838 tion Sample Processor 38	
1	6.2739.000	Wrench	e e
	For tightening co	nnectors	
	Length (mm):	68	2
1	6.2811.000	Molecular sieve	
	Molecular sieve. moisture indicato	Bottle containing 250 g. Pore size: 0.3 nm. Without or. For Rancimats and Karl Fischer instruments.	Hand Hand Hand Hand Hand Hand Hand Hand
1	6.3026.210	Exchange unit / 10 mL	
	Exchange unit wi and light protecti tion, antidiffusior Volume (mL):	th integrated data chip with 10 mL glass cylinder ion. PCTFE/PTFE flat stopcock, FEP tubing connec- n buret tip and standard amber glass reagent bottle. 10	
1	6.6056.131	tiamo light CD: 1 Lizenz	
	tiamo light with a 856 Conductivity	additional drivers for the 888 and 890 Titrando, the Module and the 867 pH Module.	
Qty.	Order no.	Description	
------	------------------	---------------------------------------------------	--
1	6.2122.0x0	Mains cable with C13 line socket IEC-60320-C13	
	Cable plug accor	ing to customer requirements.	
	Switzerland:	Type SEV 12 6.2122.020	
	Germany,:	Type CEE(7), VII 6.2122.040	
	USA,:	Type NEMA/ASA 6.2122.070	
1	8.890.8002EN	890 Titrando Manual	

Index

Numbers/Symbols

Dosimat	14
Dosino	14
Dosino	14
Stirrer	15
Ti Stand	15
Ti Stand	15
Dosimat	14
	Dosimat Dosino Dosino Stirrer Ti Stand Ti Stand Dosimat

В

-	
Balance	18
Bluetooth adapter	
Connecting	23

С

Computer
Connecting 11
Connecting
Balance 18
Barcode reader 21
Bluetooth adapter 23
Computer 11
Dosing devices 14
MSB devices 13
PC keyboard 20
Printer
Remote box 16
Stirrer or titration stand 15
Touch Control 10
USB hub 22
Connecting a barcode reader 21
Connector
Controller 8
MSB 2, 8
Polarizable electrode
USB 2, 8
Contact pins
Controller cable 6.2151.000 11
Coupling

D

Data chip 3	30
Device software	
Update	2
Dosing command	2
Dosing devices	
Connecting 1	4
Driver	
Install 1	1

E	
Electrostatic ch	arge 5
Exchange unit	26, 31, 30

G

GLP	29
Guarantee	42
Guide openings	7

Installation	
Driver 11	I
Instrument description 2	2
Instrument type 8	3

κ

Keyboard		20
----------	--	----

.ED								
	On						7,	30
	Status	7,	27,	30,	31,	30		

М

Magnetic Stirrer	
Connecting 15	5
Mains connection 8	3
Mains voltage 5	;
Maintenance 28	3
MEAS 2)
Measuring interface	3
Measuring mode 2)
MEAS 2)
Metrohm Serial Bus MSB, see also	
"MSB" 13	3
MSB	
Connecting devices 13	3
MSB.Connector 2)

0

Overview of the instrument 7

Ρ

PC keyboard	20
Pins	32
Piston jammed	30
Printer	18

Quality Management 29

R

Rem	ote	
I	nput	33
I	nterface	32
(Output	33
Rem	ote box	
(Connecting	16
F	Pin assignment	32

S

Safety instructions	4
Serial number	8
Service	4
Service Agreement 21	9
Spindle	7
Stirrer	
Connecting 1!	5
Stopcock jammed 30	0
System test 10	0

Т

Titrando system	1
Titration mode	2
KFT	2
Titration stand	15
Connecting	15
Touch Control	
Connecting	10

U

Update	
Device software	. 2
USB connector	. 2
USB hub	
Connecting	22

V

w

Warranty	42
Windows	12